2009/JUL/23 at 5:13 p.m. Doc number: M934627A001 Rev. A [NO17]

Medtronic

DLP®/Gundry®

Retrograde Coronary Sinus Perfusion Cannula with Manual-Inflate Cuff
Canule de perfusion rétrograde du sinus coronaire avec ballonnet à remplissage manuel
Retrograde Koronarsinus-Perfusionskanüle mit manuell befüllbarem Ballon
Cánula de perfusión retrógrada del seno coronario con balón de inflado manual
Canule voor retrograde perfusie van de sinus coronarius met handmatig op te blazen ballon
Cannula per perfusione retrograda del seno coronarico con palloncino a gonfiaggio manuale
Kanyle for retrograd perfusjon av sinus coronarius med manuelt oppblåsbar mansjett
Sepelpoukaman retrogradinen perfuusiokanyyli ja manuaalisesti täytettävä mansetti
Kanyl med manuellt uppblåsbar kuff för retrograd perfusion av sinus coronarius
Κάνουλα ανάδρομης έγχυσης στεφανιαίου κόλπου με χειροκίνητα διογκούμενο αεροθάλαμο (cuff)
Kanyle til retrograd perfusion af sinus coronarius med manuelt oppustelig ballon
Kaniula do wstecznej kardioplegii z możliwością ręcznego napełniania balonu
Cânula de perfusão retrógrada do seio coronário com balão de insuflação manual
Manuel Şişirilen Manşetle beraber Retrograd Koroner Sinüs Perfüzyon Kanülü
Retrográdní kanyla pro perfuzi koronárního sinu s ručně nafukovací manžetou
Kézzel felfújható mandzsettával ellátott, sinus coronariusba vezethető retrográd perfúziós kanül
Канюля для ретроградной перфузии коронарного синуса с накачиваемой вручную манжетой

Instructions for Use • Mode d'emploi • Gebrauchsanweisung • Instrucciones de uso
Gebruiksaanwijzing • Istruzioni per l'uso • Bruksanvisning • Käyttöohjeet
Bruksanvisning • Οδηγίες χρήσης • Brugsanvisning • Instrukcja użytkowania
Instruções de utilização • Kullanım Talimatları • Návod k použití • Használati útmutató
Инструкция по эксплуатации

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DLP®/Gundry®

Retrograde Coronary Sinus Perfusion Cannula with Manual-Inflate Cuff

Product Description

The cannula consists of a wirewound silicone cannula body with a beveled tip (6-Fr models have a non-wirewound cannula body). Two side holes are present near the tip. The back of the cannula body terminates in a locking female luer. A pressure monitoring line is an integral part of the cannula body, beginning at the tip and terminating in a locking female luer fitting or 3-way stopcock with a locking female luer fitting. An inflatable balloon is located at the distal beveled tip. Models 94006, 94010, and 94015 do not have a stylet. The inflation assembly is located at the back of the cannula body and contains a female slip luer and a one-way valve assembly.

Sterile, nonpyrogenic, single use.

Indications for Use

The cannula is intended for use during cardiopulmonary bypass surgery up to six hours or less for the delivery of cardioplegia retrogradely through the coronary sinus.

Contraindications

This device is not intended for use except as indicated above.

Warnings

Extreme caution should be exercised while introducing the catheter into the coronary sinus. Do not force the cannula into the coronary sinus as this may cause vessel damage.

Do not use liquid to inflate the balloon on the 6-Fr models. Due to the incompressibility of liquids, overdistention of the coronary sinus may occur. Only air should be used to inflate the balloon. On the 6-Fr models do not exceed an inflation volume of 2 cc of air as distal balloon and/or coronary sinus damage may occur.

DO NOT OVERINFLATE THE BALLOON. Overinflation may cause coronary sinus injury. Typically, a volume of 1 cc to 2 cc (6-Fr model), 1 cc to 3 cc (10-Fr model), 3 cc to 5 cc (13-Fr and 15-Fr models), or 3 cc to 8 cc (elongated balloon models) of air is sufficient. The actual volume needed may vary as a function of cannula placement and coronary sinus anatomy. (Saline may be used, **with all but the 6-Fr models**, but may require slightly less volume to effect a seal due to its incompressible nature.) Adequacy of balloon inflation and position should be confirmed manually or visually prior to cardioplegia administration and as evidenced by a pressure wave through the pressure monitoring port.

The pressure within the coronary sinus should be monitored while administering the cardioplegia. DO NOT EXCEED 50 mmHg SINUS PRESSURE.

Precautions

Note: Proper surgical procedures and technique are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Surgeons must evaluate the appropriateness of the procedure based on their own medical training and experience, and the type of surgical procedure being performed.

Do not exceed recommended balloon inflation volume or recommended sinus pressure within the coronary sinus.

Do not inflate the balloon until it has been properly placed within the coronary sinus.

Squeezing the blue sentinel balloon on the 6-Fr models will further inflate the distal balloon and increase its diameter. Caution should be used not to overextend an already inflated distal balloon by inadvertently compressing the blue sentinel balloon on the 6-Fr models.

Balloon should be completely evacuated before insertion.

This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or create a risk of contamination of the device, which could result in patient injury, illness, or death.

Adverse Effects

This device, as do all extracorporeal blood system devices, has possible side effects which include, but are not limited to, infections, blood loss, thrombus formation, embolic events, and dislocation. Vessel damage and complications at the puncture site may occur if the Instructions for Use are not followed.

Instructions for Use

1. Inspect the package and product for damage and expiration date. If undamaged and unexpired, open the package and transfer the cannula onto the sterile field, utilizing an aseptic technique.
2. Test the balloon by inflating with sterile saline or air (6-Fr models should only use air), check for balloon inflation and leaks. After confirming proper performance, deflate balloon by withdrawing the fluid or air.

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3. To facilitate cannula placement into the coronary sinus, slip the stylet luer snugly into the cannula luer. The solid stylet is preformed and should require no additional forming for proper insertion. The guidewire stylet may be bent into a “J” or “hockey stick” configuration to ease cannulation. Do not bend into an acute angle.
4. **6-Fr model 94106 only:** To facilitate the removal of the stylet after placement, it should be wetted with sterile saline. Fill a syringe with sterile saline. Insert it into the stylet handle. Flush sterile saline through the cannula. Disconnect the syringe from the stylet handle.
5. Cannulate the aorta as per routine following heparinization. If an open atrium procedure is being used, insert the cannula with direct visualization of the coronary sinus ostium. If a closed atrium or blind insertion technique is being used, place a small 1-cm purse string in mid to upper right atrium. Make a small atriotomy through the purse string. Do not attempt to force the cannula into the coronary sinus. Reposition if resistance is encountered.
6. **Blind Insertion Method #1:** With the surgeon on the patient’s right side, insert cannula through RA purse string with the curve of the cannula facing superiorly. With the left hand controlling the cannula, place the right index finger at the junction of the IVC and AV groove inferiorly. Place the tip of the cannula against the finger, then guide through the coronary sinus and position in the coronary sinus per surgeon’s preference (usually near the left atrial appendage). The stylet is slowly withdrawn. Egress of dark pulsatile blood further confirms proper position.

 Blind Insertion Method #2: With the surgeon on the patient’s left side, follow the above steps, except that the cannula is generally held in the right hand and guided into the coronary sinus with the index finger of the left hand as above.
7. Once in the proper position, clamp proximal tubing above the wirewound section (where applicable). Close the pressure monitoring line with a 3-way stopcock. Tighten purse string and loosely tie the suture around the tourniquet and cannula suture collar to secure. Ensure that adequate cannula length is within the right atrium.

 DO NOT OVERTIGHTEN SUTURE as balloon and pressure monitoring lines may be occluded.
8. Insert appropriate venous return cannula. Confirm no dislodgment of Retrograde Coronary Sinus Perfusion cannula.
9. Place patient on bypass. During infusions, make sure that either the left ventricle or aorta is vented.
10. Attach the female luer or stopcock of the pressure monitoring line to a pressure monitoring system and assure that all air is aspirated out of the line.

 Note: The distal tip may impinge on the venous wall, making aspiration difficult. Pressure on the catheter tip with a finger eliminates this problem, should it occur.
11. Unclamp the cannula body and connect to the primed cardioplegia delivery line, assuring that no air is present within the tubing or cannula. Open pressure monitoring line stopcock.
12. After crossclamping aorta, inflate distal balloon with a volume of 1 cc to 2 cc (6-Fr model), 1 cc to 3 cc (10-Fr model), 3 cc to 5 cc (13-Fr and 15-Fr models), or 3 cc to 8 cc (elongated balloon model) of air and gradually begin retrograde cardioplegia flow. Observe coronary sinus pressure. Additional volume may be used if an adequate seal has not been attained by the minimum volume levels indicated. However, caution must be exercised not to damage the sinus due to overinflation of the balloon. An adequate seal is best confirmed by an appropriate coronary sinus pressure wave. Further confirmation is afforded by observing red coloration of the coronary veins on the anterior surface of the heart and the corresponding blue coloration of the coronary arteries if blood cardioplegia is used.

 Caution: Do not exceed 50 mmHg pressure within the coronary veins, except when the heart is beating and balloon is inflated. This condition will cause a left ventricular pressure tracing to appear on the pressure monitor.
13. Cardiac arrest usually occurs within 30 seconds to 1 minute after cardioplegia is started. If arrest has not occurred, recheck catheter position and vein coloration. Flows may be increased as long as the pressure is under 50 mmHg.
14. Keep balloon inflated between infusions. Prior to removing aortic crossclamp, deflate the balloon.
15. After completion of the procedure, remove the cannula and close the atrium and/or tie the purse string suture.

 Note: In the event of cannula dislodgment, remove cannula, reinsert the stylet making certain that an airtight fit exists between cannula and stylet, otherwise air can be sucked into the venous lines potentially creating an air lock. Reinsert cannula following step 6, except that the venous cannula is now in place.

The following disclaimer of warranty applies to United States customers only:

Disclaimer of Warranty

ALTHOUGH THE RETROGRADE CORONARY SINUS PERFUSION CANNULA WITH MANUAL-INFLATE CUFF, HEREAFTER REFERRED TO AS “PRODUCT” HAS BEEN MANUFACTURED UNDER CAREFULLY CONTROLLED CONDITIONS, MEDTRONIC HAS NO CONTROL OVER THE CONDITIONS UNDER WHICH

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THIS PRODUCT IS USED. MEDTRONIC, THEREFORE DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. MEDTRONIC SHALL NOT BE LIABLE TO ANY PERSON OR ENTITY FOR ANY MEDICAL EXPENSES OR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, FAILURE OR MALFUNCTION OF THE PRODUCT, WHETHER A CLAIM FOR SUCH DAMAGES IS BASED UPON WARRANTY, CONTRACT, TORT OR OTHERWISE. NO PERSON HAS ANY AUTHORITY TO BIND MEDTRONIC TO ANY REPRESENTATION OR WARRANTY WITH RESPECT TO THE PRODUCT.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

The following disclaimer of warranty applies to customers outside the United States:

Disclaimer of Warranty

ALTHOUGH THE RETROGRADE CORONARY SINUS PERFUSION CANNULA WITH MANUAL-INFLATE CUFF, HEREAFTER REFERRED TO AS "PRODUCT" HAS BEEN CAREFULLY DESIGNED, MANUFACTURED AND TESTED PRIOR TO SALE, THE PRODUCT MAY FAIL TO PERFORM ITS INTENDED FUNCTION SATISFACTORILY FOR A VARIETY OF REASONS. THE WARNINGS CONTAINED IN THE PRODUCT LABELING PROVIDE MORE DETAILED INFORMATION AND ARE CONSIDERED AN INTEGRAL PART OF THIS DISCLAIMER OF WARRANTY. MEDTRONIC, THEREFORE, DISCLAIMS ALL WARRANTIES, BOTH EXPRESS AND IMPLIED, WITH RESPECT TO THE PRODUCT. MEDTRONIC SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES CAUSED BY ANY USE, DEFECT, OR FAILURE OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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All Rights Reserved
Medtronic of Canada Ltd.
Medtronic, Inc.
© Medtronic, Inc. 2009