MEGAGEN

MiNi Internal Fixture

www.imegagen.com

MegaGen Implant Co.,Ltd

45, Sechen-ro 7-gil, Dasa-eup, Dalseong-gun, Daegu 42921, Republic of Korea
Tel. 82-53-222-2828/ Fax. 82-53-289-3420

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● Do not re-use

● Do not resterilize

● Caution, consult accompanying documents

● Consult instructions for use

● Temperature limitation

● Do not use if package is damaged

● Manufacturer

Rx Only ● Caution: Federal law(USA) restricts this device to sale by or on the order of a licensed dentist

I F U - I F M - R 0 8 / 2 1 - 1 2 - 0 7

English

1. Product Summary
This product is a dental implant which is made of pure titanium, and is implanted into a mandible (the lower jaw bone), or a maxillary jaw bone (the upper bone) as a root of the tooth by the dentist to restore a chewing function of partially or fully edentulous patients. There are various sizes and shapes to implant into the anterior alveolar bone for different conditions of dental loss.
It is used to in the following situations and with the clinical protocols:
- Delayed loading.
- Immediate loading when good primary stability is achieved and with appropriate occlusal loading.

2. Intended Use
Dental medical device which is implanted into the anterior alveolar bone of the upper or lower jaw bone to support prosthetic device such as dental prosthesis, and to restore a chewing function of patient who has defects at the anterior tooth.

3. Indication for Use [USA, Canada]
The MiNi Internal Implant System is intended for two-stage surgical procedures in the following situations and with the following clinical protocols:
- The intended use for the 3.0 mm diameter MiNi implant is limited to the replacement of maxillary lateral incisors and mandibular incisors.
- Immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- It is intended for delayed loading.

4. Instruction for Use & Procedure
(1) Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.
1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
4) Should check the sealing condition of the packaging, which is completely sterilized by gamma-irradiation, to ensure that the product is completely sealed and remains the sterilization status.
5) Check the bone shape and quality, volume through X-ray before implantation and select the appropriate fixture size and insertion position.
6) Use ISQ measurements to check a bone quality.
7) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expiared date, and etc. after opening the package.
8) It is recommend to choose an implant fixture as wide as possible as determined by the patient's oral condition.
9) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(2) Applications
The following explains roughly how to use and operate the product. For a detailed procedure, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.
1) Drill into the alveolar bone in order to place an implant fixture into the planned location with various surgical instruments. At this time, revolutions per minute (RPM) of the drill should be adjusted according to the condition of the bone and the kinds of surgical instruments. (Approximately 800 ~ 1200rpm) Irrigate sufficiently the saline solution onto the area to prevent the necrosis of bone by heat. * Tap drill should be drilled with lower RPM (Approximately 15 ~ 30 rpm)
2) After making an appropriate hole for a fixture, confirm the kinds of fixture and size described on the label, and take out the transparent plastic container included ampule from the exterior box.
* Measure the depth of the hole formed by drilling with the depth gauge to match the length of fixture.
* Use the various indicators and gauge types to determine the hole's condition to determine whether the hole's direction or location is correct as the operator planned.
3) Clasp the convex part downwardly directed of transparent plastic container by one hand. Then hold the arrow mark on the label sticker by another hand and take a medical paper off slowly to open the plastic container. After that, take out the sealed ampule.
4) The cap on the ampule is designed as a 'Snap-on' type, which can be easily opened by one-hand and one touch. Clasp the ampule by one hand and then push the cap up by the thumb to separate the cap from the ampule at a time.
5) Check whether the fixture is fixed on a Holder (holder type) or a Titanium leveler (leveler type). For holder type, check the Mount is connected with the fixture. For leveler type, check the Mount is not connected with the fixture. In case the fixture is external type and not connected with the Mount, you should connect the Mount to the fixture.
6) Take out carefully the fixture from the ampule by connecting Handpiece Connector to the fixture or Mount.
7) Position the fixture near the drilled site paying attention during movement that the fixture does not separate from Handpiece Connector & does not become contaminated by foreign bodies.
8) A fixture is implanted into the bone to planned depth by turning (15~30rpm) a Handpiece clockwise with 40-45Ncm torque.
9) If the fixture is not completely implanted into the bone due to the hard bone density, separate Handpiece Connector from the fixture, then finish implanting by using a Ratchet Wrench.
※ When using a Ratchet Wrench, it is recommended to torque under 80N/cm. If the torque is over 80N/cm, there is a risk of surgical failure due to excessive torque which may cause the structure deformation of inside of fixture and bone damage such as bone necrosis, partially bone fracturing.
10) After implantation of the fixture, finish the operation by connecting Cover Screw to the fixture with torque under 10Ncm by using Hand Driver to protect the inside of fixture from the foreign bodies, and then suture the surgical site.
: In the case of using Healing Abutment instead of Cover Screw, use the same torque as above.
* Abnormal surgery methods or complicated surgical techniques may increase the risk of surgical failure.
(3) Post-operative care
1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
* Minimum 6 weeks healing period is required if surrounding bone is abundant and the bone quality is good. Minimum 12 weeks healing period is required in case of cancellous bone condition.
2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
* The time of osseointegration is different depending on the bone quality of the patient, but in general, 6 months(for the maxillary), 3 months(for mandibular) are required.
3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
5) Before operating the next stage of surgery, the clinician should check the fixture's progress of Osseo-integration by x-ray image or diagnosis such as percussion reaction.
6) In case of implant failure, it must be visited immediately for treatment.

5. Adverse Effects
(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, Inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

6. Contraindication
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

7. Preservation & Handling
(1) The package hermetically sealed and sterilized with gamma radiation should be stored at low humidity and room temperature condition.
(2) The product must be opened just prior to surgery.
(3) Surgical instruments verified for safety as a medical device should be used.
(4) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(5) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(6) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(7) If the package that is not opened by dentist is defective request a replacement from manufacture.

8. Caution ⚠
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
(4) When opening the fixture, open the outer wrapper with the help of an assistant. The ampoule with the fixture packaged should be opened by the surgeon and mounted on the surgical instrument
(5) Care should be taken to avoid damaging the conical connections due to excessive force when the fixture is embedded.
(6) If the opened product or package is damaged, it is difficult to maintain the sterilization state. Discard even if it is not used.
(7) if the Pick Up is not done properly, care must be taken that the Fixture is separated and penetrated into the body during the procedure.
(8) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering t[illegible] or if the fixture placement angle and the prosthetic angle are different.
(9) Too high a density when using bone graft materials may result delayed bone fusion due to poor blood supply.

[Patient Notice]
· Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect hexamedin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.

· Care After Final Treatment
- Keep your mouth clean at all times.
- Brush your teeth, interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Please do scaling once every 6 to 12 months.

※ Warnings ----------
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.
(4) This is a single-use product and should not be reused.
(5) Risks re-using single-use medical device
1) Re-use is not validated and it carries risks such as transmission of infectious diseases.
2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
3) Failure of osseo-integration may occur.
(6) The manufacturer does not have any responsibility for the product re-sterilized by users.
(7) The contaminated products must be discarded since it cause infection and surgery fail.
(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.
(9) This product can be used in countries where the product was approved by their government or related agencies. Product availability may vary by country due to registration process.
(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

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Europe representative : MDSS GmbH

Schiffgraben 41 30175 Hannover, Germany

- CE conformity marking
- Catalogue number
- Serial number
- Batch code
- Date of manufacture
- Use by date
- Sterilized using irradiation
- Do not re-use
- Do not resterilize
- Caution, consult accompanying documents
- Consult instructions for use
- Temperature limitation
- Do not use if package is damaged
- Manufacturer
- Authorised representative in the european community
- Caution: Federal law(USA) restricts this device to sale by or on the order of a licensed dentist

IFU-IFM-R08/ 21-12-07

English

1. Product Summary

This product is a dental implant which is made of pure titanium, and is implanted into a mandible (the lower jaw bone), or a maxillary jaw bone (the upper bone) as a root of the tooth by the dentist to restore a chewing function of partially or fully edentulous patients. There are various sizes and shapes to implant into the anterior alveolar bone for different conditions of dental loss.
It is used to in the following situations and with the clinical protocols:
- Delayed loading.
- Immediate loading when good primary stability is achieved and with appropriate occlusal loading.

2. Intended Use

Dental medical device which is implanted into the anterior alveolar bone of the upper or lower jaw bone to support prosthetic device such as dental prosthesis, and to restore a chewing function of patient who has defects at the anterior tooth.

3. Indication for Use [USA, Canada]

The MiNi Internal Implant System is intended for two-stage surgical procedures in the following situations and with the following clinical protocols:
- The intended use for the 3.0 mm diameter MiNi implant is limited to the replacement of maxillary lateral incisors and mandibular incisors.
- Immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge.
- It is intended for delayed loading.

4. Instruction for Use & Procedure

(1) Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.

1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
4) Should check the sealing condition of the packaging, which is completely sterilized by gamma-irradiation, to ensure that the product is completely sealed and remains the sterilization status.
5) Check the bone shape and quality, volume through X-ray before implantation and select the appropriate fixture size and insertion position.
6) Use ISQ measements to check a bone quality.
7) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expiared date, and etc. after opening the package.
8) It is recommend to choose an implant fixture as wide as possible as determined by the patient's oral condition.
9) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(2) Applications
The following explains roughly how to use and operate the product. For a detailed procedure, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.
1) Drill into the alveolar bone in order to place an implant fixture into the planned location with various surgical instruments. At this time, revolutions per minute (RPM) of the drill should be adjusted according to the condition of the bone and the kinds of surgical instruments. (Approximately 800 ~ 1200rpm) Irrigate sufficiently the saline solution onto the area to prevent the necrosis of bone by heat. * Tap drill should be drilled with lower RPM (Approximately 15 ~ 30 rpm)
2) After making an appropriate hole for a fixture, confirm the kinds of fixture and size described on the label, and take out the transparent plastic container included ample from the exterior box.
* Measure the depth of the hole formed by drilling with the depth gauge to match the length of fixture.
* Use the various indicators and gauge types to determine the hole's condition to determine whether the hole's direction or location is correct as the operator planned.
3) Clasp the convex part downwardly directed of transparent plastic container by one hand. Then hold the arrow mark on the label sticker by another hand and take a medical paper off slowly to open the plastic container. After that, take out the sealed ample.
4) The cap on the ample is designed as a 'Snap-on' type, which can be easily opened by one-hand and one touch. Clasp the ample by one hand and then push the cap up by the thumb to separate the cap from the ample at a time.
5) Check whether the fixture is fixed on a Holder (holder type) or a Titanium leveler (leveler type). For holder type, check the Mount is connected with the fixture. For leveler type, check the Mount is not connected with the fixture. In case the fixture is external type and not connected with the Mount, you should connect the Mount to the fixture.
6) Take out carefully the fixture from the ample by connecting Handpiece Connector to the fixture or Mount.
7) Position the fixture near the drilled site paying attention during movement that the fixture does not separate from Handpiece Connector & does not become contaminated by foreign bodies.
8) A fixture is implanted into the bone to planned depth by turning (15~30rpm) a Handpiece clockwise with 40~80Ncm torque.
9) If the fixture is not completely implanted into the bone due to the hard bone density, separate Handpiece Connector from the fixture, then finish implanting by using a Ratchet Wrench.
※ When using a Ratchet Wrench, it is recommended to torque under 80N/cm. If the torque is over 80N/cm, there is a risk of surgical failure due to excessive torque which may cause the structure deformation of inside of fixture and bone damage such as bone necrosis, partially bone fracturing.
10) After implantation of the fixture, finish the operation by connecting Cover Screw to the fixture with torque under 10Ncm by using Hand Driver to protect the inside of fixture from the foreign bodies, and then suture the surgical site.
* In the case of using Healing Abutment instead of Cover Screw, use the same torque as above.
* Abnormal surgery methods or complicated surgical techniques may increase the risk of surgical failure.
(3) Post-operative care
1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
* Minimum 6 weeks healing period is required if surrounding bone is abundant and the bone quality is good. Minimum 12 weeks healing period is required in case of cancellous bone condition.
2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
* The time of osseointegration is different depending on the bone quality of the patient, but in general, 6 months(for the maxillary), 3 months(for mandibular) are required.
3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
5) Before operating the next stage of surgery, the clinician should check the fixture's progress of Osseo-integration by x-ray image or diagnosis such as percussion reaction.
6) In case of implant failure, it must be visited immediately for treatment.

5. Adverse Effects

(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open sinus site may result from inattentive surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

6. Contraindication

Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

7. Preservation & Handling

(1) The package hermetically sealed and sterilized with gamma radiation should be stored at low humidity and room temperature condition.
(2) The product must be opened just prior to surgery.
(3) Surgical instruments verified for safety as a medical device should be used.
(4) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(5) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(6) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(7) If the package that is not opened by dentist is defective, request a replacement from manufacture.

8. Caution ⚠

(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
(4) When opening the fixture, open the outer wrapper with the help of an assistant. The ampoule with the fixture packaged should be opened by the surgeon and mounted on the surgical instrument
(5) Care should be taken to avoid damaging the conical connection due to excessive force when the fixture is embedded.
(6) If the opened product or package is damaged, it is difficult to maintain the sterilization state. Discard even if it is not used.
(7) If the Pick Up is not done properly, care must be taken that the Fixture is separated and penetrated into the body during the procedure.
(8) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, if the fixture placement angle and the prosthetic angle are different.
(9) Too high a density when using bone graft materials may result delayed bone fusion due to poor blood supply.

[Patient Notice]
- Bleeding after the first procedure
 - Please bite a gauze about 1 hour for checking the bleeding.
 - Cold pack is required for swelling area.
 - Take painkillers if you have pain.
 - Do not exercise excessively, take a bath.
 - Eat soft foods and avoid irritating foods for a day or two.
 - Avoid smoking and don't drink alcohol for a month.
 - Do not brush the area after the first operation, disinfect by chlorhexidine gargle.
 - Be careful with the treatment area until the final prosthetic treatment ends.

- Care After Final Treatment
 - Keep your mouth clean at all times.
 - Brush your teeth, interdental toothbrush and floss after eating.
 - Avoid hard or chewy foods or chew a lot.
 - Please visit regularly to check for abnormalities.
 - Please do scaling once every 6 to 12 months.

※ Warnings
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.
(4) This is a single-use product and should not be reused.
(5) Risks re-using single-use medical device
1) Re-use is not validated and it carries risks such as transmission of infectious diseases.
2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
3) Failure of osseo-integration may occur.
(6) The manufacturer does not have any responsibility for the product re-sterilized by users.
(7) The contaminated products must be discarded since it cause infection and surgery fail.
(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.
(9) This product can be used in countries where the product was approved by their government or related agencies. Product availability may vary by country due to registration process.
(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

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45, Sechen-ro 7-gil, Dasa-eup, Dalseong-gun, Daegu 42921, Republic of Korea

Tel. 82-53-222-2628 / Fax. 82-53-289-3420

Representative in Europe
MegaImpl LLC
220090, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

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Applications

- XPEED AnyRidge Internal Fixture
- AnyOne Internal Fixture
- AnyOne External Fixture*
- ST Internal Fixture
- BLUEDIAMOND IMPLANT
- EZ Plus Internal Fixture**
- Rescue Internal Fixture**
- Rescue External Fixture**
- ExFeel Internal Fixture**
- ExFeel External Fixture**

(*) This model is not available for sale in Canada because it is not licensed by Health Canada.

(**) This model is not available for sale in Europe and Canada because it is not licensed by Health Canada and Europe.

English

1. Product Summary

This fixture is a dental implant, which is made of pure titanium, and is implanted into a mandible (the lower jaw bone), or a maxillary jaw bone (the upper bone) as a root of the tooth through the surgical operation by the dentist to restore a chewing function of partially or fully edentulous patients. There are various sizes and shapes to implant into the anterior, posterior site for various conditions of dental loss. It is used to in the following situations and with the clinical protocols:

- Delayed loading.

– Immediate loading when good primary stability is achieved and with appropriate occlusal loading.

Larger implants (diameter 5mm or more) are dedicated for the molar region.

2. Intended Use

Dental medical device, which is implanted in the upper or lower jaw bone to support prosthetic device such as dental prosthesis, to restore a chewing function of partially or fully edentulous patients.

3. Indication for Use – Only USA

ST Internal Fixture/ BLUEDIAMOND IMPLANT

The ST Internal Implant System and AnyRidge Octa 1 Implant System are intended to be surgically placed in the maxillary or mandibular arches for the purpose of providing prosthetic support for dental restorations (Crown, bridges, and overdentures) in partially or fully edentulous individuals. It is used to restore a patient's chewing function. The implants (less than Ø6.0 mm) are dedicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading. Larger implants are dedicated for the molar region and are indicated for delayed loading.

AnyOne Internal Fixture/ XPEED AnyRidge Internal Fixture/ AnyOne External Fixture

The Xpeed AnyRidge Internal System and AnyOne Internal/External Implant System are intended to be surgically placed in the maxillary or mandibular molar areas for the purpose providing prosthetic support for dental restorations (Crown, bridges, and overdentures) in partially or fully edentulous individuals. It is used to restore a patient's chewing function. Smaller implants (less than Ø6.0 mm) are dedicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading. Larger implants are dedicated for the molar region and are indicated for delayed loading.

EZ Plus Internal Fixture**

The EZ Plus Implant Systems are intended to be surgically placed in the upper or lower jaw to support prosthetic devices, such as artificial teeth, and to restore a patient's chewing function. This may be accomplished using either a single stage procedure or a single stage surgical procedure. Large angle abutments (e.g. 25°) on small diameter implants of the EZ Plus Implant System are not commercially available for the anterior region of the mouth and not intended for use in the posterior region of the mouth due to limited strength of the implant

Rescue Interna Fixturel*/ Rescue External Fixture**

The Rescue Implant System is intended to be surgically placed in the maxillary or mandibular arches for the purpose providing prosthetic support for dental restorations (Crown, bridges, and overdentures) in partially or fully edentulous individuals. These implants are intended to be used when a better implant is not available.

ExFeel Internal Fixture* / ExFeel External Fixture**

The ExFeel Dental Implant Systems are intended to be placed in the upper or lower jaw to support prosthetic devices, such as artificial teeth, and to restore a patient's chewing function. This may be accomplished using either a two stage surgical procedure or a single stage surgical procedure.

4. Instruction for Use & Procedure

(1) Pre-operative planning

For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.

1) Should check out the patient's status through the clinical, radiography photo, and model analysis.

2) Should ask a patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.

3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, we recommend that you ask for the consultation to determine the absence of abnormal conditions before preparing for the surgery.

4) Should check the sealing condition of the packaging, which is completely sterilized by gamma radiation. Also ensure that the product is completely sealed and remains the sterilization status.

5) Check the bone shape and quality, volume through X-ray before implantation and select the appropriate fixture size and insertion position.

6) Use ISQ measurements to check a bone quality.

7) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as fissure, cavity, stain, minute dust, expired date, and etc., after opening the package.

8) It is recommended to choose an implant fixture as wide as possible as determined by the patient's oral condition.

* In posterior site, the Manufacturer recommends choosing an Implant fixture diameter wider than the crown width. If a smaller fixture diameter is chosen, There is a risk of implant fracture due to patient's heavy masticatory force or abnormal eating habits.

9) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.

(2) Applications

The following explains roughly how to use and operate the product. For a detailed surgical procedure for the implant fixture and implant's components used with the products, please refer to the Prosthesis Manual of relevant Implant System issued by MegaGen Implant Co., Ltd.

1) Drill should be executed in order that the implant fixture into the planned location with various surgical instruments. At this time, revolutions per minute (RPM) of the drill should be adjusted according to the condition of the bone and the kinds of surgical instruments. (Approximately 800 ~ 1200rpm) Irrigate sufficiently the saline solution onto the area to prevent the necrosis of bone by heat. • Tap drill should be drilled with lower RPM (Approximately 15 ~ 30 rpm)

2) After making an appropriate hole for a fixture, confirm the kinds of fixture and size described on the label, and take out the transparent plastic container included ampule from the exterior box.

• Measure the depth of the hole formed by drilling with the depth gauge to match the length of fixture.

• Use the various indicators and gauge types to determine the hole's condition to determine whether the hole's direction or location is correct as the operator planned.

3) Grasp the convex part downwardly directed of transparent plastic container by one hand. Then hold the arrow mark on the label sticker by another hand and take a medical paper off slowly to open the plastic container so that; take out the sealed ampule.

[Ampule type]

4) The cap on the ampule is designed as a 'Snap-on' type, which can be easily opened by one-hand and one touch. Clasp the ampule by one hand and then push the cap up by the thumb to separate the cap from the ampule as a line.

5) Check whether the fixture is fixed on a Holder (holder type) or a Titanium Screw (leveler type). For holder type, check the Mount is connected with the fixture. For leveler type, check the Mount is not connected with the fixture. In case the fixture is external type and not connected with the Mount, you should connect the Mount to the fixture.

6) Take out carefully the fixture from the ampule by connecting Handpiece Connector to the fixture or Mount.

[Lego-type]

4) Open the upper ampule (long block) and take out the inner ample with one hand holding the ampule (short block).

5) The cap of the inner ample is designed as a 'snap-on' type that can be easily opened with one touch, grab the ampule with one hand and slide the lid with your thumb to remove the lid at once.

6) Make sure that the fixture is a mount type or no-mount type, and then put the fixture with a dedicated handpiece connector.

(No-mount type) Carefully take out the fixture from the ampule in the vertical direction after insert the dedicated handpiece connector into the internal hex / octa part of the conical connect.

(Mount type) Carefully take out the fixture from the ampule in the horizontal direction after insert the dedicated handpiece connector into the mount part.

7) Perform the fixture near the drilled site paying attention during movement that the fixture does not separate from Handpiece Connector or does not become contaminated by foreign bodies.

8) A fixture is implanted into the bone to planned depth by turning (15~30rpm) a Handpiece clockwise with 45Ncm or more torque.

* If the fixture is not completely implanted into the bone due to the hard bone density, separate Handpiece Connector from the fixture, then finish implanting by using a Ratchet Wrench.

※ When using a Ratchet Wrench, It is recommended to torque under 80N/cm. If the torque is over 80N/cm, there is a risk of surgical failure (i.e. fracture of the bone which may cause the structure deformation of inside of fixture and bone damage such as bone necrosis, partially bone fracturing.

10) After implantation of the fixture, finish the operation by connecting Cover Screw to the fixture with torque under 10Ncm by using Hand Driver. Protect the fixture carefully from the foreign bodies, and then suture the surgical site.

(Non-submerged) Using Healing Abutment instead of Cover Screw, use the same torque as above.

* Non-submerged surgery methods or complicated surgical techniques may increase the risk of surgical failure.

(3) Post-operative care

1) After surgery, healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, etc.). Oral condition, Patient's health condition of patient and clinician's close inspection on the progress are needed.

* Minimum 6 weeks healing period is required if surrounding bone is abundant and the bone quality is acceptable. Minimum 12 weeks healing period is required in case of cancellous bone condition, placing implant only, or poor bone quality.

2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).

* The healing period, osseointegration period, different depending on the bone quality of the patient, but in general, 6 months(for the mavillary), 3 months(for the mandibular) are needed.

3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.

4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.

5) Before operating the next stage of surgery, the clinician should check the fixture's progress of Osseo-integration by x-ray image or diagnosis and determine the next procedure.

(4) In case of implant failure, it must be visited immediately for treatment.

5. Adverse Effects

1) It could result in implantation failure, like bone loss, low initial stability, fixture failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to operating anesthesia.

(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, numbness, paresthesia at the operation site may result from inattentive surgical practice or poor patient overall health.

(3) Do not use the product in the area if the local area is usually a temporary symptom, but some cases are permanent.

6. Contraindications

Contraindications should be observed when selecting patients. These include, but are not limited to:

(1) Poor bone quality and inadequate bone quantity

(2) Bad oral hygiene and pathological infection at the alveolar bone.

(3) Immune disease such as hematodyscrasia or malignancy

(4) Alcohol drug abuse, heavy smoking and a mental patient.

(5) Other periodontal and surgical operation sites of patient.

(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI), and Cobalt-Chrom-Leglerung.

7. Preservation & Handling

(1) The package hermetically sealed and sterilized with gamma radiation should be stored at 15~22℃ condition.

(2) The product must be opened just prior to surgery.

(3) Surgical instruments and drugs that are not a medical device should be used.

(4) Unpacked or damaged product should be not be used and it cannot be returned to the manufacturer or distributor.

(5) When disposing of such devices, safely dispose of them in accordance with the hospital's administrative regulations regarding the disposal of medical devices.

(6) Any serious occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

(7) If the package has torn or not opened by dentist is defective, request a replacement from manufacture.

8. Caution ⚠

(1) Be fully aware of the instructions for use prior to using Megagen Implant System.

(2) Do not use the product, If the valid date (shelf life) were expired.

(3) Do not use the product contaminated by mistake of the operating surgeon.

(4) When opening the fixture, open the outer wrapper with the help of the assistant, and then transfer the fixture package should be opened by the surgeon and mounted on the surgical instrument

(5) Care should be taken to avoid damaging the conical connection due to excessive force when the fixture is embedded.

(6) If the opened product or package is damaged, it is difficult to maintain the sterilization state. Discard even if it is not used.

(7) If the Pick Up is not done properly, care must be taken that the fixture is not contaminated into the body during the procedure.

(8) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.

(9) Do not use the product if the bone graft materials may result in delayed bone fusion due to poor blood supply.

(10) Patients aged 18 and over only.

(11) Implantation should be completed the growth of the jaw bone (Partial or complete edentate patient).

[Patient Notice]

- Caution after the first procedure
- Please bite a gauze about 1 hour for checking the hemostasis.
- Cold pack is required for swelling area.
- Take painkillers If you have pain.
- Do not exercise excessively on the day.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation.
- disinfect hexamedin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.
- Care After Final Treatment
- Keep your mouth clean at all times.
- Brush your teeth, Interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Take an X-ray photo every 6 to 12 months.

※ Warnings

(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.

(2) The product must be used in clean and pleasant environment area.

(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.

(4) This is a single-use product and should not be reused.

(5) Risks re-using single-use medical device

① There is a high probability of infectious such as transmission of infectious diseases.

② Product can be contaminated in extensive malfunctions, including the defective connection, jammed adjustment, damage, and screw loosening.

③ Failure of osseo-integration may occur.

④ The manufacturer will not have any responsibility for the product re-sterilized by users.

(6) The unsterilized products must be discarded since it cause infection and surgery fail.

(8) The manufacturer does not have any responsibility for problems (using conditions, etc.) that do not proper instructions.

(9) This product can be used in countries where the product is approved by the government related agencies. Product availability may vary by country due to registration process.

(10) users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.

(11) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in injury.

※ Link to SSCP (The summary of safety and clinical performance)

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(*) This model is not available for sale in Canada because it is not licensed by Health Canada.
(**) This model is not available for sale in Europe and Canada because it is not licensed by Health Canada and Europe.

English

1. Product Summary

This product is a dental implant, which is made of pure titanium, and is implanted into a mandible (the lower jaw bone), or a maxillary jaw bone (the upper bone) as a root of the tooth through the surgical operation by the dentist to restore a chewing function of partially or fully edentulous patients. There are various sizes and shapes to implant into the anterior, posterior site for different conditions of dental loss. It is used to in the following situations and with the clinical protocols:
- Delayed loading.
- Immediate loading when good primary stability is achieved and with appropriate occlusal loading.

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3. Indication for Use – Only USA

ST Internal Fixture/ BLUEDIAMOND IMPLANT

The ST Internal Implant System and AnyRidge Octa 1 Implant System are intended to be surgically placed in the maxillary or mandibular arches for the purpose of providing prosthetic support for dental prostheses (Crown, bridges, and overdentures) in partially or fully edentulous individuals. It is used to restore a patient's chewing function. The implants placed in the 6.0 mm) are dedicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading. Larger implants are dedicated for the molar region and are indicated for delayed loading.

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EZ Plus Internal Fixture**

The EZ Plus Implant Systems are intended to be surgically placed in the upper or lower jaw to support prosthetic devices, such as artificial teeth, and to restore a patient's chewing function. This may be accomplished using either a two stage surgical procedure or a single stage surgical procedure. Large angle abutments (e.g. 25º) on small diameter implants of the EZ Plus Implant System are not commonly used but may be used in anterior region of the mouth and not intended for use in the posterior region of the mouth due to limited strength of the implant

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The Rescue Implant System is intended to be surgically placed in the maxillary or mandibular ridge for the purpose providing prosthetic support for dental restorations (Crown, bridges, and overdentures) in partially or fully edentulous individuals. These implants are intended to be used when a smaller implant is not indicated.

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1) Should check out the patient's medical history and the clinical, radiography photo, and model analysis.

2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.

3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.

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5) Check the bone shape and quality volume through X-ray before implantation and select the appropriate fixture size and insertion position.

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8) It is recommend to choose an implant fixture as wide as possible as determined by the patient's oral condition.
* In posterior site, the Manufacturer recommends choosing an Implant fixture diameter wider than the crown diameter. If a smaller diameter is chosen, There is a risk of implant fracture due to patient's heavy masticatory force or abnormal eating habits.

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(2) Applications

The following explains roughly how to use and operate the product. For a detailed description of the surgical procedure and implant system components used with the products, please refer to the Prosthesis Manual of relevant Implant System issued by MegaGen Implant Co., Ltd.

1) Drill should be vertical to the bone for the proper implant fixture into the planned location with various surgical instruments. At this time, revolutions per minute (RPM) of the drill should be adjusted according to the condition of the bone and the kinds of surgical instruments. (Approximately 800 ~ 1200rpm) Irrigate sufficiently the saline solution onto the area to prevent the necrosis of bone by heat. * Tap drill should be drilled with lower RPM (Approximately 15 ~ 30 rpm)

2) After making an appropriate hole for a fixture, confirm the kinds of fixture and size described on the label, and take out the transparent plastic container included ampule from the exterior box.
* Measure the depth of the hole formed by drilling with the depth gauge to match the length of fixture.
* Use the various indicators and gauge types to determine the hole's condition to determine whether the hole's direction or location is correct as the operator planned.

3) Grasp the convex part downwardly directed of transparent plastic container by one hand. Then hold the arrow mark on the label sticker by another hand and take a medical paper off slowly to open the plastic container so that, take out the sealed ampule.
[Ampule type]

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5) Check whether the fixture is fixed on a Holder (holder type) or a Titanium Screw (screw type). For holder type, check the Mount is connected with the fixture. For leveler type, check the Mount is not connected with the fixture. In case the fixture is external type and not connected with the Mount, you should connect the Mount to the fixture.

6) Take out carefully the fixture from the ampule by connecting Hex-type Connector to the fixture or Mount.
[Lego-type]

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6) Make sure that the fixture is a mount type or no-mount type, and take out the fixture with the dedicated handpiece connector.
(No-mount type) Carefully take out the fixture from the ampule in the vertical direction after inset the dedicated handpiece connector into the internal hex / octa part of the conical connect.
(Mount type) Carefully take out the fixture from the ampule in the horizontal direction after insert the dedicated handpiece connector into the mount part.

7) 7) Attach the fixture near the drilled site paying attention during movement that the fixture does not separate from Handpiece Connector or does not become contaminated by foreign bodies.

8) A fixture is implanted into the bone to planned depth by turning (15~30rpm) a Handpiece clockwise with 45Ncm or more torque.
* If the fixture is not completely implanted into the bone due to the hard bone density, separate Handpiece Connector from the fixture, then finish implanting by using a Ratchet Wrench.
※ When using a Ratchet Wrench, It is recommended to torque under 80N/cm. If the torque is over 80N/cm, there is a risk of surgical failure due to the excessive force which may cause the structure deformation of inside of fixture and bone damage such as bone necrosis, partially bone fracturing.

10) After implantation of the fixture, finish the operation by connecting Cover Screw to the fixture with torque under 10Ncm by using Hand Wrench. To protect the sealing capability of the fixture, Tightening torque of Cover Screw, use the same torque as above.
* Abrupter surgery methods or complicated surgical techniques may increase the risk of surgical failure.

(3) Post-operative care

1) After surgery, healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, etc.) and location. Typically, healing period of patient and clinician's close inspection on the progress are needed.
* Minimum 6 weeks healing period is required if surrounding bone is abundant and the bone quality is considered optimum. The healing period is required in case of cancellous bone condition, placing bone graft, etc. The determination is authorized to the attending physician.

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3) The healing period is considerably different depending on the bone quality of the patient, but in general, 6 months(for the mavillary), 3 months(for the mandibular) are needed.

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5) Before operating the next stage of surgery, the clinician should check the fixture's progress of Osseo-integration by x-ray image or diagnosis and determine the loading protocol.

5. Adverse Effects

(1) It could result in implantation failure, like bone loss, low initial stability, fixture failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to operating anesthesia.

(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, numbness, paresthesia, and scarring at the site may result from inattentive surgical practice or poor patient overall health.

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(1) Poor bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Uncontrolled systemic diseases such as hematodyscrasia and disorder of nutrition
(4) Alcohol drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation (Radiation) patient.
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI), and Cobalt-Chrom-Leglering.

7. Preservation & Handling

(1) The package hermetically sealed and sterilized with gamma radiation should be stored at 15~22℃ condition.
(2) The product must be opened just prior to surgery.
(3) Surgical instruments and drills used for a medical device should be used.
(4) Scratched or damaged product should be not used it and it cannot be returned to the manufacturer or distributor.
(5) When disposing of such devices, safely dispose of them in accordance with local authorities' authorities' regulations regarding the disposal of medical devices.

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45, Sechen-ro 7-gil, Dasa-eup, Dalseong-gun, Daegu 42921, Republic of Korea

Tel. 82-53-222-2828/ Fax. 82-53-289-3420

Representative in the Republic of Belarus:
MegaImplant LLC
220090, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

- REF ● Catalogue number
- SN ● Serial number
- LOT ● Batch code
- ● Date of manufacture
- ● Use by date
- STERILE R ● Sterilized using irradiation
- ● Do not re-use
- ● Do not resterilize
- ● Non-sterile
- ● Caution, consult accompanying documents
- ● Consult instructions for use
- ● Temperature limitation
- ● Do not use if package is damaged
- ● Manufacturer

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IFU-IA-R12/ 21-12-07

⊞ Applications

- ● Abutments & Accessories for AnyOne Internal / External* System
- ● Abutments & Accessories for XPEED AnyRidge Internal System
- ● Abutments & Accessories for EZ Plus Internal System*
- ● Abutments & Accessories for Rescue Internal* / External System*
- ● Abutments & Accessories for ExFeel Internal* / External System*
- ● Abutments & Accessories for MINI Internal System
- ● Abutments & Accessories for ST Internal System
- ● Abutments & Accessories for BLUEDIAMOND IMPLANT System
- ● Abutments & Accessories for AnyOne OneStage System
- ● Accessories for Intermezzo System*
- ● Accessories for Advanced Intermezzo System

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English

1) Product Summary

To recover a patient's masticatory function, it is necessary to have a dental implant which is implanted into a patient's alveolar bone as a root of a tooth and abutment which is protruded through the gum in the oral cavity. These products include various Implant Component Systems to offer a masticatory function and esthetic result very similar to natural tooth. The prosthetic component intended for temporary placement, such as the Temporary Abutment, Temporary Cylinder, Healing Cap has a maximum intra-oral use of 180-days.

(2) Intended Use/ Indication for Use

Applicable prosthetic components (except for Abutment & Accessories for MINI Internal System and Accessories for Intermezzo System, Advanced Intermezzo System) are used to fabricate appropriate prosthetic restorations that correspond to fixture implanted in the alveolar bone. These products include various sizes and shapes to select according to its intended use and its use.

Product lines are as follows:

- Products for healing or formation of soft tissue
- Products for fabricating of final prosthesis
- Products for fabricating of provisional prosthesis
- Products for impression taking
- Products for interconnection of prosthesis

※ Abutments & Accessories for MINI Internal System

Applicable prosthetic components are used to fabricate appropriate prosthetic restorations that correspond to fixture implanted in the alveolar bone. These products include various sizes to select according to its intended use.

- Products for healing or formation of soft tissue
- Products for fabricating of final prosthesis
- Products for fabricating of provisional prosthesis
- Products for impression taking

※ Accessories for Intermezzo & Advanced Intermezzo System

Applicable prosthetic components are used to fabricate appropriate prosthetic restorations that correspond to fixture implanted in the alveolar bone. The Fixture support the final prosthesis.

- Products for healing or formation of soft tissue
- Products for fabricating of provisional prosthesis
- Products for impression taking

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1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
4) Should check the sealing condition of the packaging, which is completely sterilized by gamma-irradiation, to ensure that the product is completely sealed and remains the sterilization status.
5) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, explared date, and etc. after opening the package.
6) Check the product information such as reference code, size, manufacturing date and expirtion date before starting the surgery.
7) Check for compatibility with the opponent component before oral installation.

(4) Instruction for Use & Procedure

For a detailed procedure and the features of our dental implant system components used with the products, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.

(5) Post-operative care

1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
5) In case of implant failure, it must be visited immediately for treatment.

(6) Adverse Effects

1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, Inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.
3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

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1) Low bone quality and inadequate bone quantity
2) Bad oral hygiene and pathological infection at the alveolar bone.
3) Internal diseases such as hematodyscrasia and diabetes, or mal-nutrition
4) Alcoholic, drug abuse, heavy smoking and a mental patient.
5) Other periodontal and surgical operation prohibited patient
6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

(8) Preservation & Handling

1) Cover Screw and Healing Abutment are products sterilized by gamma-irradiation.
2) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270°F) and wrapped with 30minutes drying time before use. (Sterilizer Type : Gravity Type)
3) Use only FDA cleared sterilization accessories including wrap.
4) Products made of plastic may deform in process of high pressure steam sterilization, therefore should not exceed 160℃ or pressure 0.45MPa.
5) After applying the investment, cool the investment slowly for coagulation to avoid the deformation by the delicate shrinking or expansion.
6) It should be used surgical instruments of relevant Implant system.
7) Melting temperature for gold materials: above 1481.5℃
8) Burnout temperature for plastic materials: above 332℃
9) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)
10) Screws fractures may occur if retightened five or more times.
11) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
12) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
13) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
14) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
15) If the package that is not opened by dentist is defective, request a replacement from manufacture.
16) Sterilized products should be stored at 15℃~22℃ condition, and non-sterile products should be stored at 15℃~30℃ condition.

(9) Caution ⚠

1) Be fully aware of the instructions for use prior to using Megagen Implant System.
2) Do not use the product, if the valid date (shelf life) were expired.
3) Do not use the product contaminated by mistake of the operating surgeon.
4) If the opened product or package is damaged, it is difficult to maintain the sterilization state. Discard even if it is not used.
5) If the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.
6) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.
7) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.
8) Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient)

[Patient Notice]

· Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation,disinfect hexamedin gargle.

Be careful with the treatment area until the final prosthetic treatment ends.

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※ Warnings

(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.
(4) This is a single-use product and should not be reused.
(5) Risks re-using single-use medical device
1) Re-use is not validated and it carries risks such as transmission of infectious diseases.
2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
3) Reuse between patients is prohibited.
(6) The manufacturer does not have any responsibility for the product re-sterilized by users.
(7) The contaminated products must be discarded since it cause infection and surgery fail.
(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.
(9) This product can be used in countries where the product was approved by their government or related agencies. Product availability may vary by country due to registration process.
* If users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.
(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or Image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
(11) Small diameter implants and angled abutments are not recommended for the posterior region.
※ Link to SSCP (The summary of safety and clinical performance)

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MDSS GmbH

Schiffgraben 41 30175 Hannover, Germany

Representative in the Republic of Belarus

MegaImplant LLC

220090, Minsk, Logoisky Trakt 22a, office 320;

Telephone: +375 29 134 55 55

- CE 2460 ● CE conformity marking
- REF ● Catalogue number
- SN ● Serial number
- LOT ● Batch code
- ● Date of manufacture
- ● Use by date
- STERILE R ● Sterilized using irradiation
- ● Do not re-use
- ● Do not resterilize
- NON STERILE ● Non-sterile
- ● Caution, consult accompanying documents
- ● Consult instructions for use
- ● Temperature limitation
- ● Do not use if package is damaged
- ● Manufacturer

EC REP ● Authorised representative in the european community

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- Products for healing or formation of soft tissue
- Products for fabricating of provisional prosthesis
- Products for Impression taking

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1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.

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4) Should check the sealing condition of the packaging, which is completely sterilized by gamma-irradiation, to ensure that the product is completely sealed and remains the sterilization status.

5) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, explared date, and etc. after opening the package.

6) Check the product information such as reference code, size, manufacturing date and expirtion date before starting the surgery.

7) Check for compatibility with the opponent component before oral installation.

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(5) Post-operative care

1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.

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5) In case of implant failure, it must be visited immediately for treatment.

(6) Adverse Effects

1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.

(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, Inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.

3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

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1) Low bone quality and inadequate bone quantity

2) Bad oral hygiene and pathological infection at the alveolar bone.

3) Internal diseases such as hematodyscrasia and diabetes, or mal-nutrition

4) Alcoholic, drug abuse, heavy smoking and a mental patient.

5) Other periodontal and surgical operation prohibited patient

6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

(8) Preservation & Handling

1) Cover Screw and Healing Abutment are products sterilized by gamma-irradiation.

2) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270°F) and wrapped with 30minutes drying time before use. (Sterilizer Type : Gravity Type)

3) Use only FDA cleared sterilization accessories including wrap.

4) Products made of plastic may deform in process of high pressure steam sterilization, therefore should not exceed 160℃ or pressure 0.45MPa.

5) After applying the investment, cool the investment slowly for coagulation to avoid the deformation by the delicate shrinking or expansion.

6) It should be used surgical instruments of relevant Implant system.

7) Melting temperature for gold materials: above 1481.5℃

8) Burnout temperature for plastic materials: above 332℃

9) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)

10) Screws fractures may occur if retightened five or more times.

11) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.

12) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.

13) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.

14) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

15) If the package that is not opened by dentist is defective, request a replacement from manufacture.

16) Sterilized products should be stored at 15℃~22℃ condition, and non-sterile products should be stored at 15℃~30℃ condition.

(9) Caution ⚠

1) Be fully aware of the instructions for use prior to using Megagen Implant System.

2) Do not use the product, if the valid date (shelf life) were expired.

3) Do not use the product contaminated by mistake of the operating surgeon.

4) If the opened product or package is damaged, it is difficult to maintain the sterilization state. Discard even if it is not used.

5) If the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.

6) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.

7) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.

8) Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient)

[Patient Notice]

· Caution after the first procedure

- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation,disinfect hexamedin gargle.

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· Care After Final Treatment

- Keep your mouth clean at all times.
- Brush your teeth, interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Please do scaling once every 6 to 12 months.

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(2) The product must be used in clean area.

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1) Re-use is not validated and it carries risks such as transmission of infectious diseases.

2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.

3) Reuse between patients is prohibited.

(6) The manufacturer does not have any responsibility for the product re-sterilized by users.

(7) The contaminated products must be discarded since it cause infection and surgery fail.

(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.

(9) This product can be used in countries where the product was approved by their government or related agencies. Product availability may vary by country due to registration process.

• If users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.

(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or Image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

(11) Small diameter implants and angled abutments are not recommended for the posterior region.

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Representative in the **Republic of Belarus**
Megalmplant LLC
220068, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

For details, visit website at www.megagen.com

English

1. Product Summary
TiGEN Abutment is intended to support for the patient- customized zirconia crown made by CAD/CAM system.
Scan Abutment is used for digital impression using dental scanner. And it is allowed to design, manufacture patient- specific abutment (customized abutment) using CAD/CAM.

2. Intended Use / Indication for Use
TiGEN Abutment is used to fabricate a customized prosthesis using a CAD/CAM system.
TiGEN Abutment is used to connect with a fixture and functions as the implant prosthesis.
In Scan Abutment system, an abutment is scanned by a 3D scanner to create an image file and is sent to CAD/CAM system. And design and manufacture a customized abutment. The 3D image of the Scan Abutment is used to know the implantation position and orientation of fixture.

3. Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed prior to the surgical treatment for patient safety, evaluation, and also the product with medical and dentistry's point of view.
(1) Should check out the patient's oral status through the clinical, radiographic, and photographic examinations.
(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
(3) Should check a patient's medical history to evaluate the status of patient's health condition, the various clinical tests should be performed to determine the abnormal and also the patient's ability to withstand the surgery.
(4) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormality prior to use, including the appearance, structure, expiration date etc. after opening the package.
(5) Check all the information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(6) Check whether to the opponent component before oral installation.

4. Instruction for Use & Procedure
* For a detailed procedure and the features of our dental implant system components, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.
(1) TiGEN Abutment
(a) Preparation before use
(a) Check the specifications of the abutment to be manufactured by carefully checking the order sheet from the dental clinic.
(b) Confirm whether the product is the one intended to be used by checking the product label, and also check the product packaging and product integrity.
(2) Using the product
(a) CAD/CAM system: dental laboratories
a) Scan the patient's oral condition in dental laboratories using a plaster model which is taken from impression of the dental clinic. Abutment, or scan it directly using a dental scanner at dental clinic.
Either scanned data or STL data should be delivered to patients.
b) Based on the 3D data, design the shape of the customized prosthesis using the designed CAD program.
c) Securely mount the TiGEN Abutment into the milling machine then chase it. Follow the instructions for use of the machine from the particular manufacturer.
d)
nician.
e) Completely clean the final prosthesis to remove any foreign matters, and send the final prosthesis to the dental clinic.
(b) Dental clinic
a) Install the final prosthesis into the oral cavity of the patient. If the prosthesis does not fit the patient, return it to the CAD/CAM center, dental laboratories for re-manufacturing.
b) Install the prosthesis according to the proper surgical protocol.
(3) Insert Abutment
(a) CAD/CAM center, dental laboratories
(a) Attach the Scan Abutment to the plaster model taken from the impression of the oral cavity of the patient.
(b) Scan the entire plaster model, including the Scan Abutment.
(c) Send the scanned file to the CAD program to have a customized abutment manufactured.
(2) Dental clinic
(a) Install the Scan Abutment to the fixture implanted inside the oral cavity of the patient.
(b) Scan the oral cavity of the patient, including the Scan Abutment.
(c) Send the scanned data to the CAD/CAM center, dental laboratories equipped with a CAD/CAM system, to manufacture a customized abutment.
* The angle should be within 20° from the axis of the fixture.
* Do not use any product inappropriate for or that cannot be connected to the implant joint. For implant compatibility, refer to the catalog.

5. Post-operative care
(1) After surgery, the healing period depends on each patient's clinical condition (bone quality / bone quantity), initial stability / loading conditions, oral condition, etc.; therefore the careful attention of patient and clinician's close inspection on the progress are needed.
(2) Give advice to the patient to avoid the undue excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(4) Attach the identification tag containing the Lot number, which is enclosed in the external box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
(5) In case of implant failure, it must be visited immediately for treatment.

6. Adverse Effects
(1) Could result in implantation failure, like bone loss, low initial stability, fracture, failure, infection, nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of, but not limited to severe infection or bleeding.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, necrosis, irritation, malocclusion, and an open suture site may result from careless surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary condition, but some cases are permanent.

7. Contraindications
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity.
(2) Any condition that may interfere with placement of the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or malfunction.
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Pregnancy, inflammation, fever, infection, soft tissue ulcer, necrosis, irritation, malocclusion, and an open suture site may result from careless surgical practice or poor patient overall health.
(6) Numbness from nerve damage in the facial area is usually a temporary condition, but some cases are permanent.

8, Preservation & Handling
(1) Non-sterile supplied products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and sterilization time with 30minutes drying time before use. (Sterilizer type - Gravity Type)
(2) It is the only FDA cleared sterilization accessories including wrap.
(3) It should be used surgical instruments of relevant Implant system.
(4) If the recommended sterilizing conditions are not observed, risks such as loosening, product deformation and fracture may occur.(refer caution)
(5) Screws fractures may occur if retightened five or more times.
(6) Save advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(8) When disposing of such devices, safety dispose of them in accordance with the normal operating procedures regarding the disposal of medical devices.
(9) Any serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(10) If the package that is not opened by dentist is defective, request replacement from manufacturer.
(11) Non-sterilized products should be stored at 15℃-30℃ condition.

9, Caution ⚠
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
(4) If the Pick Up is not done properly, care must be taken that the screw components are separated and penetrated into the body during the procedure.
(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle is not the prosthetic angle as per design.
* If the product is a single-use product, therefore it should not be reused as there is no scientific data available between patients.
(6) Patients aged 18 and over who have completed the growth of the skeletal system is typical or complete essential patient.

[Patient Notice]
- Please bite the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect hexamidin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.
- Care After Final Treatment
- Keep your mouth clean at all times.
- Brush your teeth, interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Please visit scaling once every 6 to 12 months.

⚠ Warnings
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and long-term damage.
(4) This is a single-use product and should not be reused.
(5) Re-use is not validated and it carries risks such as transmission of infectious diseases which could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
(6) Reuse between patients is prohibited.
(7) The manufacturer does not have any responsibility for the product re-sterilized by users.
(8) Contaminated products must be discarded since it cause infection and surgery fail.
(9) The manufacturer does not have any responsibility for problems using competitors products or improper instruments.
(10) The product is used in conjunction with components that are approved by their government or related agencies. Product availability may vary depending due to registration process.
* If users use similar products that are not manufactured or sold by Megagen, we cannot be responsible in quality and precision compatibility.
(11) TiGEN Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for MR-related heating, migration or artifact. The safety of TiGEN Abutment of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
* Link to SSCP (The summary of safety and clinical performance)

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EC REP **MDSS GmbH**
Schiffgraben 41 30175 Hannover, Germany

Representative in the Republic of Belarus
Megaimplant LLC
220069, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

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English

1. Product Summary
TiGEN Abutment is intended to support for the patient- customized zirconia crown made by CAD/CAM system.
Scan Abutment is used for digital impression using dental scanner. And it is designed to design, manufacture patient- specific abutment (customized abutment) using CAD/CAM.

2. Intended Use / Indication for Use
TiGEN Abutment is used to fabricate a customized prosthesis using a CAD/CAM system.
TiGEN Abutment is used to connect with a fixture and functions as the implant prosthesis.
In Scan Abutment system, an abutment is scanned by a 3D scanner to create 3D image file and it sent to CAD/CAM program, and manufacture a customized abutment. The 3D image of the Scan Abutment is used to know the implantation position and orientation of fixture.

3. Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed by the dentist or dental laboratory prior to evaluation, and also the product with medical and dentistry's point of view.
(1) Should check out the patient's oral status through the clinical, radiographic, photographic analysis.
(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
(3) Should check a patient's medical history to evaluate the status of patient. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the informal condition before the implant fixture surgery.
(4) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormality by visual inspection, contamination, damage, expiration and etc. after opening the package.
(5) Check all information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(6) Check whether is the opponent component before oral installation.

4. Instruction for Use & Procedure
※ For a detailed procedure and the features of our dental implant system please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.
(1) TiGEN Abutment
(a) Preparation before use
(a) The specifications of the abutment to be manufactured by dentist carefully checking the order sheet from the dental clinic.
(b) Confirm that the product is the same as the intended to be used by checking the product label, and also check the product packaging and product.
(c) Using the product
(a) Select the patient's dental laboratories.
(a) Scan the patient's oral condition in dental laboratories using a plaster model which is taken impression in dental clinic. Scan Abutment, or scan it directly using a dental scanner at dental clinic.
(b) Send the patient's oral condition scanned between patients to the dental laboratories using the designated CAD program.
(b) Based on the 3D data, design the shape of the customized prosthesis using the designated CAD program.
(c) Securely mount the TiGEN Abutment into the milling machine then manufacture the designed CAD prosthesis according to instructions for use of the machine from the particular manufacturer.
(d)
nician.
(a) Completely clean the final prosthesis to remove any foreign matters.
(b) Send the final prosthesis to the dental clinic.
(b) Dental clinic
(1) Install the final prosthesis into the oral cavity of the patient. If the prosthesis does not fit the patient, return it to the CAD/CAM center, dental laboratories for re-manufacturing.
(b) Install the prosthesis according to the proper surgical protocol.
(c) Abutment
(1) CAD/CAM center, dental laboratories
(a) Attach the Scan Abutment to the plaster model taken from the impression of the oral cavity of the patient.
(b) Scan the entire plaster model, including the Scan Abutment.
(c) Send the scanned file to the CAD program to have a customized abutment manufactured.
(2) Dental clinic
(a) Install the Scan Abutment to the fixture implanted inside the oral cavity of the patient.
(b) Scan the oral cavity of the patient, including the Scan Abutment.
(c) Send the scanned data to the CAD/CAM center, dental laboratories equipped with a CAD/CAM system, to manufacture a customized abutment.
* The single should be within 20° from the axis of the fixture.
* Do not use any product inappropriate for or that cannot be connected to the implant joint. For implant compatibility, refer to the catalog.

5. Postoperative care
(1) After surgery, the healing period depends on each patient's clinical condition (bone quality and bone quantity), initial stability, loading conditions, oral condition, etc., therefore the careful attention to patient and clinician's close inspection on the progress are needed.
(2) Give advice to the patient not to have any excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(4) Attach the identification tag containing the Lot number, which is supplied in the external box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
(5) In any case of failure, it must be visited immediately for treatment.

6. Adverse Effects
(1) It could result in implantation failure, like bone loss, low initial stability, fracture, failure, implant nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection or hemorrhage.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, tissue toxicity, necrosis, irritation, malocclusion, and an open suture site may result from careless surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary condition, but some cases are permanent.

7. Contraindications
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity.
(2) Acute infection at the implantation site or at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or malnutrition.
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) People with general and psychosomatic disorders.
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4 / Ti-6Al-4V ELI) and cobalt chromium alloy.

8. Preservation & Handling
(1) Non-sterile products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and 30minutes drying time before use. (Sterilizer type - Gravity Type)
(2) Use only FDA cleared sterilization accessories including wrap.
(3) It should be used surgical instruments of relevant Implant system.
(4) If the recommended sterilizing conditions are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog)
(5) Screws fractures may occur if retightened five or more times.
(6) Save advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(8) When disposing of such devices, safely dispose of them in accordance with normal procedures for handling infectious materials and disposal of medical devices.
(9) Any serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(10) If the package that is not opened by dentist is defective, request replacement from manufacturer.
(11) Non-sterilized products should be stored at 15℃-30℃ condition.

9. Caution
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
(4) If the Pick Up is not done properly, care must be taken that the screw/connector is separated and penetrated into the body during the procedure.
(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle is beyond the prosthetic angle correction.
(6) If the product is a single-use product, therefore it should not be reused, and the prosthetic angle use differently between patients.
(7) Patients aged 18 and over who have completed the growth of the jaw bone is typical or complete edentate patient.
[Patient Notice]
- Gauze must be the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid alcohol and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect hexadine gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.
- Care After Final Treatment
- Keep your mouth clean at all times.
- Brush your teeth, interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Please visit scaling once every 6 to 12 months.

※ Warnings
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and long-term damage.
(4) This is a single-use product and should not be reused.
(5) Re-use is not validated and it carries risks such as transmission of infectious diseases and failure.
(6) Re-use could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
(7) Reuse between patients is prohibited.
(8) The manufacturer does not have any responsibility for the product re-sterilized by users.
(9) Non-sterilized products must be discarded since it cause infection and surgery fail.
(10) The manufacturer does not have any responsibility for problems using competitors products or improper instruments.
(11) The product is a CDRH approved product and the product was approved by their government or related agencies. Product availability may vary country due to registration process.
* If users use similar products that are not manufactured or sold by Megagen, we cannot be responsible for quality and precision compatibility.
(12) This Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for MR-related heating, image artifact, and migration. The safety and effects of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
※ Link to SSCP (The summary of safety and clinical performance)

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Tel. 82-53-222-2828 / Fax. 82-53-289-3420

Representative In the **Republic of Belarus**
Magellan Logistics
220090, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

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IFU-IAZr-R07/ 21-12-07

English

1. Product Summary
ZrGEN Abutment is intended to support for the customized zirconia crown made by CAD/CAM system.
Scan Abutment is used for digital impression using dental scanner. And it is allowed to design, manufacture customized abutment using CAD/CAM.

2. Intended Use / Indication for Use
ZrGEN Abutment is used for connecting with a fixture and functions as the implant prosthesis. It is used with a cemented zirconia-based crown.
Scan abutment is scanned by a 3D scanner to create 3D image file and it is sent to a dental CAD to design and manufacture a customized abutment. The 3D image of the Scan Abutment is used to know the implantation position and direction of a fixture.

3. Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.
(1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
(3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
(4) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expiared date, and etc. after opening the package.
(5) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.
(6) Check for compatibility with the opponent component before oral installation.

4. Instruction for Use & Procedure
※ For a detailed procedure and the features of our dental implant system components used with the products, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.
(1) ZrGEN Abutment
1) Preparation before use
(a) Check the specifications of the abutment to be manufactured by carefully checking the order sheet from the dental clinic.
(b) Confirm whether the product is intended to be used by checking the product label, and also check the product packaging and product itself.
2) Using the product
(a) CAD/CAM center, Dental laboratories
a) Scan the patient's oral condition in dental laboratories using a dental plaster model which is taken impression in dental clinic with Scan Abutment, or directly using a dental scanner at dental clinic. Obtain a scanned 3D data file and open it using dental CAD program.
b) Based on the 3D data, design the shape of the customized prosthesis using the dental CAD program.
c) Using a milling machine, shape zirconia material to manufacture the customized prosthesis part. Refer to the instructions for use of the machine from the particular manufacturer.
d) Sinter the customized zirconia prosthesis part. Refer to equipment manufacturer's product manual for the conditions of sintering.
e) A well-trained dental technician perform final touch-up on the milled and sintered zirconia prosthesis part.
f) Cement the finished zirconia prosthesis part to the ZrGEN Abutment. Use an adhesive interacting strongly to zirconia and refer to the respective manufacturer's product manual for using and caution.
g) Completely clean the final prosthesis to remove any foreign matters.
h) Send the final prosthesis to the dental clinic.
(b) Dental clinic /Hospital
a) Install the final prosthesis into the oral cavity of the patient. If the prosthesis does not fit the patient, return it to the CAD/CAM center, dental laboratories for re-manufacturing.
b) Install the prosthesis according to the proper surgical protocol.
(2) Scan Abutment
1) CAD/CAM center, dental laboratories
(a) Attach the Scan Abutment into the dental plaster model taken from the impression of the oral cavity of the patient.
(b) Scan the entire dental plaster model, including the Scan Abutment by dental scanner.
(c) Send the scanned file to the CAD program to manufacture a customized abutment.
(3) Dental clinic
(a) Install the Scan Abutment to the fixture implanted inside the oral cavity of the patient.
(b) Scan the oral cavity of the patient including the Scan Abutment.
(c) Send the scanned data to the CAD/CAM center, dental laboratories equipped with a CAD/CAM system to manufacture a customized abutment.
• The angle should be within 20° from the axis of the fixture.
• Do not use any product inappropriate for or that cannot be connected to the implant joint. For implant compatibility, refer to the catalog.

5. Post-operative care
(1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
(2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
(5) In case of implant failure, it must be visited immediately for treatment

6. Adverse Effects
(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by excessive surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient's overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

7. Contraindications
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

8. Preservation & Handling
(1) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and wrapped with 30minutes drying time before use. (Sterilizer Type : Gravity Type)
(2) Use only FDA cleared sterilization accessories including wrap.
(3) It should be used surgical instruments of relevant Implant system.
(4) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)
(5) Screws fractures may occur if retightened five or more times.
(6) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(8) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(9) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(10) If the package that is not opened by dentist is defective, request a replacement from manufacture.
(11) Non-sterilized products should be stored at 15℃-30℃ condition.

9. Caution ⚠
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
(4) If the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.
(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.
(6) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.
(7) Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient).

[Patient Notice]
- Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is suggested for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods after finishing the operation for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect with Chlorhexidin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.

- Care After Final Treatment
- Keep your mouth clean at all times.
- Brush your teeth, interdental toothbrush and floss after eating.
- Avoid hard or chewy foods or chew a lot.
- Please visit regularly to check for abnormalities.
- Please do scaling once every 6 to 12 months.

※ Warnings
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.
(4) This is a single-use product and should not be reused.
(5) Risks re-using single-use medical device
1) Re-use is not validated and it carries risks such as transmission of infectious diseases.
2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
3) Reuse between patients is prohibited.
(6) The manufacturer does not have any responsibility for the product re-used by MegaGen.
(7) The contaminated products must be discarded since it cause infection and surgery fail.
(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.
(9) This product can be used in countries where the product was approved by MegaGen and related agencies. Product availability may vary by country due to registration process.
• If users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.
(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.
※ Link to SSCP (The summary of safety and clinical performance)

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English

1. Product Summary

ZrGEN Abutment is intended to support for the customized zirconia crown made by CAD/CAM system.

Scan Abutment is used for digital impression using dental scanner. And it is allowed to design, manufacture customized abutment using CAD/CAM.

2. Intended Use / Indication for Use

ZrGEN Abutment is used for connecting with a fixture and functions as the implant prosthesis. It is used with a cemented zirconia-based crown.

Scan abutment is scanned by a 3D scanner to create 3D image file and it is sent to a dental CAD to design and manufacture a customized abutment. The 3D image of the Scan Abutment is used to know the implantation position and direction of a fixture.

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(1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.

(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.

(3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.

(4) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, explared date, and etc. after opening the package.

(5) Check the product information such as reference code, size, manufacturing date and expiration date before starting the surgery.

(6) Check for compatibility with the opponent component before oral installation.

4. Instruction for Use & Procedure

※ For a detailed procedure and the features of our dental implant system components used with the products, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.

(1) ZrGEN Abutment

1) Preparation before use

(a) Check the specifications of the abutment to be manufactured by carefully checking the order sheet from the dental clinic.

(b) Confirm whether the product is intended to be used by checking the product label, and also check the product packaging and product itself.

2) Using the product

(a) CAD/CAM center, Dental laboratories

a) Scan the patient's oral condition in dental laboratories using a dental plaster model which is taken impression in dental clinic with Scan Abutment, or directly using intraoral scanner at dental clinic. Obtain a scanned 3D data file and open it using dental CAD program.

b) Based on the 3D data, design the shape of the customized prosthesis using the dental CAD program.

c) Using a milling machine, shape zirconia material to manufacture the customized prosthesis part. Refer to the instructions for use of the machine from the particular manufacturer.

d) Sinter the customized zirconia prosthesis part. Refer to equipment manufacturer's product manual for the conditions of sintering.

e) A well-trained dental technician perform final touch-up on the milled and sintered zirconia prosthesis part.

f) Cement the finished zirconia prosthesis part to the ZrGEN Abutment. Use an adhesive interacting strongly to zirconia and refer to the respective manufacturer's product manual for using and caution.

g) Completely clean the final prosthesis to remove any foreign matters.

h) Send the final prosthesis to the dental clinic.

(b) Dental clinic /Hospital

a) Install the final prosthesis into the oral cavity of the patient. If the prosthesis does not fit the patient, return it to the CAD/CAM center, dental laboratories for re-manufacturing.

b) Install the fixture and according to the proper surgical protocol.

(2) Scan Abutment

1) CAD/CAM center, dental laboratories

(a) Attach the Scan Abutment into the dental plaster model taken from the impression of the oral cavity of the patient.

(b) Scan the entire dental plaster model, including the Scan Abutment by optical scanner.

(c) Send the scanned file to the CAD program to manufacture a customized abutment.

(2) Dental clinic

(a) Install the Scan Abutment to the fixture implanted inside the oral cavity of the patient.

(b) Scan the oral cavity of the patient including the Scan Abutment.

(c) Send the scanned data to the CAD/CAM center, dental laboratories equipped with a CAD/CAM system to manufacture a customized abutment.

* The angle should be within 20° from the axis of the fixture.
* Do not use any product inappropriate for or that cannot be connected to the implant joint. For implant compatibility, refer to the catalog.

5. Post-operative care

(1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, etc.). Therefore the surgeon should pay attention of patient and clinician's close inspection on the progress are needed.

(2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).

(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.

(4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.

(5) In case of implant failure, it must be visited immediately for treatment

6. Adverse Effects

(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.

(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient's overall health.

(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.

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(1) Low bone quality and inadequate bone quantity

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(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition

(4) Alcoholic, drug abuse, heavy smoking and a mental patient.

(5) Other periodontal and surgical operation prohibited patient

(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.

8. Preseration & Handling

(1) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and wrapped with 30minutes drying time before use. (Sterilizer Type : Gravity Type)

(2) Use only FDA cleared sterilization accessories including wrap.

(3) It should be used surgical instruments of relevant Implant system.

(4) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)

(5) Screws fractures may occur if retightened five or more times.

(6) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.

(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.

(8) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.

(9) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

(10) If the package that is not opened by dentist is defective, request a replacement from manufacture.

(11) Non-sterilized products should be stored at 15℃-30℃ condition.

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(2) Do not use the product, if the valid date (shelf life) were expired.

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(4) If the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.

(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.

(6) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.

(7) Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient).

[Patient Notice]

- Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is useful for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft food and avoid chewing food for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect henceforth gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.

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(5) Risks re-using single use medical device

(1) Re-use is not validated and it carries risks such as transmission of infectious diseases.

(2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.

(3) Reuse between patients is prohibited.

(6) The manufacturer does not have any responsibility for the product re-used by others.

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(9) This product can be used in countries where the product was approved by the national related agencies. Product availability may vary by country due to registration process.

* If users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.

(10) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been evaluated heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

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Representative in the Republic of Belarus
MegaImplant LLC
220090, Minsk, Logoisky Trakt 22a, office 320;
Telephone: +375 29 134 55 55

For details, visit website at www.imegagen.com

Applications (Meg-Rhein Abutment)

- AnyOne Internal Fixture/ AnyOne External Fixture*
- XPEED AnyRidge Internal Fixture
- MiNi Internal Fixture / MiNi Overdenture Fixture**
- AnyOne Onestage Fixture
- BLUEDIAMOND IMPLANT
- ST Internal Fixture

Applications (Meg-Loc / Meg-Ball Abutment)

- AnyOne Internal Fixture/ AnyOne External Fixture*
- XPEED AnyRidge Internal Fixture
- MiNi Internal Fixture
- AnyOne OneStage Fixture
- BLUEDIAMOND IMPLANT
- ST Internal Fixture

Applications (Meg-Magnet Abutment*)

- AnyOne Internal Fixture/ AnyOne External Fixture*
- XPEED AnyRidge Internal Fixture
- AnyOne OneStage Fixture
- BLUEDIAMOND IMPLANT
- ST Internal Fixture

(a) This model is not available for sale in Europe because it is not licensed by Europe.
(*) This model is not available for sale in Canada because it is not licensed by Health Canada.
(**) This model is not available for sale in Canada and USA because it is not licensed by Health Canada and FDA.

English

1. Product Summary
To recover a patient's masticatory function, it is necessary to have a dental implant which is implanted into a patient's alveolar bone as a root of a tooth and abutment which is protruded through the gum in the oral cavity. These products include various Implant Component Systems to offer a masticatory function and esthetic result very similar to natural tooth.

2. Intended Use
The Overdenture System is intended to be used in the upper or lower jaw and used for supporting tooth replacements to restore chewing function. Intended for fully edentulous jaw retaining a tissue supported overdenture. Applicable prosthetic components are used to fabricate appropriate prosthetic restorations that correspond to fixture implanted in the alveolar bone. These products include various sizes and shapes to select according to its intended use and its use. Dental attachment is intended to be used in fixed overdenture restorations that can be attached.

3. Indication for Use [USA, Canada]
[USA]
Meg-Ball Attachment System, Meg-Loc Abutment, Meg-Magnet Abutment is intended to be used in the upper or lower jaw and used for supporting tooth replacements to restore chewing function. Intended for fully edentulous jaw retaining a tissue supported overdenture. The abutments in combination with endosseous implants are used as the foundation for anchoring tooth replacements in either jaw. The attachments are used in fixed overdenture restorations that can be attached with a snap-in system.
[Canada]
For Partial or complete edentate patient: Permanent partial or total loss of teeth dueto aging, illness or accident
- for supporting tooth replacements to restore chewing function
- for retaining a tissue supported overdenture in fully edentulous jaw

4. Pre-operative planning
For a successful treatment, the following confirmation procedures should be performed carefully to confirm the treatment plan and evaluation, and also the product with medical and dentistry's point of view.
(1) Should check out the patient's oral status through the clinical, radiography photo, and model analysis.
(2) Should ask the patient's decision in case of various treatment methods, after providing the sufficient explanation of the merit and demerit of each method and treatment schedule.
(3) Should check a patient's medical history to evaluate the status of patient's health condition. If there is any doubt about a patient's health condition, the various clinical tests should be performed to determine the absence of abnormal conditions before preparing for the surgery.
(4) Should not use the product if the package is not completely sealed or that is torn, and also should certainly check the abnormalities such as foreign body, stain, minute dust, expired date, and etc. after opening the package.
(5) Check the product information such as reference code, size, manufacturing date and expirtion date before starting the surgery.
(6) Check for compatibility with the opponent component before oral installation.

5. Instruction for Use & Procedure
For a detailed procedure and the features of our dental implant system components used with the products, please refer to the catalog of relevant Implant System issued by MegaGen Implant Co., Ltd.

Meg-Rhein

1. How to use
The following explains roughly how to use and operate the product. For a detailed procedure, please refer to the catalog or leaflet of relevant Implant System issued by Megagen Implant Co., Ltd.
(1) Insert a Retentive Cap (Black) into a Stainless Steel Housing by using Insertion Tool.

Image	Note (Retentivity)
	Black (for laboratory use)
	Yellow (Extra Soft - 0.6kg)
	Pink (Soft - 1.2kg)
	White (Standard - 1.8kg)
	Violet (Strong - 2.7kg)

(a) The retention force is different according to the color of Retentive Cap, and is correspond to the following order.
(Black<Yellow<Pink<White<Violet)
(b) The acceptable angle for connection between two implants of the Meg-Rhein Abutment and the Retentive Cap is 20°(141CAE) and 40°(TH), and the recommended number of use is 744 times, and the period of use is 6 months.
(2) Connect the Stainless Steel Housing to the Meg-Rhein Abutment which is connected with the fixture already.
(3) Paint on the top part of Stainless Steel Housing for dental positioning.
(4) Attach the denture.
(5) After detaching the denture, make holes on the denture at the point marked for dental positioning so that the Stainless Steel Housing can be seen from the denture.
(6) Settle the Stainless Steel Housing in the holes by putting impression material (Pink Resin) to set the denture completely in the oral cavity.
(7) After impression material is hardened, detach the denture from the Meg-Rhein Abutment.
(8) The denture which is combined with the Stainless Steel Housing is separated from Meg-Rhein Abutment due to the hardening of impression material.
(9) Remove the Retentive Cap (Black) from the Stainless Steel Housing by using Removal Tool.
(10) Select an appropriate Retentive Cap.
(11) Insert the selected Retentive Cap into the Stainless Steel Housing by using Insertion Tool.
(12) Attach the final denture to the Meg-Rhein Abutment, and confirm whether the installation of denture was done appropriately by using a X-ray equipment or CT scanner.

Meg-Loc

1. How to use
(1) After confirming that the gingiva is completely formed, remove the healing abutment connected to the fixture by using a hex driver.
(2) Clean the internal thread of the fixture cleanly with saline solution or distilled water.
(3) Select Meg-Loc Abutment of appropriate height according to thickness of gingiva and fasten Meg-Loc Abutment to connector with recommended torque.
(4) When the Meg-Loc Abutment is installed, apply the impression material around the impression material, then place the impression material on the personal tray, place it in the mouth, and then make the impression.
(5) Make a work model and check the occlusion with the denture.
(6) After checking the position where the Meg-Loc Abutment is to be completed in the completed denture, make a hole where Attachment will be placed, and attach a Block-out Spacer to the Meg-Log Abutment, and insert the Attachment into the hole of the denture.
(7) After injecting Resin, place it correctly, remove the block-out spacer, and install the oral prosthesis.
(8) The recommended number of use is 744 times, and the period of use is 6 months.

*Attachment is not a manufacture, but a purchase.

Meg-Ball

1. How to use
1.1 Meg-Ball Abutment
(1) After confirming that the gingiva is completely formed, use Healing Remove the abutment.
(2) Clean the internal thread of the fixture cleanly with saline solution or distilled water.
(3) Select the appropriate Meg-Ball Abutment according to the thickness of the gingiva and use a dedicated ball driver Attach the abutment to the fixture with the recommended torque.
(4) After the abutment is installed, apply impression material around the impression material, place the impression material on the personal tray, place it in the mouth, and make the impression.
1.2 Retentive Ring
(1) Remove the retentive ring in the metal housing with tweezers when the retentive ring wears due to frequent detachment of the denture located in the patient's mouth and does not function properly.
(2) Reinsert the worn Retentive Ring with Instrument like tweezer.
1.3 Metal Housing
(1) Make impression, make working model, make denture and check occlusion.
(2) Check the position where the abutment is to be fastened to the denture, make a hole to inject the resin, fasten the housing positioner at the desired angle to the abutment, and insert the metal housing.
(3) After injecting the resin, place the denture correctly, and when the denture and the metal housing are firmly attached, remove the denture from the oral cavity.
(4) Remove the housing positioner inside the denture and install it in the mouth.
(5) The recommended number of use is 744 times, and the period of use is 6 months.
1.4 Housing Positioner
(1) After making a working model and denture, check the position where the abutment is to be fastened to the denture, make a hole to inject resin, fasten the housing positioner at the desired angle to the abutment, and insert the metal housing.
(2) After injecting the resin, place the denture correctly, and if the denture and the metal housing are firmly attached, remove the denture from the mouth.
(3) Remove the housing positioner inside the denture.

Meg-Magnet

1. How to use
1.1 Meg-Magnet Abutment
(1) After confirming that the gingiva is completely formed, use Healing Remove the abutment.
(2) Clean the internal thread of the fixture cleanly with saline solution or distilled water.
(3) Select the appropriate Meg-Magnet Abutment according to the thickness of the gingiva and combine the abutment to the fixture with the recommended torque by using a hex driver.
(4) After the abutment is installed, apply impression material around the impression material, place the impression material on the personal tray, place it in the mouth, and make the impression.
1.2 Magnet
(1) Make impression, make working model, make denture and check occlusion.
(2) Check the position where the abutment is to be fastened to the denture, make a hole to inject the resin, fasten the magnetic positioner at the abutment, and insert the magnet.
(3) After injecting the resin, place the denture correctly, and when the denture and the magnet are firmly attached, remove the denture from the oral cavity.
(4) Remove the magnetic positioner inside the denture and install it in the mouth.
(5) The recommended number of use is 744 times, and the period of use is 6 months.
1.3 Magnetic Positioner
(1) After making a working model and denture, check the position where the abutment is to be fastened to the denture, make a hole to inject resin, fasten the magnetic positioner at the abutment, and insert the magnet.
(2) After injecting the resin, place the denture correctly, and if the denture and the metal housing are firmly attached, remove the denture from the mouth.
(3) Remove the magnetic positioner inside the denture.
5. Post-operative care
(1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
(2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
(5) In case of implant failure, it must be visited immediately for treatment.
6. Adverse Effects
(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.
7. Contraindication
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.
8. Preservation & Handling
(1) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and wrapped with 30minutes drying time. (Sterilizer Type : Gravity Type)
(2) In case of Magnet, the magnetic force can be lost at high temperature, so it is prohibited to autoclave and heat above 70 ℃. Alcohol sterilization is recommended if necessary.
(2) Use only FDA cleared sterilization accessories including wrap.
(3) It should be used surgical instruments of relevant Implant system.
(4) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)
(5) Screws fractures may occur if retightened five or more times.
(6) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(8) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(9) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(10) If the package that is not opened by dentist is defective, request a replacement from manufacture.
(11) Non-sterilized products should be stored at 15℃-30℃ condition.
9. Caution ⚠
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
5) If the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.
(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.
(6) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.
(7) Remove Magnet if taking MRI.
(8) Patients aged 18 and over who have completed the growth of the jaw bone (Partial or complete edentate patient)
[Patient Notice]
- Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect hexamedin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.
- Care After Final Treatment
- If the patients use a new denture, their may have a pain area by press the gingival. This is a denture adaptation process that can take as long as one to two weeks or two to three months. During this time, you should visit your dentist frequently to treat your sore spots.
- At first the foreign body is severe, pronunciation is not good, and a lot of saliva comes out. But over time, the patients can use the denture like part of body and learn how to use it. Adapt slowly.
- It is important to change your eating habits. Having a denture makes it difficult to eat like nature teeth. Start with soft foods and slowly increase the amount. The food should be chopped unconditionally and chewed on both sides evenly.
- Dentures should be removed at night. The gingival pressed by the dentures all day must be take a rest and if the patients put the dentures back in the morning, it can make pain, but this is normal.
- When storing dentures, put them in a bowl of cold water. Don't Put it in hot water or boil. And never drop it on the floor.
- The teeth must be brushed. Brush dentures and wipe them with a marching toothbrush. Toothpaste can damage dentures because it has many abrasives. So soapy water is recommended.
- Denture is what chew by gingival. Therefore, the gingival change over time, so it is important to check the hospital every six months to one year. And it is important adjust the wrong dentures at that time.

10. Warnings
(1) Dental surgery must only be performed by dentists and surgeons with training and experience in oral surgery as well as diagnosis and preoperative planning.
(2) The product must be used in clean area.
(3) Inappropriate selection of patient and surgery can cause the implantation failure and loss of alveolar bone.
(4) This is a single-use product and should not be reused.
(5) Risks re-using single-use medical device
(1) Re-use is not validated and it carries risks such as transmission of infectious diseases.
(2) It could be resulted in extensive malfunctions, including the defective connection, jammed abutment, damage, and screw loosening.
(3) Reuse between patients is prohibited.
(6) The manufacturer does not have any responsibility for the product re-sterilized by users.
(7) The contaminated products must be discarded since it cause infection and surgery fail.
(8) The manufacturer does not have any responsibility for problems using competitors' products or improper instruments.
(9) This product can be used in countries where the product was approved by their government or related agencies. Product availability may vary by country due to registration process.
* If users use similar products that are not manufactured or sold by MegaGen, there may be a problem in quality and precision compatibility.
(10) Please note to patients that their dentures may take out when eating sticky foods.
(11) The Product should not be used if packaging is open or damaged
(12) Consumption of sticky foods should be limited to the violet retentive caps, using a minimum of five Meg-Rhein Abutments to support a single overdenture. Alternatively, the dentures may be taken out when eating sticky foods.
(13) The Dental Implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of Dental Implant in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

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	● **Manufacturer**
REF	● Catalogue number
SN	● Serial number
LOT	● Batch code
	● Date of manufacture
	● Use by date
NON STERILE	● **Non-sterile**
	● Do not re-use
STERILIZE	● Do not resterilize
	● Caution, consult accompanying documents
	● Consult instructions for use
	● Temperature limitation
	● Do not use if package is damaged
Rx Only	● For USA Only (Federal law(USA) restricts this device to sale by or on the order of a licensed dentist)

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EC REP MDSS GmbH Schiffgraben 41 30175 Hannover, Germany

Representative in the Republic of Belarus MegaImplant LLC
220090, Minsk, Logoisky Trakt 22a, office 320; Telephone: +375 29 134 55 55

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Meg-Magnet

1. How to use
1.1 Meg-Magnet Abutment
(1) After confirming that the gingiva is completely formed, use Healing Remove the abutment.
(2) Clean the internal thread of the fixture cleanly with saline solution or distilled water.
(3) Select the appropriate Meg-Magnet Abutment according to the thickness of the gingiva and combine the abutment to the fixture with the recommended torque by using a hex driver.
(4) After the abutment is installed, apply impression material around the impression material, place the impression material on the personal tray, place it in the mouth, and make the impression.
1.2 Magnet
(1) Make impression, make working model, make denture and check occlusion.
(2) Check the position where the abutment is to be fastened to the denture, make a hole to inject the resin, fasten the magnetic positioner at the abutment, and insert the magnet.
(3) After injecting the resin, place the denture correctly, and when the denture and the magnet are firmly attached, remove the denture from the oral cavity.
(4) Remove the magnetic positioner inside the denture and install it in the mouth.
(5) The recommended number of use is 744 times, and the period of use is 6 months.
1.3 Magnetic Positioner
(1) After making a working model and denture, check the position where the abutment is to be fastened to the denture, make a hole to inject resin, fasten the magnetic positioner at the abutment, and insert the magnet.
(2) After injecting the resin, place the denture correctly, and if the denture and the metal housing are firmly attached, remove the denture from the mouth.
(3) Remove the magnetic positioner inside the denture.
5. Post-operative care
(1) After surgery, the healing period depends on each patient's clinical condition (bone quality, bone quantity, initial stability, loading condition, Oral condition, etc.). Therefore the careful attention of patient and clinician's close inspection on the progress are needed.
(2) Give advice to the patient for recognizing that excessive pressure such as masticatory force is not applied to the surgical site during the healing period (Osseo-integration period).
(3) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(4) Attach the identification tag containing the Lot number, which is enclosed in the exterior box, on the patient's chart, x-ray film, etc. to make it possible to trace the product.
(5) In case of implant failure, it must be visited immediately for treatment.
6. Adverse Effects
(1) It could result in implantation failure, like bone loss, low initial stability, fixation failure, irritation and nerve damage, degradation of masticatory function by careless surgical practice. Also there is a risk of death due to severe infection in severe cases.
(2) Complications such as malaise, flushing, itching, pain, edema, bleeding, inflammation, fever, infection, soft tissue ulcer, feeling of irritation, malocclusion, and an open suture site may result from inattentive surgical practice or poor patient overall health.
(3) Numbness from nerve damage in the facial area is usually a temporary symptom, but some cases are permanent.
7. Contraindication
Contraindications should be observed when selecting patients. These include, but are not limited to:
(1) Low bone quality and inadequate bone quantity
(2) Bad oral hygiene and pathological infection at the alveolar bone.
(3) Internal diseases such as hematodyscrasia and diabetes, or mal- nutrition
(4) Alcoholic, drug abuse, heavy smoking and a mental patient.
(5) Other periodontal and surgical operation prohibited patient
(6) People with allergies and hypersensitivity for titanium (CP Ti Grade 4, Ti-6Al-4V ELI) and cobalt chromium alloy.
8. Preservation & Handling
(1) Non-sterilized products should be sterilized with high pressure steam sterilization for 15 minutes at temperatures 132℃(270℉) and wrapped with 30minutes drying time. (Sterilizer Type : Gravity Type)
In case of Magnet, the magnetic force can be lost at high temperature, so it is prohibited to autoclave and heat above 70 ℃. Alcohol sterilization is recommended if necessary.
(2) Use only FDA cleared sterilization accessories including wrap.
(3) It should be used surgical instruments of relevant Implant system.
(4) If the recommended tightening torques are not observed, risks such as loosening, product deformation and fracture may occur.(refer catalog.)
(5) Screws fractures may occur if retightened five or more times.
(6) Give advice to the patient for taking good care of the patient's oral hygiene to avoid side effects.
(7) Unpacked or damaged product should be not used and it cannot be returned to the manufacturer or distributor.
(8) When disposing of such devices, safely dispose of them in accordance with the relevant authorities' regulations regarding the disposal of medical devices.
(9) If serious incident occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.
(10) If the package that is not opened by dentist is defective, request a replacement from manufacture.
(11) Non-sterilized products should be stored at 15℃-30℃ condition.
9. Caution ⚠
(1) Be fully aware of the instructions for use prior to using Megagen Implant System.
(2) Do not use the product, if the valid date (shelf life) were expired.
(3) Do not use the product contaminated by mistake of the operating surgeon.
5) if the Pick Up is not done properly, care must be taken that the product is separated and penetrated into the body during the procedure.
(5) Care should be taken that product deformation and fracture can easily occur if the prosthesis is manufactured without considering the patient's tooth shape and load distribution, or if the fixture placement angle and the prosthetic angle are different.
(6) The product is a single-use product, therefore it should not be reused and user not to reuse products between patients.
(7) Remove Magnet if taking MRI.
(8) Patients aged 18 and over who have completed the growth of the jaw bone
(Partial or complete edentate patient)
[Patient Notice]
- Caution after the first procedure
- Please bite a gauze about 1 hour for checking the bleeding.
- Cold pack is required for swelling area.
- Take painkillers if you have pain.
- Do not exercise excessively, take a bath.
- Eat soft foods and avoid irritating foods for a day or two.
- Avoid smoking and don't drink alcohol for a month.
- Do not brush the area after the first operation, disinfect hexamedin gargle.
- Be careful with the treatment area until the final prosthetic treatment ends.
- Care After Final Treatment
- If the patients use a new denture, their may have a pain area by press the gingival. This is a denture adaptation process that can take as long as one to two weeks or two to three months. During this time, you should visit your dentist frequently to treat your sore spots.
- At first the foreign body is severe, pronunciation is not good, and a lot of saliva comes out. But over time, the patients can use the denture like part of body and learn how to use it. Adapt slowly.
- It is important to change your eating habits. Having a denture makes it difficult to eat like nature teeth. Start with soft foods and slowly increase the amount. The food should be chopped unconditionally and chewed on both sides evenly.
- Dentures should be removed at night. The gingival pressed by the dentures all day must be take a rest and if the patients put the dentures back in the morning, it can make pain, but this is normal.
- When storing dentures, put them in a bowl of cold water. Don't Put it in hot water or boil. And never drop it on the floor.
- The teeth must be brushed. Brush dentures and wipe them with a marching toothpaste. Toothpaste can damage dentures because it has many abrasives. So soapy water is recommended.
- Denture is what chew by gingival. Therefore, the gingival change over time, so it is important to check the hospital every six months to one year. And it is important adjust the wrong dentures at that time.

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	• Manufacturer
CE 2460	• CE conformity marking
REF	• Catalogue number
SN	• Serial number
LOT	• Batch code
	• Date of manufacture
	• Use by date
NON STERILE	• Non-sterile
	• Do not re-use
STERILIZE	• Do not resterilize
	• Caution, consult accompanying documents
	• Consult instructions for use
	• Temperature limitation
	• Do not use if package is damaged
Rx Only	• For USA Only (Federal law(USA) restricts this device to sale by or on the order of a licensed dentist)
EC REP	• Authorised representative in the european community

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