# Temporary Pacemaker User Manual

## Mode 6108

Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

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Product Name: Temporary Pacemaker

Mode:6108

After-sales Service Provider:
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

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# Table of Contents

Symbols ..... 1
1 Device Overview ..... 3
1-1 Description ..... 3
1-2 Package Content ..... 3
1-3 Features ..... 4
1-4 Intended Use ..... 5
1-5 Indications ..... 6
1-6 Contraindications ..... 6
1-7 Patient Cable/Analysis Cable ..... 7
2 Warnings and Precautions ..... 8
2-1 Warnings ..... 8
2-2 Precautions ..... 13
2-3 Adverse Effects ..... 16
3 User Interface ..... 17
3-1 Dial and Buttons ..... 19
3-2 Patient Cable Socket ..... 20
3-3 Indicator Lights ..... 21
3-4 Low Battery Alert ..... 23
3-5 Battery Compartment ..... 23
3-6 Hanging Ring ..... 23
3-7 Protective Cover ..... 24
3-8 Back Clip ..... 24
4 Device Operation ..... 25
4-1 Preparation ..... 25
4-2 Power On ..... 27
4-3 Power-Off ..... 28
4-4 Synchronous Pacing Mode (SSI) ..... 28
4-5 Asynchronous Pacing Mode (SOO) ..... 29
4-6 Pulse Amplitude ..... 30
4-7 Pulse Frequency ..... 30
4-8 Sensing Sensitivity ..... 31
4-9 Run-Away Protection ..... 31
4-10 Sensing Events ..... 32
4-11 Refractory Period ..... 32
4-12 Noise Reversion ..... 32

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
4-13 Battery Status ..... 33
4-14 Battery Replacement ..... 34
4-15 Cleaning and Disinfection ..... 35
4-16 Post-Cleaning Inspection ..... 36
5 Application ..... 37
5-1 Preparation ..... 37
5-2 Sensing Peak Measurement ..... 37
5-3 Pacing Threshold Test ..... 38
6 Device Specifications ..... 39
6-1 Parameter ..... 39
6-2 Compliance with Standards ..... 44
7 Support, Disclaimer and Warranty Information ..... 51
7-1 Support ..... 51
7-2 Disclaimers ..... 52
7-3 Warranty Information ..... 53

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

Humidity limit (0%-90%RH)

General warning sign

Refer to instruction manual/ booklet

MR Unsafe

Defibrillation-Proof CF-Type Applied Part

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


Material code

Patient Cable

Press both battery compartment buttons simultaneously to open

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# 1 Device Overview

## 1-1 Description

The Lifetech Cardio 6108 Temporary Pacemaker consists of a Non-implantable Pulse Generator and a Patient Cable.

The Lifetech 6108 Temporary Pacemaker is a handheld device powered by two AA (LR6) alkaline batteries. The device is directly connected to one end of the patient cable, while the other end of the patient cable is connected either directly or via an adapter cable to the pacing lead, which in turn is connected to the myocardial tissue. The pacing system consists of four components: the temporary pacemaker, the patient cable, the pacing lead, and the myocardial tissue.

The 6108 Temporary Pacemaker offers two pacing modes: SSI (synchronous pacing) and SOO (asynchronous pacing), providing single-chamber synchronous/asynchronous pacing therapy. Its pulse rate, pulse amplitude, and sensing sensitivity are adjustable.

The 6108 Temporary Pacemaker should only be operated by trained medical professionals.

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## 1-3 Features

The functional features of the 6108 Temporary Pacemaker are summarized as follows.

**Pacing Features**

| | |
|---|---|
| **Pacing Modes** | SSI, SOO. |
| **Pulse Width** | 1.8ms |
| **Amplitude** | Adjustable, range: 0.1 – 10.0V |
| **Rate** | Adjustable, range: 30 – 180 ppm |
| **Sensitivity** | Adjustable, range: 0.5 – 20.0mV |

**Safety Features**

**Interference Reversion**: When electrical interference is sensed, such as continuous noise with a cycle length of 40 milliseconds or less (e.g., 50 Hz or 60 Hz interference), the device should revert to asynchronous pacing at the set pulse rate until the interference ceases.

**Frequency Limitation (Runaway Pro taction)**:Ensures the device's output pulse rate does not exceed 220ppm.

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**Defibrillation Protection**: The device's defibrillation protection function is designed to minimize the risk of malfunction when external defibrillation voltage is applied.

**Low Battery Warning:** When the battery becomes critically low, the amber "⏻" light begins to flash and the buzzer sounds. After the low battery warning is triggered, the device should be able to continue functioning for at least 24 hours.
When the battery compartment is opened, the amber "⏻" light also begins to flash and the buzzer sounds**.**

**Backup Power:** After the battery is removed, the backup power supply can continue to provide power for at least 55 seconds, during which the amber "⏻" light flashes**.**

**Power-off Protection:** When the device is operating normally, pressing the "⏻" button will cause the "Pace," "Sense," and "⏻" indicator lights to flash in a rolling sequence. If the "⏻" button is pressed again before the rolling indicator lights disappear, the device will shut down; otherwise, the shutdown operation will be canceled.

## 1-4 Intended Use

It is for external temporary pacing of either the atrium or the ventricle.

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1) Sick sinus syndrome
2) Sinus bradycardia
3) Atrial and/or ventricular arrhythmias
4) Complete atrioventricular block
5) Asystole
6) Bradycardia accompanied by congestive heart failure
7) Patient support, management and evaluation before the implantation of implantable pacemaker
8) Support during the replacement of implantable pacemaker
9) Cardiac complications occurring during intervention or surgery
10) Support after the cardiac surgery
11) Acute myocardial infarction complicated with cardiac conduction block

## 1-6 Contraindications

There are No known contraindications when used for temporary cardiac pacing.

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## 1-7 Patient Cable

The patient cable provided as an accessory can be used to connect with the 6108 Temporary Pacemaker for temporary pacing. The cable is approximately 1860 mm in length, with one end featuring a plug and the other end an electrode connector.
The electrode connector can be connected to pacing leads or adapter cables with a connector diameter not exceeding 3.0 mm. The plug is used to connect to the patient cable socket on the top of the 6108 Temporary Pacemaker.

For extra cables, contact Lifetech Cardio (see front page for details).

| Accessory | Order Code |
|---|---|
| Patient Cable | 0202002007 |

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## 2-1 Warnings

### Daily Use and Maintenance

**Device Maintenance:** Do not disassemble or attempt to repair the 6108 temporary pacemaker yourself. Any disassembly or repair not performed by Lifetech Cardio may result in device malfunction or compromise the safety of its use.

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**Liquid Protection:** Operators may wipe the device using a cloth dampened with clean water. Do not immerse the device in any liquid. Any ingress of liquid may compromise the device's functionality and safety.

**Fire and Explosion Prevention:** Improper use of the Temporary Pacemaker may generate electrical sparks. Therefore, it is not suitable for use in environments containing flammable or explosive gases, such as anesthetic gases, oxygen, or nitrous oxide.

**Monitoring:** No diagnostic functions

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are intended.

**Not MRI-Compatible:** Magnetic Resonance Imaging (MRI) is a medical imaging technique that uses magnetic fields to generate internal views of the body. MRI scanning should not be performed on patients using the 6108 Temporary Pacemaker, as it may induce tachycardia, cause serious bodily harm, or lead to system malfunction or damage.

**Equipment Maintenance:** Do not disassemble or repair the 6108 Temporary Pacemaker yourself. Any disassembly or repair work not performed by Lifetech Cardio may cause equipment malfunction or affect the safe use of the equipment.

**Usage:** Do not connect, place in close proximity to or stacked the 6108 Temporary Pacemaker with other devices. If close or stacked use is necessary, the system should be observed and verified to ensure proper operation under the specific configuration.
Except for transducers and cables sold by the device or system manufacturer as replacement internal components, the use of unauthorized accessories, transducers, or cables may result in increased electromagnetic emissions or reduced immunity of the device or system.

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**Disposal:** Please dispose of the 6108 Temporary Pacemaker in accordance with local regulations.

## Clinical Use

**Threshold Testing:** During pacing threshold testing, loss of capture may occur, which could lead to cardiac arrest.

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**Atrial Pacing:** In cases of atrioventricular conduction disturbances or atrial fibrillation/flutter, single-chamber atrial pacing is ineffective in stimulating the heart and should not be used as a treatment.

**Pacing Cessation:** Abrupt termination of pacing therapy may lead to cardiac arrest. Before stopping pacing, it is recommended that the operator first confirm the 6108 temporary pacemaker is in synchronous mode and then adjust the pulse rate to slightly below the patient's intrinsic heart rate.

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Pacing should only be discontinued after confirming the return of the patient's own rhythm.

**Introduced Current:** Implanted electrodes form a direct electrical connection with the patient′ s heart, meaning that any introduced current may be directly conducted to the heart, potentially causing fatal arrhythmias. Additionally, introduced current may lead to inappropriate sensing by the device, resulting in suppression of pacing pulses and loss of capture. Lastly, such current may also damage the device.

**Patient Cable:** Connect the patient cable to the temporary pacemaker first, and then connect the patient cable to the pacing lead. Using patient cables other than those specified by Lifetech Cardio may increase the likelihood of inappropriate sensing or reduce the therapeutic effectiveness of pacing pulses.

**Electrostatic Discharge (ESD):** To prevent the effects of static electricity on the device and the patient, operators should discharge any static buildup from their body before using the device, such as by touching a large, grounded conductive surface.

**Implanting electrodes:** Prevent electrode contact with conductive/damp surfaces during operation.

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**Mains-Powered Equipment:** If using equipment connected to the mains power supply (220V), ensure that it has a properly configured grounding wire.

**Electrosurgical Unit (ESU):**
Maintain as much distance as possible, typically no less than 15 cm, between the ESU current input site and the electrode tip of the pacing lead.

**Ablation:**
Maintain as much distance as possible, typically no less than 15 cm, between the ablation electrode and the electrode tip of the implanted pacing lead. Direct contact between the ablation electrode and the pacing lead tip should be avoided.

**Defibrillation:**
Maintain as much distance as possible, typically no less than 15 cm, between the defibrillation pads and the electrode tip of the implanted pacing lead. After each defibrillation, check the ECG monitor and confirm that the 6108 temporary pacemaker is functioning properly. Note that when defibrillation is performed with the 6108 Temporary Pacemaker connected, part of the defibrillation energy may be shunted through the conduction path between the pacing lead tip and the myocardium, potentially reducing the effectiveness of the defibrillation therapy.

Note: Before each defibrillation, it is recommended to disconnect the 6108

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temporary pacemaker from the patient and/or implanted pacing lead whenever possible.

**RF Communication Devices:** Keep handheld RF communication devices, including antenna cables and external antennas, as far away as possible from the patient cable of the 6108 temporary pacemaker, typically no less than 30 cm, to avoid potential device damage.

**Radio Frequency Resistance:** The 6108 Temporary Pacemaker may malfunction if used near RF transmitters due to limited signal resistance.

**Electromagnetic Interference:** Strong electromagnetic interference may cause the 6108 Temporary Pacemaker to misinterpret signals, resulting in changes to its output. In such cases, it is recommended to switch the device to SOO mode and set the pacing rate slightly higher than the patient's intrinsic heart rate.

**Patient monitoring:** Continuously monitor the patient while the temporary pacemaker is in use to ensure proper functioning of the temporary pacemaker and to provide the patient with appropriate treatment.

## 2-2 Precautions

**Lead Polarity:** Because bipolar pacing leads offer greater resistance to interference compared to unipolar leads,

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it is recommended to use bipolar pacing leads whenever possible with the 6108 temporary pacemaker.

**Simultaneous Pacing:** Simultaneous pacing refers to the condition where the 6108 temporary pacemaker operates concurrently with another independent pacing device. This may cause confusion in the pacing and sensing functions of the 6108 pacemaker. The use of the 6108 temporary pacemaker for simultaneous pacing is strictly prohibited.

**Battery Replacement:** The 6108 temporary pacemaker is powered by two AA (LR6) alkaline batteries. When the low battery indicator light turns on or the battery level indicator shows empty, the batteries should be replaced immediately. The device supports at least 55 seconds of operation without batteries to allow for battery replacement.
A new set of batteries should be used for each patient. If the device is not used at all, the batteries should be replaced every two weeks to ensure sufficient power when needed.

**Environmental Precautions:** Operating the 6108 Temporary Pacemaker in environments with ab-normal atmospheric pressure may cause device damage. To ensure proper operation, the following pre-cautions should be observed during handling, transport,

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and storage:

- The device should be used clinically by trained medical professionals.
- Operators must perform proper elec-trostatic discharge procedures before handling the device.
- Batteries should be removed when the device is not in use.
- Do not subject the device to physical damage such as mechanical impact.
- Do not place the device in environments outside the specified operating and storage temperature and humidity ranges (refer to “Device Specifications” ).
- Dispose of used batteries in accordance with local environmental regulations.
- Do not place the 6108 Temporary Pacemaker where it can be accessed by the patient.
- Do not immerse the 6108 temporary pacemaker or the patient cable in water.
- Do not sterilize the device using gamma radiation, steam, or similar methods.

If there is a sudden change in the ambient temperature or humidity where

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the 6108 temporary pacemaker is located, wait until the device temperature equalizes with the environment before use.

**Electromagnetic Compatibility (EMC) Notice:** For detailed information on electromagnetic compatibility, refer to section "6.3 Electromagnetic Compatibility Statement." Users must operate the device within the electromagnetic environment specified in that section.
Potential adverse environmental factors include, but are not limited to, those listed above. As a general principle, the device should always be operated with patient health and safety as the top priority, minimizing any environmental impact that could compromise device performance.

## 2-3 Adverse Effects

Potential Adverse Effects of Using the 6108 Temporary Pacemaker Include:

1) Abrupt interruption of temporary pacing therapy before the intrinsic rhythm is restored may result in arrhythmia.


2) In SOO mode, pacing pulses may coincide with the patient's intrinsic T wave, which may lead to tachyarrhythmia.

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
# 3 User Interface

The 6108 Temporary Pacemaker interface includes the following parts:

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(1) "Pacing" Indicator Light
(2) "Sensing" Indicator Light
(3) "Amplitude" Knob
(4) "Sensing" Knob
(5) "Rate" Knob
(6) Power Button and "Low Battery" Light
(7) Patient Cable Connector
(8) Hanging Ring
(9) Back Clip
(10) Battery Compartment
(11) Protective Cover
(12) Battery Compartment key

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## 3-1 Dial and Buttons

There are a total of three knobs and three buttons on the control panel of the 6108 temporary pacemaker. Their functions are described below.

**Knob：**

"Amplitude" Knob: This knob is used to adjust the pulse amplitude. Turning the knob clockwise increases the pacing voltage, while turning it counterclockwise decreases the pacing voltage.

"Rate" Knob: This knob is used to adjust the base pacing rate. Turning the knob clockwise increases the pacing rate, while turning it counterclockwise decreases the pacing rate.

"Sensitivity" Knob: This knob is used to adjust the sensing sensitivity value. Turning the knob counterclockwise increases the sensing threshold (i.e., decreases sensitivity), while turning it clockwise decreases the sensing threshold (i.e., increases sensitivity).

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
The “⏻” button is located at the bottom right corner of the device and is used to power the device on and off. When the device is off, pressing the “⏻” button once will turn it on. There is also a “⏻” indicator light at the same location. When the battery is low, this light will flash
To turn off the device, a two-step shutdown process is required:
After the first press of the “⏻ button, the “Pace,” “Sense,” and “⏻” indicator lights will flash in a rolling sequence. To shut down the device, press the “⏻” button again before the rolling lights stop. If the second press occurs after the lights stop flashing, the device will not shut down and will instead repeat the rolling indicator sequence.

**Battery Compartment keys**
The battery compartment keys are located on both sides of the battery compartment. To open the compartment for battery installation or removal, the user should press both buttons simultaneously.


## 3-2 Patient Cable Socket

The patient cable socket is located on the top of the 6108 Temporary Pacemaker, through which the patient cable connects to the device.

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
Please connect the patient cable as follows:

1. Align the plug of the patient cable with the guide rails on the patient cable socket, then insert the patient cable into the socket of the 6108 Temporary Pacemaker.
2. Continue inserting the cable until a "click" sound is heard.
3. Gently pull on the cable portion to ensure a secure connection. The cable should not come loose.

Disconnect the patient cable as follows:

1. Hold the 6108 Temporary Pace maker firmly to keep it stable. Using two fingers, gently pinch the locking ring on the patient cable plug and pull the cable gently in the direction away from the device.
2. Continue to pull gently until the patient cable is fully disconnected from the socket on the 6108 temporary pacemaker.

Note: Do not pull on the cable body itself to disconnect the patient cable. This may not only fail to detach the cable but could also cause damage to it.

## 3-3 Indicator Lights

The "Pacing" and "Sensing" indicator lights are located side by side at the top front of the 6108 Temporary Pacemaker. In addition, the "⏻" indicator

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light is located at the "⏻" button. During the device's power-on self-test, all three indicator lights will illuminate simultaneously.

**Pacing Indicator Light**

The "Pace" indicator light is green and located on the top left of the device. When the device is operating normally, the light flashes in sync with each pacing event—flashing once for every pacing pulse delivered.

**Sensing Indicator Light:**

The "Sensing" indicator light is white and located on the top right of the device. When the device is operating normally, the light flashes in sync with sensing events, that is flashing once each time the device detects an electrical cardiac activity outside the refractory period.

**"⏻" Indicator Light**

The "⏻" indicator light is amber and located at the same position as the "⏻" button. The 6108 Temporary Pacemaker continuously monitors battery status, and when the battery level is low, this indicator will flash. Once the "⏻" light starts flashing, please check the power supply and replace both AA (LR6) alkaline batteries promptly. After replacement, verify that the "⏻" light has stopped flashing.

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## 3-4 Low Battery Alert

The 6108 Temporary Pacemaker continuously monitors battery status. When the battery level becomes critically low, the amber "⏻" indicator light begins to flash and the buzzer sounds. After the low battery warning is triggered, the device should be able to continue operating for at least 24 hours.

## 3-5 Battery Compartment

The battery compartment is located on the side of the device and is designed to hold two AA (LR6) alkaline batteries.

**To open the battery compartment:**
Press the buttons on both sides of the battery compartment simultaneously to release and open it.

**To close the battery compartment:**
After inserting the batteries according to the indicated polarity, slide the compartment back into the device. A clicking sound indicates that the battery compartment is properly secured in place.

## 3-6 Hanging Ring

The hanging ring is located at the upper part of the back of the device. Operators can use this ring to hang the device. For example, on an IV pole or

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a hospital bed rail. To use it, pull the hanging ring outward from its slot and rotate it 180 degrees until it gently clicks into place. When not in use, please return the ring to its original position.

## 3-7 Protective Cover

The protective cover is designed to prevent accidental contact with the device's control knobs.

When adjusting parameters, slide the protective cover upward along the guide rail from bottom to top until the knob area is fully exposed, then turn the desired knob. After adjustments are complete, slide the cover back down from top to bottom along the guide rail on the front of the pacemaker to fully cover the knob area.

## 3-8 Back Clip

The back clip is located on the back of the device. Operators can use it to secure the device to the patient's body for convenient use.

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
# 4 Device Operation

## 4-1 Preparation

The 6108 Temporary Pacemaker should only be operated by trained medical professionals in clinical environments such as hospital operating rooms or wards. Use the patient cable provided with the product, together with an adapter or extension cable (if required), to connect to permanent or temporary pacing leads. The 6108 temporary pacemaker, patient cable, adapter/extension cable, permanent/temporary pacing lead, and the patient's myocardium together form the pacing system. The following are the steps for connecting the pacing leads. Please ensure that all connections are secure before use.


1. Connect the patient cable to the device following the instructions provided in section 3.2 "Patient Cable Socket."

2. Connect the pacing lead with the adapter/extension cable to the patient cable via the electrode connector on the patient cable.

1) Turn the tightening screw on the electrode connector of the patient cable counterclockwise until the pin of the pacing lead connector can be inserted.

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2) Fully insert the pin of the pacing lead connector into the electrode connector of the patient cable until it cannot be pushed in any further.

3) Turn the tightening screw on the electrode connector of the patient cable clockwise to secure the pin of the pacing lead connector. At this point, gently pull the pacing lead to confirm that it cannot be pulled out.

4) Pay attention to polarity when connecting the patient cable. The electrode connector on the patient cable is marked with (+) and (–) polarity symbols. When using a bipolar pacing lead, connect the adapter/extension cable of the lead tip (distal electrode) to the (–) pole, and connect the adapter/extension cable of the lead ring (proximal electrode) to the (+) pole. When using a unipolar pacing lead, connect the cable from the cardiac electrode to the (–) pole, and connect the cable from the body surface electrode to the (+) pole.

3. Double check the connections among the patient cable, pacing lead, and the 6108 Temporary Pacemaker to ensure they are secure and safe.
When disconnect the patient cable, follow the instructions provided in section 3.2 "Patient Cable Socket" to disconnect the patient cable.

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## 4-2 Power On

Press the "⏻" button on the 6108 Temporary pacemaker to power it on.

Before powering on, ensure that all device settings meet the usage requirements.

Note: Do not operate the 6108 Temporary Pacemaker if the ambient temperature or humidity is outside the specified operating range.

If there is a sudden change in ambient temperature (e.g., bringing the device indoors from the cold in winter), allow the device to reach the ambient temperature before proceeding with use.

Note: Do not press any buttons or turn any knobs while the device is powering on, as this may cause the device to misinterpret the input as a malfunction and fail to start properly.

After the device is powered on, it operates according to the following parameters:

| | |
|---|---|
| Pacing mode： | When the sensitivity knob is set to "ASY," the device operates in SOO mode. In all other cases, it operates in SSI mode. |
| Rate： | As per knob setting |
| Amplitude： | As per knob setting |
| Pulse width： | 1.8ms |
| Sensitivity： | As per knob setting |

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Note: After adjusting the parameters, the protective cover should be closed. If the device detects an error during startup, the amber "⏻" indicator will flash and the buzzer will sound. Device startup will also be delayed.

In addition, the device performs continuous real-time self-checks during operation to ensure it is functioning properly through comprehensive diagnostics.

## 4-3 Power-Off

The 6108 Temporary Pacemaker uses a two-step shutdown mechanism to prevent accidental power-off.

To shut down the device, follow these steps:

1. Press "⏻" button once – pacing, sensing, and low battery lights flash sequentially.
2. Press the "⏻" button again before the rolling indicator lights stop flashing: the device will power off.

In any case, if the user wishes to cancel the shutdown process and keep the device running, do not press the "⏻" button again before the rolling indicator lights stop flashing.

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When the 6108 Temporary Pacemaker is operating in single-chamber synchronous pacing mode (SSI), the pacing timing is synchronized with the patient's intrinsic heart rhythm. When the patient's intrinsic heart rate is higher than the pulse rate set on the 6108 temporary pacemaker, the device's pacing is inhibited.

| First Letter | Second Letter | Third Letter |
|---|---|---|
| Pacing chamber | Sensing chamber | Response mode |
| O = no function<br>S = single chamber | O = no function<br>S = single chamber | O = Indicates pulse generator has no sensing function<br>T = Triggered<br>I = Inhibited |

To set the device to SSI mode:
Adjust the "Sensitivity" knob to any setting other than "ASY."

## 4-5 Asynchronous Pacing Mode (SOO)

In asynchronous mode, the 6108 Temporary Pacemaker delivers pacing pulses at the preset pulse rate without sensing the patient's electrocardiographic signals.

**To set the device to SOO mode:**

Turn the "Sensitivity" knob counterclockwise and set the sensitivity value to "ASY." .

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coincide with the patient's intrinsic T wave, potentially leading to life-threatening arrhythmias. Cautious clinical judgment should be exercised when using asynchronous pacing (e.g., in patients without intrinsic rhythm or those fully dependent on pacing), and appropriate pacing energy should be selected.

## 4-6 Pulse Amplitude

The user can adjust the pulse amplitude using the "Voltage" knob, with a range from 0.1 V to 10.0 V. To decrease the pulse amplitude, turn the knob counterclockwise; to increase it, turn the knob clockwise.

The pulse width of the pacing output is a fixed value of 1.8ms.

## 4-7 Pulse Frequency

The pulse rate refers to the pacing pulse delivery frequency of the 6108 Temporary Pacemaker.
In SSI pacing mode, the device begins timing from the most recent pacing or sensing event, searching for the patient's intrinsic heartbeat within the preset interval (t = 60/pulse rate in seconds).

If cardiac activity is detected within this interval, the device inhibits pacing and restarts the timing cycle from that point. If no intrinsic activity is detected, a pacing pulse is delivered at the end of

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the interval. In asynchronous pacing mode, the device does not sense the patient's intrinsic cardiac activity and delivers pacing pulses at the fixed preset rate.

**Adjust the pulse rate as follows:**
To decrease the pulse rate, turn the knob counterclockwise. To increase the base rate, turn the knob clockwise.

## 4-8 Sensing Sensitivity

Adjust sensing knob to set sensitivity:
Clockwise: increases sensitivity (lower threshold)
Counterclockwise: decreases sensitivity (higher threshold)

## 4-9 Run-Away Protection

The runaway protection function is designed to limit uncontrolled pacing output and ensure patient safety. The 6108 Temporary Pacemaker's runaway protection is composed of three independent mechanisms that monitor each other. Under this protection system, the device's output rate will not exceed 220ppm.

When the device first detects runaway output, it will automatically revert to the current set parameters. If the runaway condition persists and is detected again, the device will shut down after the buzzer sounds.

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## 4-10 Sensing Events

In synchronous pacing mode, the 6108 Temporary Pacemaker uses a standard threshold detection method to identify P/R wave peaks. When the pacing lead is properly connected and implanted in the patient's heart chamber, the device collects intracardiac ECG signals through the lead and compares them with the sensitivity threshold set on the device. When the collected EGM signal voltage exceeds the sensing sensitivity threshold, the device interprets it as a sensed P/R wave, that is, a sensing event is detected.

## 4-11 Refractory Period

The refractory period is a fixed time interval following each pacing or sensing event. For the 6108 Temporary Pacemaker, the refractory period is 250 milliseconds. Sensing events that occur during the refractory period will be ignored and will not alter the timing of output pulses.

## 4-12 Noise Reversion

The 6108 Temporary Pacemaker uses frequency detection methods to continuously monitor and analyze ECG signals in real-time for noise interference that may affect sensing accuracy. When excessive environmental noise is detected, the device automatically switches to asynchronous pacing mode and operates at the preset pacing rate. When the interference is reduced to a level that no longer affects

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sensing accuracy, interference reversion (noise reversion) will automatically stop.

## 4-13 Battery Status

In SSI pacing mode, with a base rate of 70ppm, a pulse amplitude of 5.0 V, and a sensing sensitivity of 2.0 mV, two new AA (LR6) alkaline batteries can power the 6108 temporary pacemaker for 18 days. The device continuously monitors battery status during operation and alerts the user when the battery is low.

The low battery warning is indicated by the flashing of the amber "⏻" light and an audible buzzer. If either warning occurs, the user must replace the batteries immediately.

**Low battery warnings include:**

1. "⏻" Light: Located at the lower right corner on the front of the device. This indicator flashes amber to remind the user to replace the batteries.
2. Buzzer: When the batteries are nearly depleted, the 6108 Temporary Pacemaker emits a continuous "beep-beep" sound, alerting the user to replace the batteries immediately to ensure proper device operation.

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## 4-14 Battery Replacement

Please use new batteries for pacing therapy with each new patient.

1. If the device is in continuous use, replace the batteries at least once every two weeks.
2. Replace the batteries promptly whenever any low battery indicator appears.
3. During battery replacement, the 6108 Temporary Pacemaker can operate without batteries for at least 55 seconds.
4. Remove the batteries when the device is not in use.

**Warning :**

1. Whenever possible, avoid replacing the batteries while the device is connected to the patient.
2. Do not use rechargeable batteries.
3. Using batteries that do not match the recommended model or size, or using batteries with contaminated electrode surfaces, may result in unstable output or no output from the device.

The battery compartment is located on the side of the 6108 Temporary Pacemaker and is designed to hold two AA (LR6) alkaline batteries. There are battery compartment buttons on both sides; pressing both buttons simultaneously will open the compartment.

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**Battery Replacement**

1. Open the battery compartment following the instructions in section "3.5 Battery Compartment" and remove the old batteries.
2. Insert new batteries according to the polarity markings inside the compartment.
3. Close the battery compartment as described in section "3.5 Battery Compartment."
4. If you notice any debris, liquid, or signs of rust or corrosion on the battery terminals during replacement, please contact Lifetech Cardio after-sales service for support.

## 4-15 Cleaning and Disinfection

Wipe the 6108 Temporary Pacemaker and the Patient Cable thoroughly using a damp cloth or sponge moistened with clean water.

Then, disinfect the device and patient cable by wiping them with a sponge or cloth moistened with 75% ethanol.

After each cleaning and disinfection, a qualified technician must verify that the device is functioning properly. Before each use, ensure that the device and patient cable are completely dry, especially the electrode connectors and the patient cable socket.

Note: Immersing the 6108 Temporary Pacemaker in water or disinfectant, or disinfecting the device using gamma

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radiation or high-temperature steam, may damage the device and render it inoperable.

## 4-16 Post-Cleaning Inspection

After each cleaning/disinfection, perform the following checks on the patient cable and the 6108 Temporary Pacemaker:

- Inspect the appearance of the device and patient cable to ensure there is no damage and that the surfaces are free of contamination.
- Ensure that both the device and patient cable are completely dry, especially the electrode connectors and the patient cable socket.
- Confirm that the display screen functions properly and that all buttons and knobs are operating normally.
- Ensure that the patient cable is correctly connected to the device and that the device completes its power-on self-test without any errors.

If one or more of the above checks cannot be completed successfully, please contact Lifetech Cardio after-sales service. Contact information can be found on the front page.

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# 5 Application

## 5-1 Preparation

1) Connect the electrocardiograph and/or other monitoring devices to the patient as needed based on clinical requirements.
2) Prepare the 6108 Temporary Pacemaker, power it on, and verify that the device is functioning properly and has sufficient battery power.
3) Connect the 6108 Temporary Pacemaker to the patient according to the instructions in section "4.1 Preparation."

## 5-2 Sensing Peak Measurement

1) In SSI pacing mode, set the pacing rate to at least 10 bpm lower than the patient's intrinsic rate, as described in 4.7 Pulse Rate.
2) Observe the "Sensing" indicator. If it does not flash continuously, lower the sensitivity value (increase sensitivity) as described in 4.8 until continuous flashing occurs. If it flashes continuously, raise the sensitivity value (decrease sensitivity) until it no longer flashes continuously.
3) This value represents the estimated sensing peak for the lead position.
4) Set the sensing value to half of the measured peak, as described in 4.8. For example, if the peak is 2 mV, set sensitivity to 1 mV.

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5) Confirm that the pacing rate meets clinical needs. Adjust as needed following 4.7 Pulse Rate.

## 5-3 Pacing Threshold Test

**Preparation**:

1) Set the pacing rate at least 10 ppm higher than the patient's intrinsic rate, per 4.7.

2) Following 4.6 Pulse Amplitude, start at 0.1V and gradually increase until consistent capture is observed. This value is the pacing threshold.
3) Set the output to at least twice the pacing threshold. For example, if threshold is 2V, set amplitude to ≥4V.
4) Ensure the pacing rate is appropriate; if not, adjust according to 4.7 Pulse Rate

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Service life: 8 years.

The 6108 Temporary Pacemaker has no replaceable internal parts. Any faults can be returned for factory service. Lifetech Cardio guarantees availability of repair parts for at least 10 years from delivery.

**Device Dimensions**

The physical dimensions of the 6108 Temporary Pacemaker are as follows (with the hanging ring retracted):

| | Dimensions ( mm ) | Tolerance ( mm ) |
|---|---|---|
| Length | 133 | ±10 |
| Width | 75 | ±8 |
| Thickness | 58 | ±5 |

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| | | |
|---|---|---|
| Maximum diameter of compatible adapter connectors | 3 | NA |

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**Battery**

The 6108 Temporary Pacemaker uses two AA (LR6) alkaline batteries.

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1) At an ambient temperature of 20°C, the 6108 Temporary Pacemaker is expected to operate continuously for at least 14 days using two new AA (LR6) alkaline batteries.

2) Under the same temperature condition (20 °C), using default pacing and sensing parameters (pacing rate: 70 ppm, pulse width: 1.8ms, pulse amplitude: 5.0 V, sensing sensitivity: 2.0 mV), the typical continuous operating time is up to 18 days with two new AA (LR6) alkaline batteries.

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| Battery Type | Battery Life(days) |
|---|---|
| AA (LR6) alkaline battery | Typically:18; Mnimum:14 |

**Backup Power Supply**

The 6108 Temporary Pacemaker can rely on backup power to ensure at least 55 seconds of continuous operation.

**Software Release Version**

The software release version of the 6108 Temporary Pacemaker is Version 1.

**Pacing Modes**

According to clinical requirements, the

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### Chamber

The 6108 Temporary Pacemaker can pace either the atrium or the ventricle independently.

### Runaway Protection Rate

The runaway protection mechanism of the 6108 Temporary Pacemaker is triggered at rates exceeding 220 ppm.

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| | | | |
|---|---|---|---|
| Knob value | Knob value | 1.8/ms | Knob value |

**Rate**

| Range (ppm) | Tolerance (%) | Step (ppm) |
|---|---|---|
| Min: 30<br>Max: 180 | ±5 | 10, from 30 to 50<br>5,from 50 to 80<br>10,from 80 to 100<br>20,from 100 to 180 |

**Amplitude**

| Range (V) | Tolerance (V/%) | Step (times/minute/ppm) |
|---|---|---|
| Min: 0.1<br>Max: 10 | ±0.1,from 0.1 to 0.5 V:<br>±10%,from 1.0 to 10.0 V: | 0.2,from 0.1 to 0.5<br>0.5,from 0.5–3.0:<br>1.0,from 3.0–10.0: |

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**Sensitivity**

| Range (mv) | Tolerance (%) | Step (mv) |
|---|---|---|
| Min: 0.5<br>Max: 20 | ±20, from 0.5 to 1.0mV | 0.5, from 0.5 to 3.0<br>2.0, from 3.0 to 6.0 |
| | ±10, from 1.5 to 20.0 mV: | 5.0, from 10.0 to 20.0 |

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Essential performance of the device: Based on GB 9706.1-2020 and YY

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9706.231-2023, the essential performance of this product is as follows:

1. Battery depletion indicator

2. ME equipment parameter stability

3. Pulse amplitude stability

4. Deliberate action required to change settings

5. Parameter stability at onset of the battery depletion indicator

6. Runaway protection

7. Interference reversal in the presence of sensed electromagnetic disturbance or electrical energy sources

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Power supply mode: Two AA 1.5V alkaline batteries (LR6)

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Signal output or input section: EGM input, pacing signal output, and data

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Altitude: ⩽4000m

Other: Handheld device

## 6-3 Electromagnetic Compatibility (EMC) Statement

| **Guidance and Manufacturer's Declaration - Electromagnetic Emissions** | | |
|---|---|---|
| The 6108 Temporary Pacemaker is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment. | | |
| **Emission test** | **Compliance** | **Electromagnetic environment - guide** |
| RF Emissions<br>GB 4824 | Group 1 | The Temporary Pacemaker uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF Emissions<br>GB 4824 | Class A | The 6108 Temporary Pacemaker is suitable for use in all establishments other than domestic households and those directly connected to the public low-voltage power supply network. |
| Harmonic emissions<br>GB 17625.1 | Not applicable | |
| Voltage fluctuations/flicker emissions<br>GB 17625.2 | Not applicable | |

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| Guidance and Manufacturer's Declaration - Electromagnetic Immunity | | | |
|---|---|---|---|
| The 6108 Temporary Pacemaker is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment. | | | |
| **Immunity test** | **IEC 60601 test level** | **Compliance level** | **Assessment** |
| Electrostatic discharge GB/T 17626.2 | ±6 kV contact discharge ± 8 kV air discharge | ±8 kV contact discharge<br>± 15 kV air discharge | Floors should be wood, concrete or ceramic tile. If covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical Fast Transient/Burst GB/T 17626.4 | Not Applicable | Not Applicable | / |
| Surge GB/T 17626.5 | Not Applicable | Not Applicable | / |
| Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines GB/T 17626.11 | Not Applicable | Not Applicable | / |
| Power Frequency Magnetic Field (50 Hz/60 Hz) GB/T 17626.8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment |
| Conducted RF GB/T | 3 V (rms) | 3 V (rms) | Portable and |

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| | | | |
|---|---|---|---|
| 17626.6<br>Radiated RF GB/T 17626.3 | 150 kHz to 80 MHz (excluding ISM bands)<br><br>10 V (rms)<br>150 kHz to 80 MHz (for ISM bands)<br><br>10 V/m<br>80 MHz to 2.5 GHz | 10 V (rms)<br><br>10 V/m | mobile RF communications equipment should be used no closer to any part of the 6108 Temporary Pacemaker, including cables, than the recommended separation distance calculated from the applicable formula. Recommended separation distance:d = [1.17]√Pd = [1.20]√Pd = [1.20]√P (80 MHz – 800 MHz)d = [2.30]√P (800 MHz – 2.5 GHz)Where:P = max rated output power of the transmitter in watts (W)d = recommended separation distance in meters (m)<br>The field strengths of fixed RF transmitters should be determined |

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| | | | |
|---|---|---|---|
| | | | by an electro-magnetic site survey and should be less than the appli-cable compli-ance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following sym-bol.  |

Note 1: At 80 MHz and 800 MHz, the separation distance formula for the higher frequency range applies.

Note 2: These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) ISM bands between 150 kHz and 80 MHz are: 6.765‑6.795 MHz, 13.553‑13.567 MHz, 26.957‑27.283 MHz, and 40.66‑40.70 MHz.

b) An additional factor of 10/3 is used to reduce the likelihood of interference from mobile/portable communications equipment inadvertently brought into the patient area.

c) The field strength from fixed transmitters (e.g., base stations for radio, TV, cordless phones, etc.) is not predictable. A site survey is recommended.

d) In the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and the 6108 Temporary Pacemaker**

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The 6108 Temporary Pacemaker is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The purchaser or user of the 6108 Temporary Pacemaker can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the 6108 Temporary Pacemaker, according to the maximum output power of the communication equipment, as recommended below.

| Maximum Rated Output Power of Transmitter (W) | Corresponding Separation Distances by Frequency (m) | | | |
|---|---|---|---|---|
| | 150 kHz～80 MHz (excluding ISM bands) d=[1.17]√P | 150 kHz～80 MHz (excluding ISM bands) d=[1.20] √P | 80 MHz～800 MHz d=[1.20] √P | 800 MHz～2.5 GHz d=[2.30] √P |
| 0.01 | 0.12 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.37 | 0.38 | 0.38 | 0.73 |
| 1 | 1.17 | 1.20 | 1.20 | 2.30 |
| 10 | 3.70 | 3.79 | 3.79 | 7.27 |
| 100 | 11.70 | 12.00 | 12.00 | 23.00 |

For transmitters rated at maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter. P is the maximum rated output power of the transmitter in watts (W), as specified by the manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: The ISM frequency bands between 150 kHz and 80 MHz are 6.765-6.795 MHz, 13.553-13.567 MHz, 26.957-27.283 MHz, and 40.66-40.70 MHz.

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# 7 Support, Disclaimer and Warranty Information

## 7-1 Support

Lifetech Cardio has a professional technical team to provide operators with technical support and after-sales service. Contact information is provided on the front page. Each 6108 Temporary Pacemaker has a serial number on the back; please provide this information when requesting after-sales service or technical support.

**Fault Handling**

In case of equipment failure or accidents during proper use (e.g., fluid ingress, dropping), which do not harm the operator, please replace the device with another 6108 Temporary Pacemaker to continue pacing the patient. Turn off the faulty device, attempt to restore its original condition, and restart it after waiting ten seconds. If the issue persists, please contact Lifetech Cardio after-sales service for technical support. Contact information is provided on the front page.

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**Professional Training**

The Lifetech Cardio after-sales team is dedicated to providing training on the use of the 6108 Temporary Pacemaker for operators in need.

## 7-2 Disclaimers

**6108 Temporary Pacemaker**

The 6108 Temporary Pacemaker provides temporary pacing for arrhythmias and bradycardia. However, Lifetech Cardio does not guarantee that the device will be effective for all patients.

**Patient Cable**

All cables are subject to fatigue and damage due to repeated use. Patient cables may also suffer from improper handling or repetitive stress under various use environments. Therefore, Lifetech Cardio does not guarantee that the patient cables are immune to damage from long-term repetitive use. Although the patient cables are carefully designed, manufactured, and tested, they may still be damaged due to unforeseen circumstances.

Warnings on product labels provide additional details and are part of this disclaimer. Lifetech Cardio disclaims

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all express and implied warranties related to this product. Lifetech Cardio shall not be liable for any incidental or consequential events arising from the use, defect, or failure of the patient cable, whether the claims are based on warranty, contract, or tort.

The above disclaimers and limitations are not intended to conflict with mandatory provisions of applicable laws. If any part or clause is deemed unlawful, unenforceable, or inconsistent with applicable laws by a court of competent jurisdiction, the remainder of the warranty shall remain in effect, and all rights and obligations shall be honored as if the invalid portion had not been included.

## 7-3 Warranty Information

**Overview**

The 6108 Temporary Pacemaker includes a one-year warranty from the delivery date, provided no man-made damage has occurred.

If non-artificial damage occurs during normal use, please contact Lifetech Cardio after-sales service. Contact details are on the front page.

Lifetech Cardio reserves the right to choose between the following repair options:

1) Repair or replacement of any defective parts of the device;

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equivalent functionality.

**Limitations**

This warranty applies only under the following conditions:

1. The device has not been tampered with, disassembled, or repaired by anyone other than Lifetech Cardio.
2. The malfunction was not caused by improper handling, misuse, or accident.
3. Lifetech Cardio was contacted within 60 days of device failure.

**Disclaimer**

This warranty is limited to the terms explicitly stated. Lifetech Cardio shall not be liable for any incidental or consequential events related to the use, defect, or failure of the device, regardless of whether claims are based on warranty, contract, or tort.

The above disclaimers and limitations are not intended to conflict with any mandatory provisions of applicable law. If any part or provision of this warranty statement is found to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remainder of the warranty statement shall not be affected. All rights and obligations shall remain acknowledged and in effect as if the invalid portion had not been included.

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# Lifetech Cardio Temporary Pacemaker

## Manufacturer Information

**Registrant / Manufacturer**
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

**Registrant/Production Location**
1401, Building D1, Nanshan i Park
No.1001, Xueyuan Avenue, Taoyuan Street
Nanshan District
Shenzhen, Guangdong, China

Tel: +86 755 23251894
WWW.LTCARDIO.COM

## Device information and Support

Product Name: Temporary Pacemaker
Mode: 8400, 8301, 8201, 8200, 7101, 7100
After-sales Service Provider:
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.
After-Sales Service Hotline: 400-803-5503

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

4-14 Pacing Pause ..... 57
4-15 Emergency Pacing ..... 58
4-16 EGM Waveform ..... 59
4-17 Event Markers ..... 61
4-18 Analysis ..... 61
4-19 Pacing Impedance ..... 61
4-20 Peak Amplitude ..... 62
4-21 EGM Waveform Comparison View ..... 63
4-22 Pacing Percentage ..... 64
4-23 Rapid Atrial Pacing ..... 64
4-24 Rapid Ventricle Pacing ..... 66
4-25 Bluetooth Setting ..... 67
4-26 Battery Status ..... 67
4-27 Battery Replacement ..... 69
4-28 Cleaning and Disinfecting ..... 71
4-29 Inspection ..... 72
5 Analysis Application ..... 74
5-1 Preparation ..... 74
5-2 Sensitivity Test ..... 74
5-3 Pacing Capture Threshold ..... 75
5-4 Impedance Test ..... 76
6 Pacing Parameters & EGM Data Storage and Transfer Software ..... 77
7 Specification ..... 81
7-1 Parameters ..... 81
7-2 Standards Compliance ..... 89
7-3 Electromagnetic Compatibility Statement ..... 91
7-4 Model and Specification ..... 97
8 Support, Disclaimer and Warranty ..... 98
8-1 Support ..... 98
8-2 Disclaimers ..... 99
8-3 Warranty ..... 101

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Humidity limit (0%-90%RH)

General warning sign

Refer to instruction manual/ booklet.

MR Unsafe

Act according to your local rules and do not dispose of the product with household waste.

Defibrillation-proof type CF applied part

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

# 1 Device Overview

## 1-1 Description

The Lifetech Cardio Temporary Pacemaker consists of a Non-implantable Pulse Generator and a Patient Cable. Optional Analysis Cable (models 7101 and 7100 excluded)

The Lifetech Temporary Pacemaker is a handheld device powered by two AA (LR6) alkaline batteries. The device is directly connected to one end of the patient cable/the analysis cable, while the other end of the patient cable/the analysis cable is connected either directly or via an adapter cable to the pacing lead, which in turn is connected to the myocardial tissue. The pacing system consists of four components: the temporary pacemaker, the patient cable/the analysis cable, the pacing lead, and the myocardial tissue.

The Lifetech Cardio Temporary Pacemaker provides single chamber synchronous/asynchronous pacing therapy and can also be used to measure parameters such as pacing impedance and P/R wave amplitude of the pacing lead system to help determine the implantation status of the cardiac pacing electrode lead.

The device supports either SSI (synchronous) or SOO (asynchronous)

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pacing modes. The Rate, Amplitude, and Sensitivity are adjustable. It can also measure P/R Wave Amplitude and Impedance, and can record and display the Pacing Percentage events over up to 5 days. In addition, it offers online EGM morphology waveforms, pacing/sensing event markers, and marker intervals. The Lifetech Cardio temporary pacemaker can provide high-frequency atrial/ventricular pacing.

The Lifetech Cardio Temporary Pacemaker is intended to be operated by trained professionals in a clinical setting.

## 1-2 Package Content

The product package contains:

1) Temporary Pacemaker *1
2) Patient Cable*1
3) AA 1.5V Alkaline (LR6) batteries*2
4) Carrying Case*1
5) Instructions for Use *1
6) Certificate of Conformity *1
7) Analysis Cable*1(Optional accessory)
8) Pacing Parameters & EGM Data Storage and Transfer Software (Optional)

The manufacturing date of the product can be found on the device label.

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## 1-3 Features

The functional features of the Lifetech Cardio Temporary Pacemaker are summarized as follows.

**Pacing Features**

| | |
|---|---|
| **Pacing Modes** | SSI, SOO. |
| **Pulse Width** | Adjustable, range: 0.06 – 2.0ms |
| **Amplitude** | Adjustable, range: 0.1 – 10.0V |
| **Rate** | Adjustable, range: 30 – 200 ppm |
| **Sensitivity** | Adjustable, range: 0.5 – 20.0mV |

**Functional Features**

| | |
|---|---|
| **Online EGM and Marker** | Longer than 3-second rolling EGM waveform in configuration view with pacing or sensing even markers as well as event interval.<br>Longer than 1.5-second rolling EGM waveform in the “ANALYSIS” view with pacing or sensing event markers as well as event interval. |

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| | |
|---|---|
| **EGM Waveform Comparison** | In the EGM Waveform Comparison View, two real-time updated EGM waveforms (> 1.5 s) are displayed, high-pass filtered at 0.5 Hz and 30 Hz, respectively. |
| **Impedance** | True pacing impedance measurement, in the “ANALYSIS” View. |
| **Peak Amplitude** | Beat to beat measurement of sensing P or R wave amplitude in the “ANALYSIS” View. |
| **Pacing Percentage** | Pacing percentage statistics rolling trend over 5 days in the PACING% View. |
| **Rapid Atrial Pacing ( RAP )** | Deliver Rapid Atrial Pacing according parameter settings. |

**Safety Features**

| | |
|---|---|
| **Noise Reversion** | When electrical interference is sensed, the device reverts to asynchronous pacing at the preset rate until interference stops. |

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**Rate Limiting (Run-away Protection)**

Ensures that the pacing output rate does not exceed 230 ppm.

**Emergency Mode**

Pressing the "⚠" button once activates emergency pacing.

Emergency Settings:

| | |
|---|---|
| MODE | SSI |
| Rate | 70ppm |
| Amplitude | 10.0V |
| Pulse Width | 2.0ms |
| Sensitivity | 2.0mV |

**Pacing Pause**

Holding the "PAUSE" button suspends pacing for up to 10 seconds. Pacing automatically resumes when released.

**Defibrillation Protection**

The device has defibrillation protection to minimize malfunction risk during external defibrillation.

**Low Battery indicator**

When the battery power is low, the yellow power indicator light will flash, the buzzer will sound an alert, and the screen will display the message "Low battery, please replace the battery." The

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| | |
|---|---|
| | battery level indicator will also show a low battery status. Once the device begins alerting due to low battery, it can typically continue operating for at least 24 more hours. |
| **Battery Compartment Not Closed Alert** | When the battery compartment is opened, the screen of the Temporary Pacemaker will display the message: "Please replace battery and lock its cover ASAP." |
| **Backup Power** | ≥ 80 seconds continuous operation during battery replacement, accompanied by the flashing orange power indicator and a screen prompt. |
| **Automatic Lock Feature** | In the "ANALYSIS" View, the device will automatically lock the screen after 5 minutes of inactivity. In other interfaces, the screen will automatically lock after 1 minute of inactivity. When the device is in the auto-lock state, pressing any button other than the "⚠" button or the "🔒" button will prompt a message on the screen indicating that the device is locked. |

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| | |
|---|---|
| **Shutdown Protection** | When the main unit is operating normally (not in a locked state), pressing the "⏻" button will cause the display to show a shutdown confirmation message. If the operator presses the "⏻" button again before the message disappears, the device will shut down; otherwise, the shutdown operation will be canceled. |

## 1-4 Intended Use

The Lifetech Cardio Temporary Pacemaker is intended for external temporary pacing of the atrium or ventricle. It is also intended for analysis of the pacing lead system during implantation of pacemakers and defibrillators (models 7101 and 7100 excluded).

## 1-5 Indications

Examples of specific indications for temporary pacing may include:

1) Sick sinus syndrome
2) Sinus bradycardia
3) Atrial and/or ventricular arrhythmias
4) Complete atrioventricular block
5) Asystole
6) Bradycardia accompanied

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by congestive heart failure

7) Patient support, management and evaluation before the implantation of implantable pacemaker
8) Support during the replacement of implantable pacemaker
9) Cardiac complications occurring during intervention or surgery
10) Support after the cardiac surgery
11) Acute myocardial infarction complicated with cardiac conduction block
12) High-rate burst pacing for the treatment of supraventricular tachyarrhythmias

## 1-6 Contraindications

There are no known contraindications to for the use of temporary cardiac pacing and/or pacing lead system analysis.

## 1-7 Patient Cable/Analysis Cable

The Patient Cable/Analysis Cable provided in the accessory kit can be used for connection with the Temporary Pacemaker for temporary pacing/ pacing lead system analysis.

1) Patient Cable

The Patient Cable is approximately

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1860 mm in length, with one end featuring a plug and the other end equipped with an electrode connector. The electrode connector can be used to connect electrode leads or adapter cables with a connector diameter not exceeding 3.0 mm, and the plug is designed to connect to the Patient Cable Socket located on the top of Temporary Pacemaker.

2) Analysis Cable
The Analysis Cable is approximately 1860 mm in length with two Alligator Clips and a Connector Plug. The two Alligator Clips can accept the lead connector pins with the diameter more than 1.0mm and less than 5.0 mm. The Connector Plug is designed to connect the Patient Cable Socket at the top of Temporary Pacemaker.

Please contact Lifetech Cardio Support using information provided on the title page for product inquires.

| Accessory | Order Code |
|---|---|
| Patient Cable | 0202002007 |
| Analysis Cable | 0203001024 |

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## 2-1 Warnings

### Daily Use and Maintenance

**Device Modification.** Do not modify or tamper with the Temporary Pacemaker. Any tampering could render the device unfunctional and unsafe to use.

**Maintenance:** To maintain its specified performance, the device requires regular preventative maintenance. Note: Temporary pacemakers and their accessories should not be maintained or serviced while in patient use.

**Liquid Immersion.** Cleaning can be done by wiping the device with a damp cloth or sponge. Subjecting the device to liquid immersion may damage the device and render the device unfunctional and unsafe to use.

**Flammable and Explosives.** The device is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide. Do not use the device in the presence of flammable anesthetics. Explosion is a hazard under such conditions.

**Patient Monitoring.** The pacemaker is not intended for standalone patient monitoring. Use an ECG monitor or 12-lead ECG machine for appropriate

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patient monitoring.

**Diagnostic Use.** The Temporary Pacemaker is not intended for diagnostic purposes.

**MR Unsafe**. Magnetic resonance imaging (MRI) is a medical imaging that uses magnetic fields to generate internal views of the body. Do not expose the device to a magnetic resonance (MR) environment or perform MRI scans on patients who are using temporary pacemaker. The device may not function properly due to the strong magnetic and radiofrequency fields generated by the MR scanner. MRI scan may cause serious injury, tachycardia or system failure or damage.

**Usage:** Do not connect, place in close proximity to or stacked the Temporary Pacemaker with other devices. If close or stacked use is necessary, the system should be observed and verified to ensure proper operation under the specific configuration.
Except for transducers and cables sold by the device or system manufacturer as replacement internal components, the use of unauthorized accessories, transducers, or cables may result in increased electromagnetic emissions or reduced immunity of the device or system.

**Disposal:** Please dispose of the Temporary Pacemaker in accordance with

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local regulations.

Therapy Specific

**Threshold Test.** Loss of capture during threshold test may cause asystole.

**Asynchronous Pacing:** In asynchronous pacing mode, ventricular pacing pulses may coincide with the patient's intrinsic T wave, potentially triggering life-threatening arrhythmias. Cautious clinical judgment should be exercised when using asynchronous pacing (e.g., in patients with no intrinsic rhythm or those fully dependent on pacing), and appropriate pacing energy should be applied.

**Atrial Pacing:** In AV conduction disorders or atrial fibrillation/flutter, single-chamber atrial pacing is ineffective and should not be used.

**Pacing Termination.** Abrupt pacing termination may cause asystole. Before stopping, ensure synchronous mode and gradually reduce rate below intrinsic rhythm until recovery.

During Procedure

**Introduced Current.** Implanted leads constitute direct electrical connections to a patient's heart, induced or injected current by any means through the implanted leads could cause fatal arrhythmias. In addition, induced or injected current could cause over sensing in the pacemaker and result in loss

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of pacing therapy. Furthermore, injected or induced current of significant magnitude may damage the Temporary Pacemaker and lead system.

**Patient Cable/ Analysis Cable.** Make sure the Patient Cable/ Analysis Cable is connected to the device first prior to connecting a lead to the Patient Cable/ Analysis Cable. Use of accessories and cables other than those specified or provided by Lifetech Cardio of the Temporary Pacemaker could result in increased electromagnetic emissions or decreased electromagnetic immunity of the device and result in improper operation.

**Electrostatic Discharge.** Operator shall discharge any static electricity by touching a large conductive ground surface before contacting the patient, lead and the Temporary Pacemaker.

**Handling Implanted Leads.** Operate with caution and keep the metal connectors of an implanted lead from contacting the operator or any conductive or wet surfaces.

**Line-powered Equipment.** When using 220 V AC-powered equipment, ensure proper grounding.

**Electrocautery.** Keep the current path and return plate as far away as possible from the lead/pacemaker system.

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A minimum separation of 6 inches or 15 cm is generally advised.

**Ablation.** Avoid direct electrical contact between the Temporary Pacemaker/lead system and the ablation electrodes/contacts. Keep the current path and return plate as far away as possible from the lead/pacemaker system. A minimum separation of 6 inches or 15 cm is generally advised.

**Defibrillation.** Keep the defibrillation paddles as far away as possible from the lead/pacemaker system. A minimum separation of 6 inches or 15 cm is generally advised. Monitor the patient after each defibrillation therapy and verify the Temporary Pacemaker is functioning properly. Furthermore, an implanted lead system when connected to the pacemaker may shunt defibrillation energy away from the heart tissue that is close to the lead tip, this may reduce defibrillation efficacy and in addition may damage tissue surrounding the conducting electrodes.
**NOTE:** Prior to defibrillation when possible, disconnect the lead from pacemaker to eliminate the shunt pathway within patient's body through the device and lead system.

**RF Communication Equipment.**
Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than

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30 cm (12 inches) to any part of the TPM, including cables specified by Lifetech Cardio. Otherwise, degradation of the performance of this equipment could result.

**RF Immunity.** This device has been tested for radiated RF immunity only at selected frequencies, and use nearby of emitters at other frequencies could result in improper operation

**Electromagnetic Interference (EMI).** Strong electromagnetic interference may cause abnormal sensing and alter desired pacing function. It is recommended to use the asynchronous pacing mode at an appropriate base rate. Note that, the real-time waveform feature of the device can help to reveal an operating environment with strong EMI.

**Rapid Atrial Pacing (RAP).** RAP is for atrial use only. Before enabling RAP, be sure that the atrial leads are connected to the atrium, not the ventricle. High-rate burst pacing is contraindicated in the ventricle; it may result in life-threatening arrhythmias. RAP may result in tachycardia, acceleration of existing tachycardia, or fibrillation. Apply high rates under careful patient monitoring and control. Monitor the patient's ECG and blood pressure, and ensure that defibrillation equipment is immediately available.

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**Patient monitoring:** Continuously monitor the patient while the Temporary Pacemaker is in use to ensure proper functioning of the Temporary Pacemaker and to provide the patient with appropriate treatment.

## 2-2 Precautions

**Lead Polarity.** A bipolar lead system has greater immunity to interference than a unipolar lead system. The Temporary Pacemaker is optimized for use with a bipolar lead system, use only bipolar lead system with the device whenever possible.

**Emergency Pacing.** The emergency pacing mode is a preset of pacing and sensing configurations. In any functional state, even when device is turned off, pressing the Emergency button immediately configures the device for synchronous pacing (SSI) and delivers 2.0 ms pulse wide 10.0 V pacing stimulus at 70ppm. Care should be taken on the invocation of the emergency pacing mode, as this mode overwrites all programmed pacing and sensing settings to emergency settings. To exit emergency pacing and return to other pacing settings, manually adjust the relevant parameters as needed.

**Pacing Pause.** Press and hold the pause button that mutes the output of

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the pacing stimulus, and maintains programmed pacing and sensing configurations. The longest pause duration is 10 seconds. This mode shall be used with care, as pacing is not delivered to the patient during the pause.

**Concomitant Pacing.** Concomitant pacing happens when the Temporary Pacemaker is used actively and simultaneously with other independent pacing instrumentations on the same patient. Concomitant pacing could result in malfunctions of the pacing and sensing functions of the device, and the use should be strictly avoided.

**Battery Status.** The device supports more than 80 seconds of extended operation during battery replacement, care shall be taken to always check the battery status using the battery status indicators on the device (Refer to Battery Status) daily when device is in continuous operation.

**Battery Replacement.** Use two quality size AA 1.5V alkaline (LR6) batteries. Replace battery immediately when the battery indicator icon on the LCD shows empty and the Power LED starts blinking. Batteries replacement is recommended for every new patient and/or for every two weeks.

**Environmental Precautions.** Overstress beyond the normal environmen-

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tal condition could damage the Temporary Pacemaker system. To ensure the device functioning properly within its manufactured specifications, care should be taken according to the following environmental precautions:

- Device should be operated by trained professionals in a clinical environment.
- Safe electrostatic discharge practices should be followed.
- Store the device with batteries removed from Battery Compartment.
- Do not subject the device to mechanical damage, for example, as a result of mechanical shock or inappropriate handling.
- Do not store or operate the device outside the specified temperature and humidity range (see "Device Specifications" for details).
- Dispose of used batteries in accordance with local environmental regulations.
- Do not place the device where the patient may interact with it.
- Do not subject the device to liquid immersion; avoid spills on the de-

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vice, the connectors and the Patient Cable.

- Do not sterilize the device with gamma radiation or steam.

- If the ambient temperature or humidity changes suddenly, allow the pacemaker to reach ambient conditions before use.

For electromagnetic compatibility information, refer to Section 7.3 "Electromagnetic Compatibility Declaration." The device must be operated in the electromagnetic environment specified in Section 7.3.
Potential adverse environmental factors include, but are not limited to, the conditions listed above. In principle, the device should be operated with priority given to patient health and safety, minimizing the impact of environmental conditions on device performance.

## 2-3 Adverse Effects

The use of temporary pacing may lead to potential adverse effects. Some examples are listed below.

1) Abrupt termination of pacing could cause asystole before the intrinsic rhythm is re-established.

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
with T waves may induce tachyarrhythmias.

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
# 3 User Interface

The user interface of the Temporary Pacemaker consists of the following parts.

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(1) LCD screen
(2) “PACE”/“SENSE”/“POWER” LED indicators
(3) “PAUSE” button
(4) “ ⚠ ” button (Emergency button)
(5) “FUNC” button
(6) “PULSE WIDTH” button
(7) “AMPL” button
(8) Dial
(9) “ZOOM” button
(10) “⏻” button (Power button)
(11) “RATE” button
(12) “SENS” button
(13) “🔒” button (Lock button)
(14) Patient Cable Socket
(15) Hang Ring
(16) Rotary Latch
(17) Battery Compartment
(18) Battery Compartment Key

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# 3-1 Dial and Buttons

There are 1 dial for configuration adjustment and 10 buttons each for a dedicated function.

Dial

User can use the dial in each of the four configuration views, i.e. Rate configuration view, Amplitude configuration view, Pulse Width configuration view, Sensitivity configuration view, ANALYSIS view, RAP Configuration View, RAP Configuration View and Setting View to adjust the value of a particular setting.
Note that, the dial is only active in the 8 views.

Buttons

There is a total of 10 buttons grouped in the middle of the front side of the device.
Each of the 10 buttons is dedicated to one particular function:

"🔓" button (Lock button)

The "🔓" button (Lock button) is located in the center of the dial. User can toggle the lock state of all but-

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tons, except the Lock, and Emergency buttons which are always unlocked, by pressing and releasing the Lock button. Note that, the device automatically locks down buttons after 5 minutes of inactivity in the ANALYSIS View, and locks down buttons after 1 minute of inactivity in other views.

“⏻” button (Power button)

The “⏻” button (Power button) is the first one from the right amongst the top row of buttons. This button turns the device on or off.

The Power button is supervised by the button lock.

User can power on the device via a single press of the Power button.

User can power off the device through a two-step process, this process is implemented to protect the user from accidentally powering off the device:

First single button press invokes the shutdown confirmation prompt on the LCD screen and initiate the two-step shutdown process, and


Second single button press within the time when the shutdown confirmation prompt exists

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
to confirm the prompt, then the device will be turned off.
If you press the power button again after the power off prompt disappears, the device will not turn off, but the device screen will appear the power off prompt again.

“ZOOM” button

The “ZOOM” button is the second one from the right amongst the top row of buttons.
User can use this button to adjusts the vertical scale of the dynamic EGM waveforms in the Rate Configuration View, Amplitude Configuration View, Pulse Width Configuration view, Sensitivity Configuration View, ANALYSIS View, EGM Waveform Comparison View. Each press of the zoom button rotates through the zoom settings in the following fashion, 1x -> 2x -> 4x -> 8x -> 16x -> 1x -> 2x.

“PAUSE” button


The “PAUSE” button is the third one from the right amongst the top row of buttons. The Pause button is supervised by the button lock. User can mute the pacing output by pressing and holding down the pause button. The longest pause duration invoked by a single button press is limited to 10 seconds by design. Pausing the

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pacing output does not alter the ongoing timing of the Pacing and Sensing sequence.

“ ⚠ ” button (Emergency button)

The “ ⚠ ” button (Emergency button) is the first one from the left amongst the top row of buttons. This button is protected by a mechanical guard ring, from unintended activation.

Pressing the Emergency button immediately powers on the device if the device is off, and configures and operates the device, unaffected by the button lock, at the emergency presets:

Pacing Mode: SSI
Rate: 70 ppm
Amplitude: 10.0V
Pulse Width: 2.0ms
Sensitivity: 2.0mV

Note that, the emergency mode over-writes existing pacing and sensing configurations with the emergency preset. Adjusting any pulse parameter or sensing parameter will exit the emergency mode. When exiting the emergency pacing state, the parameter settings prior to the emergency pacing state cannot be automatically restored. Please manually adjust each parameter according to clinical needs.

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Configuration Buttons

The 4 configuration buttons are placed in an arc immediately surrounding the dial. These 4 buttons are supervised by the button lock.

From left to right, the buttons are:

“PULSE WIDTH” button: By pressing the “PULSE WIDTH” button, user can toggle between the home view and the Pulse Width configuration view. The value of Pulse Width can be adjusted by using the dial.

“AMPL” button: By pressing the “AMPL” button, user can toggle between the home view and the Amplitude configuration view. The value of Amplitude can be adjusted by using the dial.

“RATE” button: By pressing the “RATE” button, user can toggle between the home view and the Rate configuration view. The value of Rate can be adjusted by using the dial.

“SENS” button: By pressing the “SENS” button, user can toggle between the home view and the Sensitivity configuration view. The value of Sensitivity can be adjusted by using the dial.

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

The “FUNC” button is the triangle shaped button at the left within the button panel. The “FUNC” button is supervised by the button lock. User can use this button to rotate through the Home View -> ANALYSIS View -> EGM Waveform Comparison View-> PACING% View - >RAP View- >RVP View -> Setting View.


## 3-2 Patient Cable Socket

The Patient Cable Socket is located at the top of the Temporary Pacemaker, through which the patient cable/ analysis cable connects to the device.

To connect a Patient Cable/Analysis Cable to the device

1. Align the Connector Plug to the Patient Cable Socket, rotate the Connector Plug until the guides in the plug and socket are aligned and the plug can be further inserted into the socket.

2. Push the Connector Plug into the Patient Cable Socket fully until the connector clicks and latches into position.

3. Gently pull back on the Connector Plug to verify that the Connector Plug is securely latched.

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To disconnect a Patient Cable/Analysis Cable from the device

1. Hold the Temporary Pacemaker firmly in place. Using two fingers, gently squeeze the locking ring on the patient cable plug and pull the cable gently in the opposite direction of the device

2. Continue to gently pull the patient cable until it is fully disconnected from the device's connector.

Note: Do not pull out the Patient Cable/ Analysis Cable by holding the body part of the Patient Cable/Analysis Cable. This makes it impossible to pull out the Patient Cable/Analysis Cable and may cause cable damage.

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

The LCD screen is in the upper portion at the front side of the device. It presents the operator with pacing and sensing settings and analysis information of a connected pacing lead system. In addition, it informs the operator with symbolic and verbal prompts should user attention needed during device operation.


Startup View

When first powered on, the device shows a 3-second welcome screen. During this time, the device performs a set of build-in self-tests to ensure the device is functioning properly prior to use. When the device passes all power-on self-tests, the LCD then transits into the home screen. If the output voltage is low during the self-test, the LCD screen will display a error message "Code 0104: Low output amplitude". When the output voltage returns to normal, the message will disappear automatically and the device will resume normal operation. If other self-test errors occur during the self-test, the device will display the corresponding error code and automatically shut down after ten seconds.

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
## Home View


The home view is the default view. User can find a concise summary of the current pacing and sensing configurations. This view is laid out into 4 equal quadrants. The upper left, upper right, lower left and lower right quadrants are designated for current Rate, Amplitude, Pulse Width, and Sensitivity settings respectively.

## Information Bar

An information bar is located at the bottom of the screen to display information including, button lock status, battery level, pacing mode , emergency status and useful operation hints. This information bar is consistent across all views during device operation.


## Button Lock Indicator

When the button lock is active, the device displays a lock icon in the middle of the information bar on the LCD screen.

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
## Pacing Mode Indicator

Pacing mode is shown at the left in the information bar on the LCD screen.

## Battery status Indicator

Battery status is shown at the right inside the information bar on the LCD screen.


## Prompts

In the center of the screen, user may be prompted by the LCD via symbolic and verbal messages from time to time during normal device operation. The following prompts are designed in the Tempoaray Pacemaker to facilatate device operation.

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
7. Battery Cover Warning

Configuration Views

There are 4 configuration views designed in a unified layout. The 4 configuration views are designated to Rate, Amplitude, Pulse Width and Sensitivity respectively.

In a configuration view, user can find,

in the top left corner, the title of the configuration;

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
in the top center, an indicator bar for current settings and setting range;

in the top right corner, a numerical value for the setting;


at the bottom of the configuration, the information bar;

at the center of the view, a dynamic scolling waveform showing the online EGM morphology, markers of the device's Pacing or Sensing event channel and interval between either Pacing or Sensing events.

Rate Configuration View

Pressing the "RATE" button enters Rate Configuration view. User can configurate Rate in this view by rotating the dial.


Amplitude Configuration View

Pressing the "AMPL" button enters Amplitude Configuration view. User can configurate Amplitude in this view by rotating the dial.

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

## Pulse Width Configuration View

Pressing the “PULSE WIDTH” button enters Pacing Pulse Configuration view. User can configurate Pacing Pulse in this view by rotating the dial.


## Sensitivity Configuration View

Pressing the “SENS” button enters Sensitivity Configuration view. User can configurate Sensitivity in this view by rotating the dial.

## ANALYSIS View

User can obtain the impedance of the connected lead system, sensed P or R wave amplitude and Rate at the top of the ANALYSIS View. In addition, user can find in the view a longer than 1.5-second EGM dynamic scrolling waveform showing the real time EGM, pacing or sensing event as well as event interval.


NOTE: When the patient’s intrinsic rate is higher than the setting Rate, the device stops pacing and the impedance could not be measured. The ANALYSIS View indicates PACING REQUIRED to notify users.

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
To continue, adjust the Rate to deliver pulse and obtain the impedance of the connected lead system.


## EGM Waveform Comparison View

In this View, the user can observe more than 1.5 s of continuously updated EGM comparison waveforms. F1 represents the 0.5 Hz high-pass filtered waveform, and F2 represents the 30 Hz high-pass filtered waveform. Pacing and sensing event markers, as well as event intervals, are displayed.

## PACING% View

This view shows the user a 2D chart of Pacing Percentage during the past 5 days since last device power-on. The horizontal axis is time, where -D1 indicates the previous day, -D2 indicates two days prior, and so forth. The pacing percentage (%) refers to the percentage of pacing events in the total number of pacing and sensing events, which is calculated every three hours.

## RVP Configuration View


Pressing “FUNC” button enters this view. User can adjust the Rate / Amplitude parameters of RVP by pressing Rate/ AMPL Button and with

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
the use of the Dial. Meanwhile, the current pacing parameters are displayed in the bracket on the right side. Before and after delivering RVP, the device will deliver pulses according to the Rate and Amlitude parameters in the bracket on the right side .

RVP Popup Window

Hold “FUNC” button to delivery RVP and RVP Popup window will be displayed, showing the current Rate and Amplitude parameters during Rapid Atrial Pacing, the pacing mode is SOO, and the Rate or Amplitude parameter can be adjusted by rotating the Dial.

RAP Configuration View


Pressing “FUNC” button enters this view. User can adjust the Rate / Amplitude parameters of RAP by pressing Rate/ AMPL Button and with the use of the Dial. Meanwhile, the current pacing parameters are displayed in the bracket on the right side. Before and after delivering RAP, the TPM will deliver pulses according to the Rate and Amlitude parameters in the bracket on the right side .

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
Hold “FUNC” button to delivery RAP and RAP Popup window will be displayed, showing the current Rate and Amplitude parameters during Rapid Atrial Pacing, the pacing mode is SOO, and the Rate or Amplitude parameter can be adjusted by rotating the Dial.

Settings View

Pressing “FUNC” button enters this view. In this View, the Bluetooth connection function with the Pacing Parameters & EGM Data Storage and Transfer Software can be turned ON or OFF.

## 3-4 LEDs

There is a row of three LED indicators on the Temporary Pacemaker. They are centered along the top edge at the front side of the device. From left to right, they are “PACE” LED, “SENSE” LED and “POWER” LED indicators respectively.

During power-on self-tests, all three LEDs light up briefly to indicate device startup, and are dormant most of the time during nominal device operation.

“PACE” LED

This green colored Pace LED indicator is at the left. During normal operation,

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the pace LED flashes in synchrony to each pacing stimulus.

“SENSE” LED

This white colored Sense LED indicator is at the center. During normal device operation, it flashes in synchrony to each valid sense event outside the refractory window.

“POWER” LED

This orange colored power LED indicator is at the right. The Temporary Pacemaker continuously monitors the status of the primary battery, and flashes the power LED to notify the user if the batteries need to be replaced.

After startup, this indicator flashes if the remaining energy in the installed battery is low. User should replace the old batteries with 2 fresh size AA 1.5V Alkaline (LR6) batteries upon onset of the low battery indicator and ensure that the orange power LED is dormant.

## 3-5 Low Battery Indicator

The Temporary Pacemaker beeps and its orange-colored power LED flashes to notify user when the battery can only support the device for less than a day and needs immediate replacement. On the LCD screen, a Battery

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
Replacement message prompts the user until the batteries are replaced. Once the low battery warning first appears, the device can typically continue operating for approximately 24 hours.

## 3-6 Battery Compartment

The battery compartment is located at the base of the back side of the device. The compartment accepts two size AA 1.5V Alkaline (LR6) batteries. The compartment cover is secured by a rotary latch. In addition, the cover is tethered to the device body via a silicon cable.

To open the compartment:

1) Locate the key at the right side of the compartment and use it to unlock the rotary latch. On the right side of the battery cover, locate the battery compartment key slot. Gently pry the lower end of the battery compartment key upward and outward at the same time, then pull it backward slightly until the key is removed.
2) Use the key to rotate the latch of the battery cover until the arrow on the latch points to the unlock icon.
3) Hold the lower end of the battery cover, lift it upward and outward simultaneously until the battery cover is fully opened.

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To close the compartment

1) Align the battery cover with the battery compartment, gently push the battery cover upward to fully insert the silicone wire on it into the device body. Then press the battery cover downward lightly to fully place it back into the battery compartment.
2) Use the key to rotate the latch of the battery cover until the arrow on the latch points to the lock icon to lock the battery cover.
3) Push the key and its silicone wire into the battery compartment key slot, then press downward gently until the end of the key also enters the battery compartment key slot.

## 3-7 Hang Ring

The hang ring is located at the top of the back side of the device. User can use the hang ring to secure the device to, for example, a bed railing or an IV hanger. To use the hang ring, fold the ring over the top of the device until it clicks into position. When not in use, the ring should be kept in its recessed position against the back of the device.

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# 4 Device Operation

## 4-1 Preparation

The Temporary Pacemaker is intended to be setup and operated by trained professionals in a clinical environment. Use the supplied Patient Cable/ Analysis Cable with adapters/extension as needed to connect temporary pacing leads. The pacing lead system comprises the Temporary Pacemaker, the Patient Cable/Analysis Cable, implanted leads/temporary pacing leads, and the patient's myocardium. Follow the steps below and ensure secure connections before use.


1. Connect the patient cable to the device per Section 3.2 "Patient Cable Socket."
2. Connect the electrode lead (with adapter/extension if used) to the patient cable's electrode connector.


1) Turn the connector's clamping screw counterclockwise until the lead pin can be inserted.

2) Fully insert the lead pin into the connector until it cannot go further.

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3) Turn the clamping screw clockwise to secure the pin. Gently pull the lead to confirm it does not come loose.
4) Pay attention to polarity when connecting the patient cable. The electrode connector on the patient cable is marked with (+) and (–) polarity symbols. When using a bipolar pacing lead, connect the adapter/extension cable of the lead tip (distal electrode) to the (–) pole, and connect the adapter/extension cable of the lead ring (proximal electrode) to the (+) pole. When using a unipolar pacing lead, connect the cable from the cardiac electrode to the (–) pole, and connect the cable from the body surface electrode to the (+) pole

3. Connect the electrode lead (with adapter/extension if used) to the Analysis Cable's electrode connector.

   Pay attention to the polarity when connecting the Analysis Cable. The alligator clips of the analysis cable include a red (+) pole and a black (-) pole.

   1) If using a bipolar lead, connect the Adapter/Connector corresponding to the tip lead to the black (-) alligator clip, and the Adapter/Connector corresponding to the ring lead to the red (+) alligator clip

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2) If using a unipolar lead, connect the lead corresponding to the lead from the heart to the black (-) alligator clip, and the lead corresponding to the lead from the body surface to the red (+) alligator.

4. Recheck all connections between the device, Patient Cable/Analysis Cable, and leads to ensure secure and safe connections.

To disconnect the cable, follow Section 3.2 instructions.

## 4-2 Power-On

Pressing either the Power button or the Emergency button powers on the device immediately. A 3-second welcome screen is shown.

Note that do not operate outside the specified environmental conditions. If the device undergoes a sudden temperature change, allow it to reach ambient temperature before use.

Note that do not press any buttons during startup to avoid false key-fault detection.

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After the device is powered on, the LCD screen transitions into the home view and the device begins normal operation with the following default configurations.

Pacing Mode: SSI
Rate: 70 ppm
Amplitude: 5.0V
Pulse Width: 1.5ms
Sensitivity: 2.0mV

If the device detects an error during startup, it will delay the power-on process. In certain cases, the operator may see an error message (refer to Chapter 7 for how to handle error messages). Additionally, the device performs continuous real-time self-checks during operation to ensure proper functionality.

## 4-3 Power-Off the device

The Temporary Pacemaker employs a two-step power down process to safeguard the operator from accidentally turning the device off.

To power off the device:

1) Ensure the buttons are unlocked. If the buttons are locked, press the "🔓" button to unlock them.

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3) Press the Power button once again confirms the power-off: the device immediately begins shutting down.

In any case, if the operator wishes to exit the shutdown prompt screen and continue device operation, either of the following two methods can be used:

1) Take no action until the prompt disappears.

2) Press any key other than the "⏻" button and the " ⚠" button to exit the shutdown prom

## 4-4 Lock/Unlock

There are 10 buttons and 1 dial on the device. The device automatic locks all buttons and dial if not operated within 1 minutes (extended to 5 minutes in ANALYSIS View). Users can lock or unlock the buttons manually using the round Lock button located in the center of the dial.

Button Lock Status Indication

When the buttons are locked, a lock icon will appear in the information bar

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at the bottom of the display. While the keypad is locked, pressing any button (except the “ ⚠ ” button) will bring up a lock message in the center of the screen indicating that the buttons are locked. To unlock, press the “🔒” button.

The button lock supervises the following list of 6 buttons

- “⏻” button (Power button)
- “PAUSE” button
- “ZOOM” button
- “RATE” button
- “AMPL” button
- “PULSE WIDTH” button
- “SENS” button
- “FUNC” button
- Dial

The following 1 button are not supervised by the button lock.

- “ ⚠ ” button (Emergency button)

The “🔒” button in the center of the dial is used to toggle between locked and unlocked states.

## 4-5 Pacing Stimulus

Pulse Amplitude
User can set the Amplitude in discrete steps between 0.1V and 10.0V.

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To adjust Amplitude

1) Make sure the buttons are unlocked

2) Single press the Amplitude configuration button to get into the configuration view

3) To increase Amplitude, rotate the dial clock-wise, and to decrease Amplitude, rotate the dial counter-clock wise.

4) To exit out of configuration, press the Amplitude configuration button again to return to the “Home” View.

Pulse Width
User can set the Pulse Width in discrete steps between 60us and 2.0ms.

To adjust Pulse Width

1) Make sure the buttons are unlocked

2) Single press the Pulse Width configuration button to get into the configuration view

3) To increase Pulse Width, rotate the dial clock wise, and to decrease Pulse Width, rotate the dial counter-clock wise.

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4) To exit out of configuration, press the Pulse Width configuration button again to return to the “Home” View.

## 4-6 Synchronous Mode (SSI)

By default, the Temporary Pacemaker functions in the single chamber synchronous pacing mode. In this mode, the Pacing timing can synchronize to the patient’s intrinsic rhythm, and pacing is inhibited in the presence of patient’s intrinsic rhythm at a rate higher than the set Rate.

| First alphabet | Second alphabet | Third alphabet |
|---|---|---|
| Chamber paced | Chamber sensed | Response to sensing |
| O=None<br>S=Single | O=None<br>S=Single | O= Indicates pulse generator has no sensing function<br>T = Triggered<br>I=Inhibited |

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
2) Single press the Sensitivity configuration button to get into the configuration view.

3) Adjust the Sensitivity setting, using the dial, to a target numerical value other than “ASY”.

4) To exit out of configuration, simplify press the Sensitivity configuration button again to return to the “Home” View.

## 4-7 Asynchronous Mode (SOO)

In the asynchronous mode, the Temporary Pacemaker delivers pacing stimuli at a set fixed Rate independent of any intrinsic activity.

Warning. In the asynchronous pacing mode, pacing stimulus in the ventricle may collide with patient’s intrinsic T-waves, as a result, this could lead to onset of fatal arrhythmias. User should only utilize this mode of pacing with appropriate clinical and emergency supports.

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2) Single press the Sensitivity configuration button to get into the configuration view

3) Decrease or increase sensing sensitivity, by rotating the dial counter clock-wise, until Sensitivity threshold value becomes “ASY”.

4) To exit out of configuration, press the Sensitivity configuration button again to return to the “Home” View.

## 4-8 Rate

The pulse rate refers to the pacing pulse delivery frequency of the Temporary Pacemaker.
In SSI mode, the device begins timing from the most recent pacing or sensing event, searching for the patient’s intrinsic heartbeat within the preset interval (t = 60/pulse rate in seconds).
If cardiac activity is detected within this interval, the device inhibits pacing and restarts the timing cycle from that point. If no intrinsic activity is detected, a pacing pulse is delivered at the end of the interval. In asynchronous pacing mode, the device does not sense the patient’s intrinsic cardiac activity and delivers pacing pulses at the fixed preset rate.

User can configure the Rate in discrete steps between 30 ppm and 200 ppm.

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To configure the Rate,


1) Make sure the buttons are unlocked

2) Single press the Rate configuration button to get into the configuration view

5) To increase Rate, rotate the dial clock wise, and to decrease Rate, rotate the dial counter-clock wise.

3) To exit out of configuration, press the Rate configuration button again to return to the “Home” View .

## 4-9 Run-Away Protection

Run-Away protection is a standard pacemaker feature to protect the patient from being paced at a dangerously high rate due to malfunctions within the Pacing timing generator.

The Temporary Pacemaker implements the run-away protection through a three-way independent cross-check. With this rate protection in place, when single-fault induced run-away pacing occurs, the Rate is guaranteed to be limited to < 220 ppm.

Furthermore, upon the initial detection of a run-away malfunction, the device

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first attempts to recover normal function by resetting programmed configurations to power-on values (Section 6 Specification). If run-away condition still persists after the recovery attempt, the device first display an error message for 10 seconds and then powers off.

## 4-10 Sensing

In synchronous pacing mode, the Temporary Pacemaker uses a standard threshold detection method to identify P/R wave peaks. When the pacing lead is properly connected and implanted in the patient’s heart chamber, the device collects intracardiac ECG signals through the lead and compares them with the sensitivity threshold set on the device. When the collected EGM signal voltage exceeds the sensing sensitivity threshold, the device interprets it as a sensed P/R wave, that is, a sensing event is detected

## 4-11 Refractory period

The refractory period is initiated by a pacing or sensing event for a fixed duration of 250 ms. During the refractory period, sensing events are ignored and does not change the pacing timing.

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## 4-12 Noise reversion

The Temporary Pacemaker uses a rate based method to continuously monitor the signal integrity. In the presence of excess noise, the device automatically reverts to the asynchronous mode at a user-set Rate. The reversion is effective until the noise condition diminishes and its impact on the sensing channel is below the sensing sensitivity.

## 4-13 Sensing Sensitivity

1) Make sure the buttons are unlocked

2) Single press the Sensitivity configuration button to get into the configuration view


3) To increase or decrease sensing sensitivity, rotate the dial clock wise, and to decrease or increase sensing sensitivity, rotate the dial counter-clock wise.

4) To exit out of configuration, press the Sensitivity configuration button again to return to the “Home” View.

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User can mute the pacing output for up to 10 seconds at a time using the Pacing Pause feature. This feature can be used for briefly examining patient's intrinsic rhythm without changing existing device configurations.

When pacing is paused, user is notified via a prompt on the LCD screen and no other user input is possible prior to the termination of the pause. OSO mode refers to the mode that has no pacing pulse and only sensing event. Pause is in effect only when the pause button is pressed and held down.


<u>To pause pacing</u>

1) Make sure the buttons are unlocked

2) Press and hold the pause button down to mute the pacing output for up to 10 seconds.

3) Release the button to resume pacing output

## 4-15 Emergency Pacing

User can invoke the emergency pacing mode via a single button press, even when the device is powered off. The emergency pacing mode is implemented as a preset of pacing and sensing configurations.

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Pressing the emergency button immediately configures the device to operate with the following emergency presets:

Pacing Mode: SSI
Rate: 70 ppm
Amplitude: 10.0V
Pulse Width: 2.0ms
Sensitivity: 2.0mV

Note that, the emergency button is not supervised by the button lock, however it is protected, by a mechanical guard-ring around the button, from accidental activation.

Note that, invoking emergency pacing over-writes existing pacing and sensing configurations with the emergency preset. User can change these configurations at any time.

## 4-16 EGM Waveform

The Temporary Pacemaker features a real-time display of a longer than 3-second rolling EGM waveform in the four configuration views, and a longer than 1.5-second rolling EGM waveform in the “ANALYSIS” View.

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In the four parameter configuration views and the “ANALYSIS” View, the ECG amplitude scaling can be adjusted using the “Zoom” Button.
In the four parameter configuration views and the “ANALYSIS” View, the Temporary Pacemaker displays a Single-channel EGM waveform with a 0.5 Hz high-pass filter.

To view the longer than 3-second rolling waveform in a configuration view

1) Make sure the buttons are unlocked

2) Single press any of the four configuration buttons, i.e. for Rate, Amplitude, Pulse Width, and Sensitivity, to get into a configuration view.

3) To exit out of a configuration view, press the corresponding configuration button to return to the “Home” View.

To view the longer than 1.5-second waveform in the “ANALYSIS” View

1) Ensure the buttons are unlocked

2) Single press the “FUNC” button to get into the “ANALYSIS” View.

3) To exist out of the “ANALYSIS” View, simply press the “FUNC” button twice to rotate through the

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"PACING%" View and then return to the "Home" View .

To change the vertical zoom of the waveform

Press the "ZOOM" button to rotate through the different vertical scales of the waveform.

## 4-17 Event Markers

The Pacing and Sensing event markers as well as event interval could be shown in the four configuration views, the "ANALYSIS" View and the "EGM Waveform Comparison" View.

## 4-18 Analysis

To enter the "ANALYSIS" View:

1) Ensure the buttons are unlocked.
2) From the "Home" View, press the "FUNC" button once to get into the Analysis interface.

## 4-19 Pacing Impedance


The Temporary Pacemaker measures the Pacing impedance using actual Pacing pulses. Entering the "ANALYSIS" View starts the Pacing Impedance measurement, during which, the

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Amplitude and Pulse Width are temporarily changed to equal to 5.0V and 1.5ms respectively. When the device exits the “ANALYSIS” View, the Amplitude and Pulse Width will automatically return to the settings before entering the “ANALYSIS” View.

To measure pacing impedance, it is necessary to first enable the temporary pacemaker to deliver pacing pulses. In the SSI mode, if the pacing rate is lower than the patient’s intrinsic heart rate, the device will not deliver pacing pulses, and thus cannot measure pacing impedance. In this case, the prompt "Pacing Required" will be displayed next to the impedance item on the “ANALYSIS” View. If impedance measurement is needed under such circumstances, set the pacing pulse rate to a value higher than the patient’s intrinsic heart rate to allow the device to deliver pacing pulses, and then the impedance can be measured.


In the “ANALYSIS” View, Pacing Impedance measurement of a connected lead system is shown in the upper left on the LCD. The measurement range is from 200 ohms to 4000 ohms.

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
The Temporary Pacemaker measures the peak amplitude of a P or R wave after each sense.
To ensure the accuracy and repeatability of the measurement, please maintain 5 consecutive beats or more of sensing events.

In the “ANALYSIS” View, peak amplitude of a patient’s intrinsic waveform is displayed in the upper left on the LCD. The measurement range is 0.5 to 20.0mV with a 0.1mV resolution.

## 4-21 EGM Waveform Comparison View

In the “EGM Waveform Comparison” View, the Temporary Pacemaker displays two real-time updating EGM waveforms (>1.5 s) (high-pass filtered at 0.5 Hz and 30 Hz). Use the “ZOOM” button to adjust the EGM amplitude scaling.

In the “EGM Waveform Comparison” View, the Temporary Pacemaker displays Dual-channel EGM waveforms: F1 (0.5 Hz high-pass) and F2 (30 Hz high-pass).

Follow the steps below to observe the real-time updated EGM waveforms in the “EGM Waveform Comparison” View:


1. Confirm that the buttons are unlocked.

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
## 4-22 Pacing Percentage

From the “ANALYSIS” View, pressing the “FUNC” button again brings the device to the “PACING%” View. The graph displays a total of 5-days-worth of data after the last startup in the form of a 2D chart. The horizontal axis is time with each grid denoting one day: -D1 indicates the previous day, -D2 indicates the day before the previous day, and so on. The vertical axis is the Pacing percentage from 0 to 100 percent. The Pacing percentage data is cleared with zeros at each power-off.

## 4-23 Rapid Atrial Pacing


Press the "FUNC" button to switch to the related RAP feature views. The pacing mode is displayed in the left corner of the view and the battery level is displayed in the lower-right corner.

Before entering the “RAP Configuration” View, the RAP Prompting View will be shown firstly, indicating that only the atrium can perform Rapid

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

Atrial Pacing. Please make sure that the device is connected to the atrial lead. Press the "FUNC" button again to enter the “RAP Configuration” View.

When entering the RAP Configuration View, the device continues to operate according to the previously setting mode and parameters. User can adjust the Rate / Amplitude parameters of RAP by pressing Rate/ AMPL Button and with the use of the Dial. Meanwhile, the current pacing parameters are displayed in the bracket on the right side. Before and after delivering RAP, the device will deliver pulses according to the Rate and Amlitude parameters in the bracket on the right side.

After setting the RAP parameters, Hold “FUNC” button to delivery RAP and RAP Popup window will be displayed, showing the current Rate and Amplitude parameters during Rapid Atrial Pacing, the pacing mode is SOO, and the Rate or Amplitude parameter can be adjusted by rotating the Dial. Release “FUNC” button to stop delivering RAP. Hold “FUNC” button can deliver RAP continuously for up to 120 seconds. Hold “FUNC” button again if continuous RAP is required.

If the device continues to deliver RAP even after the “FUNC” button is released, press the “⏻” button (Power

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
button) to shup down the device to stop RAP. If RAP continues to be delivered, remove the batteries from the TPM. Return the TPM for service.


## 4-24 Rapid Ventricle Pacing

Press the "FUNC" button to switch to the related RVP feature views. The pacing mode is displayed in the left corner of the view and the battery level is displayed in the lower-right corner.

When entering the RVP Configuration View, the device continues to operate according to the previously setting mode and parameters. User can adjust the Rate / Amplitude parameters of RVP by pressing Rate/ AMPL Button and with the use of the Dial. Meanwhile, the current pacing parameters are displayed in the bracket on the right side. Before and after delivering RVP, the device will deliver pulses according to the Rate and Amlitude parameters in the bracket on the right side

After setting the RVP parameters, Hold the “FUNC” button to delivery RVP and RVP Popup window will be displayed, showing the current Rate and Amplitude parameters during Rapid Atrial Pacing, the pacing mode is SOO, and the Rate or Amplitude parameter can be

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
adjusted by rotating the Dial. Release “FUNC” button to stop delivering RVP. The interface will briefly display RVP Stopped, and the device will emit a short prompt tone. Hold “FUNC” button can deliver RVP continuously for up to 45 seconds. Hold “FUNC” button again if continuous RVP is required.

If the device continues to deliver RVP even after the “FUNC” button is released, press the “ ⏻ ” button (Power button) to shup down the device to stop RVP. If RVP continues to be delivered, remove the batteries from the device. Return the device for service

## 4-25 Bluetooth Setting

Press the “FUNC” Button to switch to the Bluetooth Setting View. Use the Dial to make a selection and press the “Zoom” Button to confirm. If no Bluetooth connection is established within 15 minutes, the device automatically disables the Bluetooth function.

## 4-26 Battery Status

The Temporary Pacemaker can operate continuously on a fresh pair of size-AA 1.5V alkaline (LR6) batteries

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
for 17 days at power-on default settings. The device monitors battery status during operation, and notifies user should batteries be replaced.
The low battery notifications presented to the user through means, include, visual: LCD status icon, LCD prompt, blinking LED, and, audio: beeper.

Note that, the beeper comes on when it is absolutely necessary to replace the batteries immediately to ensure continuing operation.

Note that, batteries should be replaced when any of the low battery notifications comes on.

<u>The low battery indications include:</u>

1) Battery status icon: the low battery status icon is shown at the right in the information bar at the bottom of the screen. The status icon shows graphically the battery status throughout the operation span of installed batteries, from full to empty.

2) Low battery prompt: user is prompted the low battery notification in the center of the screen when the battery is low.

3) Power LED: This orange colored power LED indicator is at the right amongst a row of 3 LEDs along the

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top edge at the front side of the device. The power LED flashes to notify the user when the batteries need to be replaced.

4) Low battery beeper: The TPM repeats a double beep pattern when it is absolutely necessary to replace the batteries immediately to ensure uninterrupted operation.

## 4-27 Battery Replacement

New batteries should be installed when the device is used on each new patient or for every two weeks, and/or when a low battery status indicator is on.

The Temporary Pacemaker can continue to operate no less than 80 seconds through internal backup power during battery replacement.

Warning: Battery Replacement is not recommended when the device is connected to a patient.

Warning: Do not use rechargeable batteries. Use only two new size AA 1.5 V alkaline (LR6) batteries.


Warning: Using a battery with a physical size not as recommended, or with contaminants on the battery terminals

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can result in unstable or no pacing output.

Note that, remove the batteries when the device is in storage.

The battery compartment is located at the base of the back cover of the device. The compartment accepts two AA sized (LR6) batteries. The compartment cover is secured by a rotary latch. In addition, the compartment cover is tethered to the device body via a silicon cable.

Open the battery compartment:

1) Locate the key at the right side of the compartment and use it to unlock the rotary latch.
2) Slightly pull down the battery compartment, then life the cover up.

Replace batteries:

1) Remove the used batteries from the compartment and dispose the batteries properly according to local environmental guidelines. If contaminations and corrosions are found on the contacts inside the compartment, please contact Lifetech Cardio Support.
2) Inside the battery compartment, the TPM provides visual guidance on installation polarity. Follow the visual cues to load the new batteries in place.

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
1) Position the cover into the compartment and align the guide rails, slide the cover upward until the cover is fully seated.
2) Use the key right beside the compartment to turn the rotary latch to the lock position.

Battery Cover Warning:

When the battery voltage is too low and/or the battery compartment cover is unlocked, a Battery Cover Warning message prompts the user to replace the battery and lock the battery compartment cover as soon as possible.

Note that, the battery cover shall be locked to switch to battery power.

## 4-28 Cleaning and Disinfecting

**Cleaning**

Clean the device and Patient Cable/Analysis Cable by a thorough wipe down using a damp cloth or sponge moistened with water. Prior to use, user should ensure that the device and Patient Cable/Analysis Cable are thoroughly dried, especially inside the connectors and connector blocks. Have a trained person to perform in-

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spection and verify the device and cable are functioning properly after each cleaning.

**Disinfecting**

Disinfect the device and Patient Cable/Analysis Cable by using a cloth or sponge moistened with 75% ethyl alcohol.

After disinfection, user should ensure that the device and Patient Cable/Analysis Cable are thoroughly dried, especially inside the connectors and connector blocks and have a trained person to perform inspection and verify the device and cable are functioning properly.

Note that, subjecting the device to liquid/disinfectant immersion and sterilization using gamma radiation or steam may damage the device and render the device unfunctional and unsafe to use.

## 4-29 Inspection

Have a trained person to perform inspection and verify the device and cable are functioning properly after each cleaning and disinfection:

- Perform inspection to ensure there is no physical damage or contamination on the device and cable.

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- Ensure that the device and Patient Cable/Analysis Cable are thoroughly dried, especially inside the connectors and connector blocks.

- Verify that all the displays, buttons and dial work and function properly.

- Verify that the cable and all connections are properly connected and ensure the device passes the self-test when powered on prior to use.

If any inspection indicates improper operation, the device shall be returned to Lifetech Cardio for repair.

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# 5 Analysis Application

## 5-1 Preparation

1) Ensure patient is under continuous electrocardiogram and hemodynamic monitoring.
2) Prepare the Temporary Pacemaker, power it on and confirm that the device is functioning properly with sufficient battery power.
3) Connect the Temporary Pacemaker to the patient in accordance with the instructions (4.1 Preparation)

## 5-2 Sensitivity Test

1) Set the Rate 10 ppm lower than patient's intrinsic rate in SSI mode, according to the instructions (4.8 Rate).
2) Observe whether the "SENSE" LED flashes continuously. If it fails to flash continuously, follow the instructions (4.13 Sensing Sensitivity) to decrease the sensing sensitivity setting value (thereby increasing the sensing sensitivity) until the "SENSE" LED is observed to flash continuously and/or sensing markers appear on the EGM waveform. At this point, confirm that the device can continuously sense the P/R waves.

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3) Press the “FUNC” button to switch the ANALYSIS View. Read the peak amplitude of the sensed P/R waves in the upper left on the LCD.
4) Set the Sensitivity value at half the measured peak amplitude of the sensed P/R waves . For example, if the sensed peak amplitude is 2.0mV, the recommended Sensitivity is 1.0mV.
5) Confirm that the device senses the P/R waves correctly by reviewing the EGM waveform.
6) Confirm that the current rate meets the clinical requirements. If adjustment is needed, follow the instructions (4.8 Rate) to make the corresponding settings.

## 5-3 Pacing Capture Threshold

1) Increase the Rate 10 ppm higher than the patient’s intrinsic rate in the Rate view, according to the instruction (4.8 Rate).
2) Keep the Pulse Width at 0.42ms, according to the Pulse Width section of “4.5 Pacing Stimulus”.
3) Increase the Amplitude in the Amplitude view until depolarization occurs with every Pacing spike, according to the Pulse Amplitude section of “4.5 Pacing Stimulus”. View the patient’s cardiograph to make sure that patient's heart rate increase to the setting Rate value. The Amplitude value is the capture threshold.

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4) Set the Amplitude to at least 2 times the determined capture threshold, according to the Pulse Amplitude section of “4.5 Pacing Stimulus”. For example, if the determined capture threshold is 2.0V, then recommended Amplitude is at least 4.0V.
5) Confirm that the device pacing function is normal by reviewing the EGM waveform.
6) Confirm that the current rate meets the clinical requirements. If adjustment is needed, follow the instructions (4.8 Rate) to make the corresponding settings.

## 5-4 Impedance Test

1) Press the “FUNC” button to switch the ANALYSIS View. Read the Impedance in the upper left on the LCD.
2) If display “Pacing Required”, increase the pulse rate. Observe the current pulse rate in the upper right corner of the interface until it is approximately 10 ppm higher than the patient's intrinsic rate.
3) Observe the pacing LED and EGM waveform, ensure pacing occurs from time to time.
4) Read the Impedance value.
5) Confirm that the current rate meets the clinical requirements. If adjustment is needed, follow the instructions (4.8 Rate) to make the corresponding settings.

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# 6 Pacing Parameters & EGM Data Storage and Transfer Software

## 6-1 Overview

The Pacing Parameters & EGM Data Storage and Transfer Software, when used in conjunction with the Lifetech Cardio Temporary Pacemaker, displays the pacemaker parameters, waveforms, and pacing percentage trend charts.

## 6-2 Software Download

Please download and install the software from the official website of Shenzhen Lifetech Cardio Medical Electronics Co., Ltd. before use.
**Download:** www.ltcardio.com
**System Requirements:**
Processor frequency: ≥ 2.3 GHz
Memory: ≥ 8 GB
Storage: ≥ 256 GB
Operating system: Windows 10 or compatible versions

## 6-3 End User License Agreement

The End User License Agreement (EULA) defines the terms and conditions for using the software. Users should read the agreement carefully. Accepting the agreement indicates acceptance of all terms. If the user

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does not agree, click the “Reject” Button or close the pop-up window

## 6-4 Create an ID

The “ID” Window is used to record a patient ID, which is used as the folder name for data storage. If no ID is entered, the software creates a folder name using the current system time. Personal or private information must not be used as the ID.

## 6-5 Connect

Press the “Connect” Button to scan for nearby Lifetech Cardio Temporary Pacemaker devices. Select the device from the pop-up window. A PIN code is required for first-time pairing. After a successful connection, press “Disconnect” to end the Bluetooth connection.

## 6-6 Display

**Bluetooth Connection Status**

The Bluetooth connection status is shown at the lower left of the main display view.

Status messages include Bluetooth Disconnected, Bluetooth Connection Failed, Start Scanning, Device Scan Finished, and Bluetooth Connected.

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Displayed parameters include Rate, Amplitude, Impedance, P/R Wave Amplitude, Sensitivity, and Pulse Width.

**Waveform Display**
The software plots received data as EGM waveforms, pacing markers, sensing markers, and marker intervals.

**Zoom Factor Display**
Supported zoom factors are 1X, 2X, 4X, 8X, and 16X

**Sweep Speed Display**
The sweep speed is 25 mm/s, and “X:200 ms/Div” indicates 200 ms per grid.

**EGM Waveform Comparison**

The EGM waveform comparison View is displayed automatically.

**PACING% View**
The pacing percentage View is displayed automatically.

## 6-7 Freeze

Press the “Freeze” Button to freeze the software display. The user can measure the time between two feature points. Right-click to delete a measurement. Press the “Unfreeze” Button to return to the non-frozen

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state.

**Information Bar**
The freeze time is displayed as “---” by default and updates to the system time after freezing. Event descriptions can be edited by double-clicking and confirmed by pressing Enter.

## 6-8 Settings

The storage path and center name (up to 20 characters) can be configured. Saved screenshots and recordings are stored in the selected path.

## 6-9 Save Screenshot

Screenshots are saved in PDF and PNG formats when sufficient storage space is available.

## 6-10 Browse Screenshots

Up to 50 screenshots can be browsed under the current ID.

## 6-11 Screen Recording

Displayed waveforms, markers, and marker interval information are saved in segments during recording.

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# 7 Specification

## 7-1 Parameters

Unless otherwise specified, all parameters below refer to measurements taken at 20 °C ± 2 °C, connected to a 500 Ω (±1%) load.

Service life: 8 years.

The Temporary Pacemaker has no replaceable internal parts. Any faults can be returned for factory service. Lifetech Cardio guarantees availability of repair parts for at least 10 years from delivery.

**Device Dimensions**

The physical dimensions of the Temporary Pacemaker are as follows (with the hanging ring retracted):

| | Dimensions ( mm ) | Tolerance ( mm ) |
|---|---|---|
| Length | 186 | ±10 |
| Width | 86 | ±8 |
| Thickness | 33 | ±5 |

**Patient Cable Dimensions**

The dimensions of the Patient Cable accessory are as follows:

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| | Dimensions ( mm ) | Tolerance ( mm ) |
|---|---|---|
| Length | 1860 | ±50 |
| Maximum diameter of compatible adapter connectors | 3 | NA |

**Analysis Cable**

The dimensions of the Analysis Cable accessory(optional) are as follows:

| | Dimensions ( mm ) | Tolerance ( mm ) |
|---|---|---|
| Length | 1860 | ±50 |
| Size Diameter of the opening of the Alligator Clip(mm) | >5.0 | NA |

**Storage Conditions (Batteries removed)**

| Temperature °C | Altitude (m) | Humidity (%) |
|---|---|---|
| Min: -20<br>Max: 50 | ⩽4000m | < 90 |

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**Battery**

The Temporary Pacemaker uses two AA (LR6) alkaline batteries.

1) At an ambient temperature of 20°C, the Temporary Pacemaker is expected to operate continuously for at least 14 days using two new AA (LR6) alkaline batteries.

2) Under the same temperature condition (20°C), using default pacing and sensing parameters (pacing rate: 70 ppm, pulse width: 1.8ms, pulse amplitude: 5.0 V, sensing sensitivity: 2.0 mV), the typical continuous operating time is up to 20 days with two new AA (LR6) alkaline batteries.

3) Different battery brands and models may affect performance. It is recommended to use AA (LR6) alkaline batteries such as those from Nanfu.

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| Battery Type | Battery Life(days) |
| --- | --- |
| AA (LR6) alkaline battery | Typically:20; Mnimum:14 |

### Backup Power Supply

The Temporary Pacemaker can rely on backup power to ensure at least 80 seconds of continuous operation.

### Software Release Version

The software release version of the Temporary Pacemaker is Version 3.

The released version of the Pacing Parameters & EGM Data Storage and Transfer Software for the Lifetech Cardio Temporary Pacemaker is Version 1.

### Pacing Modes

According to clinical requirements, the Temporary Pacemaker can operate in either Synchronous Pacing Mode (SSI) or Asynchronous Pacing Mode (SOO).

### Chamber

The Temporary Pacemaker can pace either the atrium or the ventricle independently.

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### Runaway Protection Rate

The runaway protection mechanism of the Temporary Pacemaker is triggered at rates exceeding 220 ppm.

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| | |
|---|---|
| Pulse Width | 2.0ms |
| Sensitivity | 2.0mV |

**Rate**

| Range (ppm) | Toler-ance (%) | Step Size (ppm) |
|---|---|---|
| Min: 30<br>Max: 200 | +/- 5 | 5, from 30 to 50<br>2, from 50 to 100<br>5, from 100 to 200 |

**Amplitude**

| Range (V) | Toler-ance (V or %) | Step Size (V) |
|---|---|---|
| Min:0.1<br>Max: 10.0 | +/-0.1V, from 0.1 to 0.9 V.<br>+/-10%, from 1.0 to 10.0 V. | 0.1, from 0.1 to 1.0<br>0.2, from 1.0 to 3.0<br>0.5, from 3.0 to 10.0 |

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**Sensitivity**

| Range (mV) | Tolerance (%) | Step Size (mV) |
|---|---|---|
| Min: 0.5<br>Max: 20 | +/-20, from 0.5 to 1.0mV.<br>+/-10, from 1.2 to 20.0mV. | 0.1, from 0.5 to 1.0<br>0.2, from 1.0 to 3.0<br>0.5, from 3.0 to 10.0<br>1.0, from 10.0 to 20.0 |

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## 7-2 Standards Compliance

Safety Standard: GB 9706.1-2020 Medical Electrical Equipment – Part 1: General Requirements for basic safety and essential performance.

EMC Standard: YY 9706.102-2021 Medical Electrical Equipment – Part 1-2: General requirements for basic safety and essential performance- Collateral Standard: Electromagnetic Compatibility – Requirements and Tests;

Product Standard: YY 9706.231-2023 Medical Electrical Equipment – Part 2-31: Particular Requirements for the Safety and essential performance of external cardiac pacemakers with internal power source.

Essential performance of the device: Based on GB 9706.1-2020 and YY 9706.231-2023, the essential performance of this product is as follows:

1. Battery depletion indicator

2. ME equipment parameter stability

3. Pulse amplitude stability

4. Disabling runaway rate protection

5. Deliberate action required to change settings

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8. Interference reversion in the presence of sensed electromagnetic disturbance or electrical energy sources

Safety Features of the Device:

Classification by electric shock protection type: Internal power supply

Classification by degree of protection against electric shock: CF type application

Classification by degree of protection against liquid ingress: IPX0

Classification by safety when used with flammable anesthetic gases mixed with air or with flammable anesthetic gases mixed with oxygen or nitrous oxide: Not applicable

Rated voltage and frequency: internal power supply DC 3V

Input power: Two AA 1.5V alkaline batteries (LR6)

Application section with protection against defibrillation discharge: CF type application

Signal output or input section: EGM input, pacing signal output, and data output

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# 7-3 Electromagnetic Compatibility Statement

| **Guidance and manufacturer's declaration - electromagnetic emission** | | |
|---|---|---|
| The Temporary Pacemaker is intended to be used in the electromagnetic environment specified below, and the purchaser or user shall ensure it is used in this electromagnetic environment: | | |
| **Emission test** | **Compliance** | **Electromagnetic environment - guide** |
| RF Emissions<br>GB 4824 | Group 1 | The Temporary Pacemaker uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. |
| RF Emissions<br>GB 4824 | Class A | The Temporary Pacemaker is suitable for use in all establishments other than domestic households and those directly connected to the public low-voltage power supply network. |
| Harmonic emissions<br>GB 17625.1 | Not applicable | |
| Voltage fluctuations/flicker emissions<br>GB 17625.2 | Not applicable | |

| **Guidance and Manufacturer's Declaration - Electromagnetic Immunity** | | | |
|---|---|---|---|
| The Temporary Pacemaker is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment. | | | |
| **Immunity test** | **IEC 60601 test level** | **Compliance level** | **Assessment** |

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| | | | |
|---|---|---|---|
| Electrostatic discharge GB/T 17626.2 | ± 6 kV contact discharge ± 8 kV air discharge | ± 8 kV contact discharge<br>± 15 kV air discharge | Floors should be wood, concrete or ceramic tile. If covered with synthetic material, the relative humidity should be at least 30%. |
| Electrical Fast Transient/Burst GB/T 17626.4 | Not Applicable | Not Applicable | / |
| Surge GB/T 17626.5 | Not Applicable | Not Applicable | / |
| Voltage Dips, Short Interruptions and Voltage Variations on Power Supply Input Lines GB/T 17626.11 | Not Applicable | Not Applicable | / |
| Power Frequency Magnetic Field (50 Hz/60 Hz) GB/T 17626.8 | 3 A/m | 3 A/m | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment |
| Conducted RF GB/T 17626.6<br>Radiated RF GB/T 17626.3 | 3 V (rms)<br><br>150 kHz to 80 MHz (excluding | 3 V (rms) | Portable and mobile RF communica- |

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| | | | |
|---|---|---|---|
| | ISM bands)<br><br>10 V (rms)<br>150 kHz to 80 MHz (for ISM bands)<br><br>10 V/m<br>80 MHz to 2.5 GHz | 10 V (rms)<br><br>10 V/m | tions equipment should be used no closer to any part of the 6108 Temporary Pacemaker, including cables, than the recommended separation distance calculated from the applicable formula. Recommended separation distance:d = [1.17]√Pd = [1.20]√Pd = [1.20]√P (80 MHz – 800 MHz)d = [2.30]√P (800 MHz – 2.5 GHz)Where:P = max rated output power of the transmitter in watts (W)d = recommended separation distance in |

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| | | | |
|---|---|---|---|
| | | | meters (m) The field strengths of fixed RF transmitters should be determined by an electromagnetic site survey and should be less than the applicable compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol.  |
| Note 1: At 80 MHz and 800 MHz, the separation distance formula for the higher frequency range applies.<br>Note 2: These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. | | | |
| a) ISM bands between 150 kHz and 80 MHz are: 6.765‑6.795 MHz, 13.553‑13.567 MHz, 26.957‑27.283 MHz, and 40.66‑40.70 MHz.<br>b) An additional factor of 10/3 is used to reduce the likelihood of interference from mobile/portable communications equipment inadvertently brought into the patient area. | | | |

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c) The field strength from fixed transmitters (e.g., base stations for radio, TV, cordless phones, etc.) is not predictable. A site survey is recommended.
d) In the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended Separation Distances Between Portable and Mobile RF Communication Equipment and the Temporary Pacemaker**

The Temporary Pacemaker is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The purchaser or user of the Temporary Pacemaker can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the Temporary Pacemaker, according to the maximum output power of the communication equipment, as recommended below.

| Maximum Rated Output Power of Transmitter (W) | Corresponding Separation Distances by Frequency (m) | | | |
|---|---|---|---|---|
| | 150 kHz～80 MHz (excluding ISM bands) d=[1.17]√P | 150 kHz～80 MHz (excluding ISM bands) d=[1.20]√P | 80 MHz～800 MHz d=[1.20]√P | 800 MHz～2.5 GHz d=[2.30]√P |
| 0.01 | 0.12 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.37 | 0.38 | 0.38 | 0.73 |
| 1 | 1.17 | 1.20 | 1.20 | 2.30 |
| 10 | 3.70 | 3.79 | 3.79 | 7.27 |
| 100 | 11.70 | 12.00 | 12.00 | 23.00 |

For transmitters rated at maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter. P is the maximum rated output power of the transmitter in watts (W), as specified by the manufacturer.
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency

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range applies.

Note 2: The ISM frequency bands between 150 kHz and 80 MHz are 6.765‑6.795 MHz, 13.553‑13.567 MHz, 26.957‑27.283 MHz, and 40.66‑40.70 MHz.

Note 3: An additional factor of 10/3 is used when calculating the recommended separation distance for transmitters in the ISM bands between 150 kHz and 80 MHz, and in the frequency range 80 MHz to 2.5 GHz, to reduce the likelihood of interference when mobile/portable communication devices are brought into the patient area inadvertently.

Note 4: These guidelines may not apply in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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## 7-4 Model and Specification

Lifetech Cardio Temporary Pacemaker have six models with different feature configuration described in section 1-3. Lifetech Cardio Temporary Pacemaker Model 8400,8301,8200,8201,7101,7100 have same Pacing Features and Safety Features, and the configuration differences in Functional Features are shown as in the table below:

| **Functional Features** | 8400 | 8301 | 8201 | 8200 | 7101 | 7100 |
|---|---|---|---|---|---|---|
| Online EGM and Marker | √ | √ | √ | √ | √ | √ |
| RAP | √ | √ | √ | √ | √ | √ |
| Pacing Percentage | √ | √ | — | — | √ | — |
| Impedance | √ | √ | √ | √ | — | — |
| Peak Amplitude | √ | √ | √ | — | — | — |
| RVP | √ | — | — | — | — | — |
| EGM Waveform Comparison | √ | — | — | — | — | — |
| Pacing Parameters & EGM Data Storage and Transfer Software | √ | — | — | — | — | — |

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# 8 Support, Disclaimer and Warranty

## 8-1 Support

Lifetech Cardio has a professional technical team to provide operators with technical support and after-sales service. Contact information is provided on the front page. Each Temporary Pacemaker has a serial number on the back; please provide this information when requesting after-sales service or technical support.

**Fault Handling**

In case of equipment failure or accidents during proper use (e.g., fluid ingress, dropping), which do not harm the operator, please replace the device with another Temporary Pacemaker to continue pacing the patient. Turn off the faulty device, attempt to restore its original condition, and restart it after waiting ten seconds. If the issue persists, please contact Lifetech Cardio after-sales service for technical support. Contact information is provided on the front page.

**Device Repair**

Do not attempt to repair the device yourself. For repair, please contact

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Lifetech Cardio after-sales service. Contact details are provided on the front page.

**Professional Training**

The Lifetech Cardio after-sales team is dedicated to providing training on the use of the Temporary Pacemaker for operators in need.

## 8-2 Disclaimers

**Temporary Pacemaker**

The Temporary Pacemaker is a temporary cardiac pacing and implant analysis device that offers specified anti-bradyarrhythmia pacing, sensing and analysis functions to support physicians in treating patients with bradyarrhythmia and heart block conductions. However, Lifetech does not guaranty the use of this device will effectively manage cardiac conditions for all patients.

**Patient Cable/Analysis Cable**

All cables and wires are subject to fatigue and damage throughout use, and Lifetech Cardio's Patient Cables are not immune to improper handling and repeated stress in varying use environments. Therefore, Lifetech Cardio makes no warranty that the cables

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will be free from malfunctions throughout use.

Although the Lifetech Cardio Patient Cable/Analysis Cable has been carefully designed, manufactured and tested prior to sales, the product may not perform its design functions satisfactorily for various reasons. The warnings on the product label provide more detailed information and are an integral part of the disclaimer. Therefore, Lifetech Cardio denies all express or implied warranties with respect to this product. Lifetech Cardio will not be liable for any incidental and consequential events arising out of or in connection with the use, defects or malfunctions of the Patient Cable/Analysis Cable, whether or not such claims are based on warranties, contracts or civil torts.

The above disclaimers are not aiming to and shall not be deemed to be inconsistent with the mandatory provisions of the applicable laws. If any part or provision of the limited warranty statement is judged to be unlawful, unenforceable or inconsistent with the applicable laws by the court with effective jurisdiction, the validity of the remainder shall not be affected, and all rights and obligations shall be recognized and fulfilled as if the limited warranty statement does not contain the individual part that has been found to be invalid.

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# 8-3 Warranty

**General**

The Temporary Pacemaker is warranted against defects in material and workmanship for a period of one (1) year from the date of delivery.

If you believe the device has a defect in material or workmanship or has failed during normal use within the warranty period, please contact Lifetech Cardio Support according to information provided on the title page for further assistance.

Lifetech Cardio can, at its sole discretion:

1) Repair the defective device in a Lifetech Cardio designated facility, or

2) Exchange with the customer an equivalent functioning product.

**Limitation**

This warranty is honored only when the following conditions are met:

1) The device is not tampered with, modified or repaired by personnel who is not qualified by Lifetech Cardio.

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2) The defect is not a result from customer's misuse, mishandling, and accident.

3) The defect is first communicated to, and the device is returned to Lifetech Cardio within 60 days of the discovery of the defect.

**Disclaimer**

The warranty statement is limited to the terms expressly stated, no other warranties are expressed or implied. Lifetech Cardio will not be liable for any incidental or consequential events, arising out of or in connection with the use, defects or malfunctions of the device, whether or not such claims are based on warranties, contracts or civil torts

The above disclaimers and limitations shall not be deemed to be inconsistent with the mandatory provisions of the applicable laws. If any part or provision of the warranty statement is judged to be unlawful, unenforceable or inconsistent with the applicable laws by the court with effective jurisdiction, the validity of the remainder shall not be affected, and all rights and obligations shall be recognized and fulfilled as if the warranty statement does not contain the individual part that has been found to be invalid.
---END---

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