Dual-Chamber Temporary Pacemaker

User Manual

Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

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Dual-Chamber Temporary Pacemaker

Manufacturer Information

Registrant /Manufacturer
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

Registrant/Production Location
1401, Building D1, Nanshan i Park
No.1001, Xueyuan Avenue, Taoyuan Street
Nanshan District
Shenzhen, Guangdong, China

Tel: +86 755 23251894
WWW.LTCARDIO.COM

Device Information and Support

Product Name:
Dual-Chamber Temporary Pacemaker

Model: 9102, 9101, 9100

~~Registration Certificate Number:~~
~~国械注准 20233120942~~

~~Product Technical Requirements Number:~~
~~国械注准 20233120942~~

~~Medical Device Production License Number:~~
~~粤食药监械生产许 20214401 号~~

After-sales Service Provider:
Shenzhen Lifetech Cardio Medical Electronics Co., Ltd.

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Document Number: LTC/TS/117I-01C Revision Date:2025-06-23 Version: 02.01

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4-8 Sensing Sensitivity 42

4-9 Atrioventricular (AV) Interval 43

4-10 Post-Ventricular Atrial Refractory Period (PVARP) 44

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

4-12 Pacing Mode Selection 45
4-13 Noise Reversion 46
4-14 Rate Limiting (Runaway Protection) 46
4-15 Emergency Pacing 47
4-16 Pacing Pause 48
4-17 Ventricular Safety Pacing 48
4-18 Wenckebach Response 49
4-19 Mode Switch 49
4-20 Intracardiac Electrogram (EGM) 49
4-21 Pacing Percentage 50
4-22 AV Conduction Time 50
4-23 Lead Impedance Monitoring (Configuration & Alerts) 51
4-24 Lead Impedance Trends 53
4-25 Rapid Atrial Pacing (RAP) 53
4-26 Battery Status 55
4-27 Replacing Batteries 56
4-28 Cleaning and Disinfection 57
4-29 Post-Cleaning/Disinfection Checks 58
5 Device Specifications 59
5-1 Parameter 59
5-2 Standards Compliance 67
5-3 EMC Declaration 69
5-4 Product Models 76
6 Support, Disclaimer, and Warranty 77
6-1 Support 77
6-2 Disclaimer 78
6-3 Warranty Information 79

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
Humidity limit (0%-90%RH)

General warning sign

Refer to instruction manual/ booklet

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
SN	Serial number
	Dul-Chamber Temporary Pacemaker
	Defibrillation-Proof CF-Type Applied Part
LOT	Batch Code
PN	Material code
	Patient Cable
	Press both battery compartment buttons simultaneously to open

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1 Device Overview

1-1 Description

The Lifetech Cardio Dual-Chamber Temporary Pacemaker consists of a Non-implantable Pulse Generator and a Patient Cable.

The device connects to one end of the patient cable, and the other end of the patient cable connects directly or via an adapter to pacing electrode leads. The pacing system consists of the dual-chamber temporary pacemaker, the patient cable, the electrode leads, and myocardial tissue.

The Lifetech Cardio Dual-Chamber Temporary Pacemaker is primarily used to treat bradycardia. It provides nine selectable pacing modes, including DDD, DOO, DDI, DAI, DVI, AOO, AAI, VOO, and VVI.

The pacing rate, atrial and ventricular pacing pulse amplitude, and sensing sensitivity are adjustable. In addition, the device can provide high-frequency pacing pulses up to 1000 bpm in asynchronous mode for the treatment of atrial tachyarrhythmias.

Dual-Chamber Temporary Pacemaker should only be used by trained professionals.

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1-2 Package Content

The product package contains:
1) Dual-Chamber Temporary Pacemaker *1
2) Patient Cable *1
3) AA (LR6) alkaline batteries *2
4) Dual-Chamber Temporary Pacemaker Carrying Case *1
5) Instructions for Use of the Dual-Chamber Temporary Pacemaker *1
6) Certificate of Conformity *1
The manufacturing date of the product can be found on the device label.

1-3 Product Functions

The functional features of the Lifetech Cardio Dual-Chamber Temporary Pacemaker are summarized as follows.

Pacing Function

Pacing Modes	DDD, DOO, DDI, DAI, DVI, AOO, AAI, VOO, VVI.
Pulse Width	0.8 ms
Amplitude	Adjustable, range: 0.1 – 20.0V
Rate	Adjustable, range: 30 – 200 ppm
Sensitivity	Adjustable, range: 0.5 – 20.0mV

The basic functions of the device are: (1) pacing and (2) sensing.
Note: During electromagnetic interference, a temporary loss of pacing output is permissible, but the duration

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must be less than 5 seconds.

Other Features

Intracardiac Electrogram (EGM)
The device can continuously display atrial and ventricular intracardiac electrograms (EGM) of more than 3 seconds, with pacing and sensing events annotated.

Pacing Percentage
The device can display a trend chart of pacing percentages over 5 days.

Atrioventricular (AV) Conduction Time
The device can display the patient's intrinsic AV conduction time trend and ventricular pacing percentage over 5 days in dual-chamber mode.

Lead Impedance Monitoring Alerts
The device provides monitoring and alerts for high and low lead impedance.

Lead Impedance Trend
The device can separately display atrial and ventricular impedance trends over 5 days.

Rapid Atrial Pacing (RAP)
The device can perform rapid atrial pacing based on preset parameters, with adjustable pulse rate and/or amplitude, for up to 120 seconds continuously.

Auto-adjustment Function
PVARP, maximum tracking rate, and AV interval can automatically adjust according to the basic pacing rate.

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Safety Features

Noise Reversion

When electrical interference is sensed, the device reverts to asynchronous pacing at the preset rate until interference stops.

Rate Limiting (Runaway Protection)

In the event of a single fault, the device ensures that the pacing output rate does not exceed 230 ppm.

Emergency Pacing

Pressing the "⚠" button once activates emergency pacing. Parameters: Mode: DOO; Rate: 80 ppm; Amplitude: 10.0 V; Width: 0.8 ms; Sensitivity: ASY.

Pacing Pause

Holding the "PAUSE" button suspends pacing for up to 10 seconds. Pacing automatically resumes when released.

Defibrillation Protection

The device has defibrillation protection to minimize malfunction risk during external defibrillation.

Low Battery Indicator

When the battery power is low, the "Power" indicator light will flash, the buzzer will sound an alert, and the screen will display the message "Low battery, please replace the battery." The battery level indicator will also show a low battery status. Once the device begins alerting due to low battery, it can typically continue operating for at least 24 more

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hours.

Battery Compartment Not Closed Alert

When the battery compartment is opened, the screen of the Dual-Chamber Temporary Pacemaker will display the message: "Please replace battery and lock its cover ASAP."

Backup Power

After the battery is removed, the backup power supply can provide continuous power for no less than 50 seconds, accompanied by the flashing "Power" indicator and a screen prompt.

Auto-Lock Screen

In the intracardiac EGM real-time display interface, the device will automatically lock the screen after 5 minutes of inactivity. In other interfaces, the screen will automatically lock after 1 minute of inactivity. When the device is in the auto-lock state, pressing any key other than the "⚠" key or the "🔒" key will prompt a message on the screen indicating that the device is locked.

Shutdown Protection

When the main unit is operating normally (not in a locked state), pressing the "⏻" key will cause the display to show a shutdown confirmation message.

If the operator presses the "⏻" key again before the message disappears, the device will shut down; otherwise, the shutdown operation will be canceled.

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This function is designed to prevent inappropriate inhibition of ventricular pacing pulses when the ventricular lead senses non-ventricular events (such as crosstalk or noise). It is active when the pacemaker operates in DVI, DDI, and DDD pacing modes.

Wenckebach Response

In DDD mode, when the intrinsic atrial rate exceeds the maximum tracking rate, the pacemaker will continue to extend the SAV interval until the P-wave falls within the PVARP, making it undetectable. Since the P-wave is not sensed, the SAV interval is not initiated, and no ventricular pacing pulse is delivered.

Mode Switching

Upon detecting atrial tachycardia, the device temporarily switches from the DDD pacing mode to the non-tracking DDI pacing mode to prevent rapid ventricular pacing caused by the elevated atrial rate. The previously set pacing mode will be restored once the atrial tachycardia terminates.

Atrial Tracking Functions

Atrial tracking is only accessible or applicable when the temporary pacemaker is set to DDD pacing mode. When atrial tracking is enabled, each sensed event on the atrial lead not only inhibits the scheduled atrial pacing pulse but also triggers the atrioventricular interval. When atrial tracking is disabled (DDI pacing mode), atrial sensing will inhibit atrial pacing but will not trigger the atrioventricular interval; ventricular pacing will occur at the selected rate.

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1-4 Intended Use

It is for single or dual-chamber external temporary pacing with Temporary Cardiac Pacing leads.

1-5 Indications

Examples of indications for the Lifetech Cardio Dual-Chamber Temporary Pacemaker include:

1) Sick sinus syndrome
2) Sinus bradycardia
3) Atrial and/or ventricular arrhythmias
4) Complete atrioventricular block
5) Asystole
6) Bradycardia accompanied by congestive heart failure
7) Patient support, management and evaluation before the implantation of implantable pacemaker
8) Support during the replacement of implantable pacemaker
9) Cardiac complications occurring during intervention or surgery
10) Support after the cardiac surgery
11) Acute myocardial infarction complicated with cardiac conduction block
12) High-rate burst pacing for the treatment of supraventricular tachyarrhythmias

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1-6 Contraindications

This product is intended for temporary cardiac pacing and currently has no known contraindications.

1-7 Patient Cable

The Patient Cable provided in the accessory kit can be used for connection with the Dual-Chamber Temporary Pacemaker for temporary pacing. The cable is approximately 1860 mm in length, with one end featuring a plug and the other end equipped with two electrode connectors.
The electrode connectors can be used to connect electrode leads or adapter cables with a connector diameter not exceeding 3.0 mm, and the plug is designed to connect to the patient cable socket located on the top of the Dual-Chamber Temporary Pacemaker.

To obtain a Patient Cable, please contact Lifetech Cardio. Contact information can be found on the front page.

Accessory	Order Code
Patient Cable	0203001001

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2-1 Warnings

Daily Use and Maintenance

Device Modification: Do not disassemble or attempt to repair the Lifetech Cardio Dual-Chamber Temporary Pacemaker yourself. Any disassembly or repair not performed by Lifetech Cardio may cause device malfunction or compromise safety.

Maintenance: To maintain its specified performance, the device requires regular preventative maintenance. Note: Temporary pacemakers and their accessories should not be maintained or serviced while in patient use.

Fluid Protection: Operators may wipe the device with a cloth dampened with clean water. Do not immerse the device in any liquid. Any liquid ingress may damage functionality and safety.

Fire/Explosion Hazard: Improper use of the Lifetech Cardio dual-chamber temporary pacemaker may generate sparks and must not be used in environments containing flammable or explosive gases, such as anesthetic gases, oxygen, or nitrous oxide.

Monitoring: The pacemaker is not intended for standalone patient monitoring. Use an ECG monitor or 12-lead ECG machine for appropriate patient

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monitoring.

Diagnostic Capability: The device has no diagnostic functions.

Not MRI-Compatible: Magnetic Resonance Imaging (MRI) is a medical imaging technique that uses magnetic fields to generate internal views of the body. Do not perform MRI scans on patients using the Lifetech Cardio Dual-Chamber Temporary Pacemaker. MRI scanning may induce tachycardia, cause serious bodily harm, or result in system malfunction or damage.

Usage: Do not place the pacemaker in close proximity to or stacked with other equipment. If unavoidable, verify normal operation in that configuration. Using accessories, transducers, or cables other than those provided by the manufacturer may increase emissions or reduce immunity.

Disposal: Dispose of the device according to local regulations.

Clinical Use

Asynchronous Pacing: In asynchronous pacing mode, ventricular pacing pulses may coincide with the patient's intrinsic T-wave, potentially triggering life-threatening arrhythmias. Clinical judgment must be exercised when using asynchronous pacing (e.g., in patients without intrinsic rhythm or those

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fully pacemaker-dependent), and appropriate pacing energy should be used. Asynchronous pacing is contraindicated when intrinsic cardiac rhythm is present.

Atrial Pacing: In AV conduction disorders or atrial fibrillation/flutter, single-chamber atrial pacing is ineffective and should not be used.

Atrial Sensing: During rapid atrial arrhythmias (e.g., AF/AFL), atrial sensing modes that trigger ventricular responses are contraindicated.

Pacing Discontinuation: Abrupt pacing termination may cause asystole. Before stopping, ensure synchronous mode and gradually reduce rate below intrinsic rhythm until recovery.

Introduced Current: The implanted leads form a direct electrical connection to the patient's heart, so any introduced current may be conducted directly to the heart, potentially causing fatal arrhythmias. Additionally, introduced current may lead to device oversensing, resulting in inappropriate inhibition of pacing pulses and loss of capture. Finally, introduced current may also damage the device.

Patient Cable: Always connect the cable to the device before attaching leads. Non-original cables may increase mis-sensing or reduce pacing efficacy.

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Electrostatic Discharge (ESD): To prevent the effects of static electricity on the device and the patient, the operator should discharge any static electricity from their body before using the device, such as by touching a large, grounded conductive surface.

Implanted Leads: Prevent electrode tips from contacting conductive or wet surfaces.

Mains-powered Equipment: When using 220V AC-powered equipment, ensure proper grounding.

Electrosurgical Unit: Keep the current input end of the electrosurgical unit as far away as possible from the electrode tip of the pacing lead, typically no less than 15 cm.

Ablation: Maintain as much distance as possible between the ablation electrode and the electrode tip of the implanted pacing lead, typically no less than 15 cm. Avoid direct contact between the ablation electrode and the electrode tip of the implanted pacing lead.

Defibrillation: Keep the defibrillation electrode pads as far away as possible from the electrode tip of the implanted pacing lead, typically no less than 15 cm. After each defibrillation, monitor the ECG and confirm that the Lifetech Cardio Dual-Chamber Temporary Pacemaker is functioning

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properly.
Please note that when defibrillation is performed while the Lifetech Cardio Dual-Chamber Temporary Pacemaker is connected, part of the defibrillation energy may be shunted through the electrical pathway between the implanted pacing lead and the myocardium, potentially reducing the therapeutic effectiveness of defibrillation.
Note: Whenever possible, disconnect the Lifetech Cardio Dual-Chamber Temporary Pacemaker from the patient and/or implanted lead prior to defibrillation.

RF Communication Devices: Keep handheld RF communication equipment, including antenna cables and external antennas, as far away as possible from the patient cable of the Lifetech Cardio Dual-Chamber Temporary Pacemaker. A minimum distance of 30 cm is generally recommended. Failure to maintain this distance may result in device malfunction or damage.

RF Immunity: The Lifetech Cardio Dual-Chamber Temporary Pacemaker is only resistant to specific designated RF signals. Therefore, operating the device in proximity to RF transmitters may result in abnormal device performance.

Electromagnetic Interference (EMI): Strong electromagnetic interference may cause the Lifetech Cardio Dual-Chamber Temporary Pacemaker to misinterpret signals, resulting in al-

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tered output. In such cases, it is recommended to switch the device to SOO mode with a pacing rate slightly higher than the patient's intrinsic heart rate.
Note: The real-time endocardial electrogram (EGM) display function can assist the operator in determining whether strong electromagnetic interference is present in the environment.

Patient monitoring: Continuously monitor the patient while the temporary pacemaker is in use to ensure proper functioning of the temporary pacemaker and to provide the patient with appropriate treatment.

2-2 Precautions

Lead Polarity: Bipolar pacing leads offer greater resistance to interference compared to unipolar leads. Therefore, when using the Lifetech Cardio Dual-Chamber Temporary Pacemaker, bipolar pacing leads are recommended whenever possible.

Emergency Pacing: In emergency pacing mode, the Lifetech Cardio Dual-Chamber Temporary Pacemaker operates with the following preset parameters: DOO mode, pulse width of 0.8 milliseconds, pulse amplitude of 10.0 volts, and pacing rate of 80 beats per minute. Pressing the " ⚠ " button

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will immediately activate the emergency pacing mode, regardless of whether the device is powered on or off. Use of the emergency pacing function should be exercised with caution. Once triggered, the device will switch to the emergency pacing parameters. To exit emergency pacing and return to other pacing settings, manually adjust the relevant parameters as needed.

Pacing Pause: By continuously pressing the "PAUSE" button, pacing pulse delivery can be suspended for up to 10 seconds. The pacing pause does not alter the sensing or pacing pulse parameters of the Lifetech Cardio Dual-Chamber Temporary Pacemaker. Please use the pacing pause function with caution, as the patient will be without pacing support during the pause.

Synchronous Pacing: Synchronous pacing refers to the operation of the Lifetech Cardio Dual-Chamber Temporary Pacemaker simultaneously with another independent pacing device. This condition can lead to confusion in the pacing and sensing functions of the Lifetech Cardio pacemaker. Synchronous pacing with the Lifetech Cardio Dual-Chamber Temporary Pacemaker is strictly prohibited.

Dual-Chamber Pacing Mode:
Dual-chamber pacing modes should

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only be selected when both atrial and ventricular channels are properly connected to the heart. If a dual-chamber mode is used while one of the channels is not connected to the heart, the open channel may pick up unintended noise. Such noise can be misinterpreted as sensed cardiac activity, potentially resulting in events including but not limited to asynchronous ventricular pacing, unintended pacing inhibition, or ventricular tachycardia.

Dual-Chamber Modes: In DDI and DDD pacing modes, the ventricular sensing amplifier may detect the atrial pacing pulse. To avoid this issue, it may be necessary to reduce the atrial pacing output, decrease the ventricular sensing sensitivity, and/or reposition the pacing leads.

Retrograde Conduction: If retrograde P-waves are sensed outside the rate-dependent automatic Post-Ventricular Atrial Refractory Period (PVARP), manually increase the PVARP until the retrograde waves fall within the refractory period. Failure to do so may result in pacemaker-mediated tachycardia (PMT).

Atrioventricular (AV) interval: Manually programming a long AV interval may cause ventricular pacing during the vulnerable period of ventricular repolarization, potentially promoting

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ventricular arrhythmias in unstable patients.

Battery status: The Lifetech Cardio dual-chamber temporary pacemaker supports at least 50 seconds of battery replacement time (no-battery operation duration). Note: During continuous device operation, please check the battery level indicator once daily to ensure sufficient power.

Battery replacement: The Lifetech Cardio dual-chamber temporary pacemaker is powered by two AA (LR6) alkaline batteries. When the low battery indicator light is on or the battery level indicator shows empty, the batteries must be replaced immediately. Please use new batteries for each patient. If the device is not in use at all, replace the batteries every two weeks to ensure sufficient power when the device is used.

Environmental Precautions: Exposure to atmospheric pressure beyond the specified range may damage the Lifetech Cardio dual-chamber temporary pacemaker. To ensure proper operation during use, transport, and storage, please observe the following:

- The device must be used in clinical settings by trained professionals.

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electrostatic discharge procedures before handling the device.

- Remove the batteries when the device is not in use.

- Do not subject the device to physical damage such as mechanical impact.

- Do not store or operate the device outside the specified temperature and humidity range (see “Device Specifications” for details).

- Dispose of used batteries in accordance with local environmental regulations.

- Do not place the pacemaker within the patient's reach.

- Do not immerse the pacemaker or patient cables in water.

- Do not disinfect the device using gamma radiation, steam, or similar methods.

- If the ambient temperature or humidity changes suddenly, allow the pacemaker to reach ambient conditions before use.

For electromagnetic compatibility information, refer to Section 6.3 “Electromagnetic Compatibility Declaration.” The device must be operated in the

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electromagnetic environment specified in Section 6.3.
Potential adverse environmental factors include, but are not limited to, the conditions listed above. In principle, the device should be operated with priority given to patient health and safety, minimizing the impact of environmental conditions on device performance.

2-3 Adverse Effects

Possible side effects associated with the use of the Lifetech Cardio Dual-Chamber Temporary Pacemaker may include:


1) Abrupt interruption before intrinsic rhythm recovery may cause arrhythmias;

2) In VOO/DOO modes, pulses overlapping with T waves may induce tachyarrhythmias.

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3 User Interface

The Lifetech Cardio Dual-Chamber Temporary Pacemaker user interface consists of the following components:

(1) Patient cable socket
(2) Indicator lights
(3) Display screen
(4) “PAUSE” button
(5) “ ⚠ ” button
(6) “FUNC” button
(7) “A-AMPL” button
(8) “RAP” button
(9) Button (◀)
(10) “SETUP” button
(11) “⏻” button
(12) “RATE” button
(13) “V-AMPL” button
(14) Button (▲)
(15) Button (▶)
(16) Button (🔒)
(17) Button (▼)
(18) Hanging ring
(19) Battery compartment key
(20) Battery compartment

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
3-1 Button

There are a total of fourteen buttons on the control panel of the Lifetech Cardio Dual-Chamber Temporary Pacemaker. The functions of each button are described below.

"⏻" Button: The"⏻" button is located at the bottom right corner of the device screen and is used to turn the device on or off. When the device is powered off, pressing the "⏻" button once will turn it on. When the device is on and needs to be turned off, ensure the keypad is unlocked, then press the "⏻" button once to display the shutdown confirmation message on the screen. While the shutdown message is still visible, press the "⏻" button again to confirm and complete the shutdown. If the message disappears before the second press, the device will not shut down; instead, the shutdown confirmation message will reappear.

"SETUP" Button: Pressing the "SETUP" button allows the user to switch to the applicable parameter configuration interfaces, including the manual parameter configuration interface, mode selection interface, and lead impedance monitoring configuration interface.


"PAUSE" Button:
Pressing the "PAUSE" button will bring up a pacing pause prompt window "Pacing Paused." Holding the "PAUSE" button for more than 1 second will stop the delivery of pacing

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pulses. The pacemaker will switch to OOO mode, which can be maintained for a maximum duration of 10 seconds.
After releasing the "PAUSE" button, the prompt window will close, and the pacemaker will automatically resume its original operating mode.


" ⚠ " Button

Pressing the " ⚠ " button activates the emergency pacing mode, in which pacing pulses are delivered in DOO mode with the following parameters: pacing rate of 80 bpm, pulse amplitude of 10.0 V, and pulse width of 0.8ms.

"FUNC" Button


Pressing the "FUNC" button switches to the applicable function interface, including the "Intracardiac Electrogram (EGM)" interface, "Pacing Ratio" interface, "Atrioventricular Conduction Time" interface, "Atrial Impedance Trend" interface, and "Ventricular Impedance Trend" interface.

"RATE" Button


Pressing the "RATE" button opens the "RATE" parameter adjustment interface. Use the "◀" and "▶" buttons to increase or decrease the value.

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"A-AMPL" Button

Pressing the "A-AMPL" button opens the "Atrial Pulse Amplitude" parameter adjustment interface. Use the "◀" and "▶" buttons to adjust the value. Press the button again to return to the main interface.

"V-AMPL" Button

Pressing the "V-AMPL" button opens the "Ventricular Pulse Amplitude" parameter adjustment interface. Use the "◀" and "▶" buttons to adjust the value. Press the button again to return to the main interface.

"RAP" Button: Press the "RAP" button to enter the Rapid Atrial Pacing (RAP) function interface. Use the "▲" and "▼" buttons to toggle between parameters, which will be highlighted in reverse color when selected. Then use the "◀" and "▶" buttons to adjust the parameter values. Long-press the "RAP" button to initiate rapid atrial pacing according to the configured parameters.

"🔓" Button:

The "🔓" button is located in the center of the directional keys and is also referred to as the Lock/Unlock button. The operator can use this button to lock all key functions except the "⚠" button. The "⚠" button cannot be locked.

Note: The screen will automatically lock if there is no operation for 1 minute.

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

Battery Compartment Key

The battery compartment keys are located on both sides of the device. The user can press both buttons simultaneously to eject the battery compartment for battery installation or removal.

3-2 Patient Cable Connector

The patient cable connector is located at the top of the Lifetech Cardio Dual-Chamber Temporary Pacemaker, through which the patient cable connects to the device.


Please connect the patient cable as follows:

1. Align the plug of the patient cable with the guide rails on the connector, then insert the cable into the socket of the Lifetech Cardio Dual-Chamber Temporary Pacemaker.

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2. Continue inserting the patient cable until a "click" sound is heard .
3. Gently pull the cable portion to confirm a secure connection; the cable should not come out.


To disconnect the patient cable, follow these steps:


1. Hold the Dual-Chamber Temporary Pacemaker firmly in place. Using two fingers, gently squeeze the locking ring on the patient cable plug and pull the cable gently in the opposite direction of the device.
2. Continue to gently pull the patient cable until it is fully disconnected from the device's connector.

Note: Do not pull the patient cable by its wire. This will not disconnect the cable and may cause damage.

3-3 Display Screen


The display screen is located at the top of the front panel of the Dual-Chamber Temporary Pacemaker. It shows the current pacing and sensing status of the device, as well as the analysis results of the pacing leads. It also provides prompt messages to the operator in the form of icons and text.


Startup Screen: When the Lifetech Cardio dual-chamber temporary pacemaker is powered on, a startup screen will be displayed for approximately 3 seconds. After the device completes its self-check, the screen will transition to the main display.

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
If any self-check errors occur during this process, an error code will be displayed, and the device will automatically shut down after ten seconds.


Main Display Screen: The main screen shows the default parameters upon startup. From top to bottom, it displays the "RATE", "A-AMPL", and "V-AMPL" parameters.


Information Bar

The information bar is located at the bottom of the screen. It displays information such as key lock status, battery level, pacing mode, emergency pacing, and more. This bar is visible across most operational screens.

Pacing Mode


The pacing mode is shown in the lower left corner of the information bar.

Lock Indicator


When the device keys are locked, a lock icon is displayed in the center of the information bar.

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Parameter Setting Indicator

When the device parameters are set to automatic mode, "AUTO" is shown in the information bar; when set to manual mode, "MANUAL" is displayed.


Battery Level Indicator

The battery level indicator is located in the lower right corner of the information bar.

Prompt Messages

When necessary, the Lifetech Cardio Dual-Chamber Temporary Pacemaker displays important text or graphic prompt messages in the center of the screen during operation. The following are built-in system prompts:

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

5. Battery Replacement Reminder

6. Two-Step Shutdown

7. Battery Cover Not Properly Closed Warning

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



"Pacing Mode Selection" Interface

Press the "SETUP" key to enter the "Pacing Mode Selection" interface. Use the "◀" and "▶" keys to scroll through the available modes. The main pacing modes include: DDD, DVI, DOO, and VVI. Press the "🔒" key to confirm and activate the selected mode.


"Lead Impedance Monitoring Configuration" Interface

Press the "SETUP" key to enter the "Lead Impedance Monitoring Configuration" interface. In this interface, you can toggle impedance monitoring on/off for atrial and ventricular leads and adjust the lower and upper impedance limits. Use the "▲" and "▼" keys to select the desired parameter, then use the "◀" and "▶" keys to modify its value.

"Rapid Atrial Pacing (RAP) Prompt" Interface

Press the "RAP" key to enter the "Rapid Atrial Pacing Prompt" interface.


"Rapid Atrial Pacing (RAP) Configuration" Interface

Press the "RAP" key to enter the "Rapid Atrial Pacing Configuration" interface. In this interface, the operator can select the Rate/Amplitude parameter using the "▲" and "▼" keys and adjust its value using the "◀" and "▶" keys. The current background pacing parameters are displayed in parentheses on the right side of the adjustable

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
parameter. Before initiating and after completing the RAP function, the Lifetech Cardio Dual-Chamber Temporary Pacemaker will deliver pulses according to the Rate and Amplitude specified in the parentheses.

"Rapid Atrial Pacing" Pop-up


When the "RAP" key is pressed and held to initiate rapid atrial pacing, a "Rapid Atrial Pacing" pop-up will appear, displaying the current Rate and Amplitude settings.

During the implementation of rapid atrial pacing, the pacing mode will be set to AOO. The operator can use the "▲", "▼" keys to select the Rate/Amplitude parameters, then press "◀", "▶"to adjust the parameters.

"Intracardiac Electrogram (EGM)" Interface

Press the "FUNC" key to switch to the "Intracardiac Electrogram (EGM)" interface. The operator can read the real-time updated atrial and ventricular EGM signals, with markers for pacing/sensing events. Press the "▲" or "▼" key to toggle between atrial or ventricular EGM; once selected, use the"◀"or "▼"key to adjust the waveform magnification.


"Pacing Percentage" Interface

This interface displays a 2D graph showing the percentage of atrial and ventricular pacing events over the past

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
five days since the last device power-on. The upper half shows the atrial pacing percentage, and the lower half shows the ventricular pacing percentage. The X-axis represents time (D1 = 1 day ago, D2 = 2 days ago, etc.). The pacing percentage is calculated as the number of pacing events divided by the total number of pacing and sensing events, updated every 3 hours.

"AV Conduction Time" Interface

This 2D graph displays the atrioventricular (AV) conduction time (ms) and the ventricular pacing percentage (%) over the past five days in dual-chamber pacing mode. AV conduction time is defined as the interval between atrial pacing/sensing and ventricular sensing. The ventricular pacing percentage is the number of ventricular pacing events divided by the total number of pacing and sensing events. The X-axis shows time (D1 = 1 day ago, D2 = 2 days ago, etc.); the left Y-axis shows AV conduction time (ms), and the right Y-axis shows ventricular pacing percentage (%).


"Atrial Impedance Trend" Interface

This 2D graph displays the atrial impedance trend over the past five days since the last power-on. The X-axis represents time (D1 = 1 day ago, D2 = 2 days ago, etc.), and the Y-axis

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shows the impedance trend during the five-day period.

"Ventricular Impedance Trend" Interface
This 2D graph displays the ventricular impedance trend over the past five days since the last power-on. The X-axis represents time (D1 = 1 day ago, D2 = 2 days ago, etc.), and the Y-axis shows the impedance trend during the five-day period.

3-4 Indicator Lights

There are five indicator lights arranged in a row at the top front panel of the Lifetech Cardio Dual-Chamber Temporary Pacemaker.

"Atrial Pacing" Indicator Light
The "Atrial Pacing" indicator light is green. When the device is operating normally, this indicator flashes in sync with atrial pacing events, that is, each time the device delivers an atrial pacing pulse, the indicator flashes once.

"Atrial Sensing" Indicator Light
The "Atrial Sensing" indicator light is white. When the device is operating normally, this indicator flashes in sync with atrial sensing events, that is, each time the device detects an atrial activity outside of the refractory period, the indicator flashes once.

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"Ventricular Pacing" Indicator Light
The "Ventricular Pacing" indicator light is green. When the device is operating normally, this indicator flashes in sync with ventricular pacing events, that is, each time the device delivers a ventricular pacing pulse, the indicator flashes once.

"Ventricular Sensing" Indicator Light
The "Ventricular Sensing" indicator light is white. When the device is operating normally, this indicator flashes in sync with ventricular sensing events, that is, each time the device detects a ventricular activity outside the refractory period, the indicator flashes once.

"Power" Indicator Light
The "Power" indicator light is orange and located at the top right of the device. The Lifetech Cardio Dual-Chamber Temporary Pacemaker continuously monitors battery power, and this indicator will flash when the battery is low. When the "Power" indicator flashes, immediately check the power supply and replace both AA (LR6) alkaline batteries with new ones. After replacement, verify that the "Power" indicator has stopped flashing.

3-5 Low Battery Warning

The Lifetech Cardio Dual-Chamber Temporary Pacemaker continuously

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monitors battery power. When the battery level becomes low, the orange "Power" indicator starts flashing and the buzzer sounds. Once the low battery warning first appears, the device can typically continue operating for approximately 24 hours.

3-6 Battery Compartment

The battery compartment is located at the bottom of the device and holds two AA alkaline (LR6) batteries.

Please follow the steps below to open the battery compartment:
Press the battery compartment buttons on both sides simultaneously to open the battery compartment.

Please follow the steps below to close the battery compartment:
After inserting the batteries into the compartment, push it back into the device until a "click" sound is heard, indicating that the compartment is securely in place.

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

The hanging ring is located at the upper part of the back of the device. Operators can use this ring to suspend the device, such as hanging it on an IV pole or bed rail. To use it, pull the ring outward from the slot and rotate it 180 degrees until it gently clicks into place. When not in use, please return the hanging ring to its original position.

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4 Device Operation

4-1 Preparation

The Lifetech Cardio Dual-Chamber Temporary Pacemaker should only be operated by trained professionals in a clinical setting. Use the supplied patient cable with adapters as needed to connect temporary pacing leads. The pacing lead system comprises the pacemaker, patient cable, adapter/extension, temporary pacing leads, and the patient's myocardium. Follow the steps below and ensure secure connections before use.

1. Connect the patient cable to the device per Section 3.2 "Patient Cable Socket."

2. Connect the electrode lead (with adapter/extension if used) to the patient cable's electrode connector.

1) Turn the connector's clamping screw counterclockwise until the lead pin can be inserted.

2) Fully insert the lead pin into the connector until it cannot go further.

3) Turn the clamping screw clockwise to secure the pin. Gently pull the lead to confirm it does not come loose.

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4) Observe chamber and polarity when connecting. The blue “A” connector is for the atrium; the gray “V” connector is for the ventricle. Polarity (+/−) is labeled on the connectors. When using a bipolar pacing lead, connect the adapter/lead wire corresponding to the lead tip (tip electrode) to the negative (-) terminal, and connect the adapter/lead wire corresponding to the ring electrode to the positive (+) terminal. When using a unipolar pacing lead, connect the lead coming from the heart to the negative (-) terminal, and connect the lead coming from the body surface to the positive (+) terminal.

3. Recheck all connections between the device, patient cable, and leads to ensure secure and safe connections.

To disconnect the cable, follow Section 3.2 instructions.

4-2 Power On

Press the power button or the “ ⚠ ” button to power on. A startup screen appears for 3 seconds.

Note: Do not operate outside the specified environmental conditions. If the device undergoes a sudden temperature change, allow it to reach ambient temperature before use.
Note: Do not press any buttons during startup to avoid false key-fault detection.

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If the device detects an error during startup, it will delay the power-on process. In certain cases, the operator may see an error code (refer to Chapter 7 for how to handle error codes). Additionally, the device performs continuous real-time self-checks during operation to ensure proper functionality.

4-3 **Power-Off**

After shutdown, all settings revert to power-on defaults. A two-step confirmation prevents accidental shutdown.

Please follow the steps below to power off the device:

1. Confirm that the buttons are unlocked. If the buttons are locked, press the “🔓” button to unlock them.
2. Press the “⏻” key once to display the shutdown prompt.
3. Press the “⏻” key again before the prompt disappears to turn off the device.

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In any case, if the operator wishes to exit the shutdown prompt screen and continue device operation, either of the following two methods can be used:

1. Take no action until the prompt disappears.
2. Press any key other than the “⏻” key and the “ ⚠ ” key to exit the shutdown prom.

4-4 Lock/Unlock

The Lifetech Cardio Dual-Chamber Temporary Pacemaker is equipped with fourteen keys. Except in the "Intracardiac Electrogram (EGM)" interface, the device will automatically lock the keys if there is no user operation for 1 minute. In the "EGM" interface, the keys will automatically lock after 5 minutes of inactivity. Use the center round button to lock/unlock. The “🔒” key cannot be locked.

Keypad and Lock Indicator

When the keys are locked, a lock icon will appear in the information bar at the bottom of the display. While the keypad is locked, pressing any key (except the “ ⚠ ” key) will bring up a lock message in the center of the screen indicating that the keys are locked. To unlock, press the “🔒” key.

During operation, the “⏻”, “SETUP”, “PAUSE”, “FUNC”, “RAP”, “RATE”, “A-AMPL”, “V-AMPL”, “▲” , “▼”, “◀”, and “▶” keys can all be locked using the “🔒” key. The “ ⚠ ” key cannot be locked.

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The “🔓” key is used to toggle between locked and unlocked states.

4-5 Pulse Amplitude

Atrial and ventricular amplitudes range 0.1 – 20.0V and can be set independently.

Please refer to the descriptions of the “A-AMPL” key and “V-AMPL” key in Section 3.1 to set the atrial pulse amplitude and ventricular pulse amplitude.

4-6 Pulse Width


Pacing pulse width is fixed at 0.8ms.

4-7 Pacing Rate

Base rate 30–200 ppm; atrial and ventricular rates are identical and not set separately.


I Adjust via the “RATE” key (see Section 3.1)

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
The atrial sensing sensitivity and ventricular sensing sensitivity both range from 0.5 mV to 20 mV. The atrial and ventricular sensing sensitivities can be set independently to different values.


Please refer to the description of the "SETUP" key in Section 3.1 and the "Manual Parameter Configuration" interface in Section 3.3 of the display screen to set the atrial sensing sensitivity and ventricular sensing sensitivity.

4-9 Atrioventricular (AV) Interval

Unless manually set, the atrioventricular (AV) interval parameter is automatically adjusted according to the basic rate. When set manually, the paced AV interval (PAV) can be adjusted within the range of 20 to 300 milliseconds.


Please refer to the "Set" button in "3.1 Buttons" and the manual parameter configuration interface in "3.3 Display" to set the AV interval. When manually setting the AV interval parameters, "AUTO" in the information prompt bar on the main interface will change to "MANUAL".

In both automatic and manual settings, the SAV (sensed AV interval) varies with the PAV (paced AV interval) within a defined range. When PAV is between 80 ms and 280 ms, SAV = PAV

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
- 30 ms; when PAV is less than 80 ms, SAV = 50 ms; and when PAV is greater than 280 ms, SAV = 250 ms.


4-10 Post-Ventricular Atrial Refractory Period (PVARP)

Unless manually configured, the Post-Ventricular Atrial Refractory Period (PVARP) is automatically adjusted according to the basic pacing rate. When manually configured, the PVARP can be set within a range of 150 to 500 milliseconds.

Please refer to the description of the "SET" key in section 3.1 "Keys" and the "Manual Parameter Configuration" interface in section 3.3 "Display" to set the Post-Ventricular Atrial Refractory Period (PVARP). When the PVARP parameter is set manually, the "AUTO" in the information prompt bar on the main screen will change to "MANUAL."


Note: When manually setting the atrioventricular delay (AV interval) and the Post-Ventricular Atrial Refractory Period (PVARP), the device may display "Prompt Message 8" as shown in section 3.3 "Display." If this message appears, please readjust the AV interval and PVARP. The manually adjusted parameters must meet the following

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
4-11 Atrial Tracking (A-TRG)


Atrial Tracking (A TRG) is only accessible and applicable when the dual-chamber temporary pacemaker is configured for dual-chamber sensing and pacing. When applicable, A TRG is enabled by default.

Please refer to the description of the "SET" key in section 3.1 and the "Manual Parameter Configuration" interface in section 3.3 "Display" to configure the Atrial Tracking (A TRG) ON/OFF settings.

4-12 Pacing Mode Selection

Press the "SET" key to enter the "Pacing Mode Selection" interface. Use the "◀" and "▶" keys to select the desired pacing mode. Available main pacing modes include: DDD, DVI, DOO, VVI. Press the "🔓" key to confirm the selected pacing mode, which will then take effect, and a pop-up message will appear for confirmation.

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4-13 Noise Reversion

When continuous noise with cycle ≤ 40 ms (e.g., 50/60 Hz) is sensed, the device reverts to asynchronous pacing at the preset rate until the noise stops.

4-14 Rate Limiting (Runaway Protection)


The runaway protection function is designed to limit excessive pacing output and ensure patient safety. The Lifetech Cardio Dual-Chamber Temporary Pacemaker incorporates three independent protection mechanisms that

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cross-check each other. Under this mechanism, the pacing rate will not exceed 230 beats per minute.

When the device first detects runaway output, it will automatically restore the default parameters. If the runaway condition persists, upon a second detection, the device will display an error code for 10 seconds and then shut down.

4-15 Emergency Pacing


Pressing the " ⚠ " key puts the device into emergency pacing mode; this key is not locked by the " 🔒 " key and has a physical anti-misoperation design.

In emergency pacing mode, the Lifetech Cardio Dual-Chamber Temporary Pacemaker operates with the following parameters: DOO mode, pulse width 0.8ms, pulse amplitude 10.0 V, pacing rate 80 bpm, and sensing sensitivity set to SAY.

Regardless of the device's power state, pressing the " ⚠ " key will directly activate emergency pacing mode. Use this function with caution. Once triggered, relevant parameters are set to emergency pacing defaults. To change any parameters under emergency pacing mode, please adjust them individually.

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4-16 Pacing Pause

The operator can hold down the "PAUSE" key to suspend pacing pulses for no more than 10 seconds. The pause function is effective only while the key is being held. This function facilitates efficient observation of the patient's intrinsic rhythm. It does not alter the sensing or pacing pulse parameters of the Lifetech Cardio Dual-Chamber Temporary Pacemaker.

When the operator presses this key, the display will show a pacing pause message. The pacing mode will change to OOO mode, during which all other buttons are disabled.

To pause pacing, follow these steps:
1. Confirm that the button is unlocked.
2. Press and hold the Pause button to pause pacing pulses.
3. Release the button to end the Pause.

4-17 Ventricular Safety Pacing

This function is intended to prevent inappropriate inhibition of ventricular pacing pulse delivery when non-ventricular events (such as noise or interference) are sensed on the ventricular lead. It operates in dual-chamber and

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synchronized pacing modes (DVI, DDI, and DDD).

If a ventricular sense occurs outside the 30 ms blanking period but within 105 ms after atrial pacing: if A-V > 105 ms, VSP is delivered at 105 ms after atrial pacing; if A-V < 105 ms, the ventricular pace is delivered within the A-V interval.


4-18 Wenckebach Response

In DDD mode, if the intrinsic atrial rate exceeds the maximum tracking rate, the SAV is prolonged until the P-wave falls within the PVARP, preventing sensing and subsequent ventricular pacing.

4-19 Mode Switch

Upon detection of atrial tachycardia, the device temporarily switches from DDD pacing mode to non-tracking DDI pacing mode to prevent rapid ventricular pacing caused by high atrial rate, and automatically returns to the previously set pacing mode once the atrial tachycardia terminates.

4-20 Intracardiac Electrogram (EGM)


By pressing the “FUNC” key, the oper ator can switch to the "Intracardiac Electrogram (EGM)" interface to view real-time atrial and

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ventricular EGM segments no longer than 3 seconds . The waveform will also display pacing and sensing event markers.


Please refer to the “FUNC” key described in section “3.1 Keys” and the “EGM” interface described in section “3.3 Display” to access this function interface. In the “EGM” interface, use the “▲” “▼” keys to select the atrial or ventricular EGM, and use the “◀” “▶” keys to adjust the zoom level of the electrogram.

4-21 Pacing Percentage

Pressing the “FUNC” key allows switching to the "Pacing Percentage" interface. This interface displays a 2D chart showing the pacing percentage (%) recorded by the Lifetech Cardio Dual-Chamber Temporary Pacemaker over the past five days since the last power-on. The horizontal axis represents time, with D1 indicating 1 day ago, D2 indicating 2 days ago, and so on. Values reset upon device shutdown.

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
The operator can switch to the "AV Conduction Time" interface by pressing the “FUNC” key. This interface displays a 2D chart showing the atrioventricular (AV) conduc-

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
tion time (ms) and ventricular pacing percentage (%) under dual-chamber pacing mode, recorded over the past five days since the last power-on of the Lifetech Cardio dual-chamber temporary pacemaker. The horizontal axis represents time, with D1 indicating 1 day ago, D2 indicating 2 days ago, and so on. The left vertical axis shows the AV conduction time (ms) over the five days, while the right vertical axis displays the ventricular pacing percentage (%) over the same period.

4-23 Lead Impedance Monitoring


Lead Impedance Monitoring Configuration
The operator can switch to the "Lead Impedance Monitoring Configuration" interface by pressing the "SETUP" key. Atrial and ventricular lead monitoring functions are turned off (OFF) by default when the device is powered on. To enable lead monitoring, the operator can press the "▲" and "▼" keys to select atrial or ventricular lead impedance monitoring in the "Lead Impedance Monitoring Configuration" interface. Then press the "◀" or "▶" key to turn ON the atrial/ventricular lead impedance monitoring function.


Note: Impedance monitoring can only been abled when the voltage (pulse amplitude) is ≥ 5.0 V. If the voltage requirement is not met and the user attempts to enable the lead impedance monitoring function, Lifetech Cardio Dual-Chamber Temporary Pacemaker will display a pop-up warning in the "Lead Impedance Monitoring Configuration" interface. Refer to "Prompt Message" 9 in Section 3.3. The operator may press

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the "A-AMPL" or "V-AMPL" key to switch to the main screen and adjust the corresponding pulse amplitude parameter.


Upon powering on, the device sets the default low impedance alert threshold to 200 Ω and the high impedance alert threshold to 2000 Ω. To modify the lower or upper threshold for lead impedance monitoring, the operator can use the "▲" and "▼" keys in the "Lead Impedance Monitoring Configuration" interface to select the low or high impedance threshold, respectively. and then press the "◀" or "▶" key to adjust the threshold value.

The available parameter options for lead imped ance monitoring are listed below.

Optional parameter	Default	Options range
Atrial lead monitoring	OFF	OFF/ON
Ventricular lead monitoring	OFF	OFF/ON
Low-impedance lower limit (Ω)	200	100/200/300 /OFF
High-impedance upper limit (Ω)	2000	1000/2000/3 000/OFF


Lead Impedance Monitoring Alerts


When the Lifetech Cardio Dual-Chamber Temporary Pacemaker detects low or high impedance, the detected abnormal lead impedance status will be displayed in the message bar of the screen. When applicable, this includes alerts such as "Atrial Lead High Impedance" or "Atrial Lead Low Impedance," or "Ventricular Lead High Impedance" or "Ventricular Lead Low Impedance."

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4-24 Lead Impedance Trends


When the atrial and ventricular lead impedance monitoring functions are turned ON, the user can press the “FUNC” key to switch to the "Atrial Impedance Trend" screen. This screen displays a 2D graph showing the trend of atrial lead impedance over the past five days since the most recent device startup. The horizontal axis represents time, where D1 indicates 1 day ago, D2 indicates 2 days ago, and so on.

The vertical axis displays the impedance values of the atrial lead over the five-day period. When the ventricular lead impedance monitoring function is turned ON, the user can press the “FUNC” key to switch to the "Ventricular Impedance Trend" screen. This screen displays a 2D graph showing the trend of ventricular lead impedance over the past five days since the most recent device startup. The horizontal axis represents time, where D1 indicates 1 day ago, D2 indicates 2 days ago, and so on. The vertical axis displays the impedance values of the ventricular lead over the five-day period.


4-25 Rapid Atrial Pacing (RAP)


Pressing the “RAP” key switches to the interface related to the “Rapid Atrial Pacing” (RAP) function, where the bottom left corner displays the pacing mode and the bottom right corner shows the battery level.

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Before entering the "Rapid Atrial Pacing Con figuration" screen, the device will first display a "RAP Prompt" screen, indicating that only the atrium supports high-frequency pacing, and prompting the user to ensure that the Lifetech Cardio dual-chamber temporary pacemaker is connected to an atrial lead. Pressing the "RAP" key again confirms the prompt and allows access to the "Rapid Atrial Pacing Configuration" screen.

Upon entering the "Rapid Atrial Pacing Configuration" screen, the Lifetech Cardio Dual-Chamber Temporary Pacemaker continues to operate according to the previously set mode and parameters. Within the "Rapid Atrial Pacing Configuration" screen, the Rate and Amplitude parameters for RAP can be adjusted. The operator may use the "▲" and "▼" keys to select the Rate or Amplitude parameter, and then use the "◀" and "▶" keys to adjust the parameter value. The current background pacing parameters are displayed in parentheses to the right of each adjustable parameter, and both before initiating and after terminating the RAP, the pacemaker will deliver pulses using the Rate and Amplitude values shown in the parentheses.

After completing the parameter adjustments, pressing and holding the "RAP" key will initiate Rapid Atrial Pacing using the Lifetech Cardio Dual-Chamber Temporary Pacemaker, and a "Rapid Atrial Pacing" pop-up window will appear, indicating that the pacing mode is AOO. During the RAP process, the operator may use the "▲" and "▼" keys to select the Rate or Amplitude parameter, and use the "◀" and "▶" keys to adjust the values. Releasing the "RAP" key will terminate the Rapid Atrial Pacing. Holding down the "RAP" key continuously allows a maximum of 120 seconds of RAP, to

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continue RAP beyond this duration, please press and hold the “RAP” key again.
If the Lifetech Cardio Dual-Chamber Temporary Pacemaker continues to deliver Rapid Atrial Pacing after the “RAP” key has been released, power off the device using the “Power” key to terminate the Rapid Atrial Pacing. If the pacemaker continues Rapid Atrial Pacing despite powering off, remove the battery from the pacemaker and return the device for maintenance.

To exit the "Rapid Atrial Pacing" function interface, press the “RAP” key to return to the main interface.

4-26 Battery Status

Under the DDD pacing mode with a base rate of 80 bpm, pulse width of 0.8 ms, pulse amplitude of 5.0 V, and sensing sensitivity (atrium: 0.5 mV, ventricle: 2.0 mV), two new AA (LR6) alkaline batteries can power the Lifetech Cardio dual-chamber temporary pacemaker for 14 days.
During operation, the device continuously monitors battery status and alerts the user to replace the batteries when the voltage is low.
When the battery is low, the Lifetech Cardio Dual-Chamber Temporary Pacemaker notifies the user via a flashing "Power" indicator light and audible beeping.
When either indication appears, the user must replace the batteries immediately.

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1. "Power" indicator light: Located at the top right of the device. This indicator flashes to alert the user to replace the battery.
2. Audible alert: The Lifetech Cardio Dual-Chamber Temporary Pacemaker emits continuous "beep-beep" sounds when the battery is nearly depleted, prompting the user to replace the battery immediately to ensure proper operation.

Note: After the low battery warning appears, the electrical parameters listed in section "5.1 Specifications" remain within their tolerance range under a load of 500Ω±5Ω and an ambient temperature of 20℃±2℃.

4-27 Replacing Batteries

Please use a new battery for each new patient pacing therapy.

1. Replace the battery at least every two weeks with continued device use.
2. Replace the battery promptly whenever the low-battery indicator appears.
3. During battery replacement, the Dual-Chamber Temporary Pacemaker can operate continuously without batteries for at least 50 seconds.
4. Remove the battery when the device is not in use.

Warning:
1. Avoid replacing the battery while the device is connected to a patient unless necessary.
2. Do not use rechargeable batteries.

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3. Using batteries that do not meet the recommended model or size, or batteries with contaminants near the electrodes, may result in unstable or no device output.
The battery compartment is located on the bottom of the Dual-Chamber Temporary Pacemaker and can accommodate two AA alkaline (LR6) batteries. There are battery compartment buttons on both sides of the device. Press both buttons simultaneously to open the battery compartment.
Note: When the battery compartment is opened, the device battery is electrically isolated and cannot power the device properly. The battery compartment must be closed for normal use.

Battery Compartment Closure Warning

The battery compartment of the Dual-Chamber Temporary Pacemaker should be completely closed when in operation. If the battery compartment is not completely closed, the device's "Power" indicator light flashes and the buzzer continues to sound until the battery compartment is completely closed.

4-28 Cleaning and Disinfection

Thoroughly wipe the Lifetech Cardio dual-chamber temporary pacemaker and patient cables using a damp cloth or sponge moistened with clean water.
Disinfect the device and patient cables using a sponge or cloth moistened with 75% ethanol.
After each cleaning/disinfection, have a qualified technician verify proper function. Ensure

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the device and cable are completely dry—especially connectors and the socket—before use.

Note: Soaking the Dual-Chamber Temporary Pacemaker in water or disinfectant, or using gamma rays or high-temperature steam to disinfect the device may damage the device and cause it to malfunction.

4-29 Post-Cleaning/Disinfection Checks

After each cleaning/disinfection, perform the following inspections on the patient cables and the Lifetech Cardio dual-chamber temporary pacemaker.

- Inspect the device and patient cable to ensure they are not damaged and free of surface contamination.
- Verify that the device and patient cable are completely dry, especially the electrode connector and the inside of the patient cable receptacle.
- Verify that the screen display and all buttons and knobs are functioning properly.
- Verify that the patient cable is properly connected to the device and that the device self-test displays no errors.

If one or more of the above inspections fail, please contact Lifetech Cardio after-sales service. For contact information, refer to the front page.

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5 Device Specifications

5-1 Parameter

Unless otherwise specified, all parameters of the Lifetech Cardio Dual-Chamber Temporary Pacemaker refer to data measured at 20 °C ± 2 °C with a 500-ohm (±1%) load.

Service life: 6 years.

The device contains no irreplaceable parts and can be returned to the factory for repair in the event of any failure. Lifetech Cardio guarantees to provide replacement parts for at least ten years after delivery.

Device Dimensions

The physical dimensions of the Lifetech Cardio Dual-Chamber Temporary Pacemaker are as follows (with the suspension ring retracted):

	Dimensions (mm)	Tolerance (mm)
Length	190	±10
Width	86	±8
Thickness	37	±5

Patient Cable Dimensions

The dimensions of the Patient Cable accessory are as follows:

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	Dimensions (mm)	Tolerance (mm)
Length	1860	±50
diameter of compatible adapter connectors	>3	NA

Storage Conditions (without batteries)

Temperature °C	Atmospheric Pressure (hPa)	Humidity (%)
Min: -20 Max: 50	700-1060	< 90

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Battery

The Lifetech Cardio Dual-Chamber Temporary Pacemaker uses two AA (LR6) alkaline batteries.

1) At 20 °C, the device can operate continuously for at least 10 days on two new AA (LR6) alkaline batteries.

2) Under the same temperature, using default pacing/sensing parameters (80 bpm, 0.8ms width, 5.0V amplitude, atrial sensitivity 0.5 mV, ventricular sensitivity 2.0 mV), the device can typically last 14 days.

3) Battery life may vary by brand; AA (LR6) alkaline batteries such as Nanfu are recommended.

Battery Type	Battery Life(days)
AA (LR6) alkaline battery	Typically:14; Mnimum:10

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The backup power ensures at least 50 seconds of continuous operation.

Software Release Version

The software release version of the Lifetech Cardio Dual-Chamber Temporary Pacemaker is Version 1.

Pacing Modes

Available pacing modes include: DDD, DOO, DDI, DAI, DVI, AOO, AAI, VOO, and VVI.

Chamber

Lifetech Cardio Dual-Chamber Temporary Pacemaker can pace the atria or ventricles independently, or perform dual-chamber pacing of the atria and ventricles.

Run-Away Protection

The runaway protection trigger frequency is greater than 230 ppm.

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Note: When PAV is between 80-280ms, SAV = PAV – 30ms; if PAV < 80ms, SAV = 50ms; if PAV > 280ms, SAV = 250 ms.

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Rate

	Range (ppm)	Toler-ance (%)	Step Size (ppm)
Basic Rate	Min: 30 Max: 200	+/- 5	5, from 30 to 60 4, from 60 to 120 10, from 120 to 200
RAP Rate	Min: 80 Max: 1000	+/- 5	30, from 80 to 200 5, from 200 to 250 10, from 250 to 360 40, from 360 to 800 50,>800

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Pulse width: 0.8ms±10% (fixed)

Maximum Tracking Frequency

Range (ppm)	Tol-er-ance (%)	Step (ppm)
Min: 110 Max:230	+/- 5	Changes with pulse fre-quency

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		or 30%, which-ever is greater, from 0.5 to 1.0mV	from 0.5 to 1.0 0.5, From 1.0 to 3.0 1.0, From 3.0 to 20.0
		±20%, from 1.5 to 20.0mV	

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Medical Electrical Equipment – Part 1: General Requirements for basic safety and essential performance.

EMC Standard: YY 9706.102-2021 Medical Electrical Equipment – Part 1-2: General requirements for basic safety and essential performance- Collateral Standard: Electromagnetic Compatibility – Requirements and Tests;

Product Standard: YY 9706.231-2023 Medical Electrical Equipment – Part 2-31: Particular Requirements for the Safety and essential performance of external cardiac pacemakers with internal power source.

Essential performance of the device: Based on GB 9706.1-2020 and YY 9706.231-2023, the essential performance of this product is as follows:

1. Battery depletion indicator

2. ME equipment parameter stability

3. Pulse amplitude stability

4. Disabling runaway rate protection

5. Deliberate action required to change settings

6. Parameter stability at onset of the battery depletion indicator

7. Runaway protection

8. Interference reversion in the presence of sensed electromagnetic disturbance or electrical energy sources

9.Limit at which the ventricle is paced in response to sensed atrial activity

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Safety Features of the Device:

Classification by electric shock protection type: Internal power supply

Classification by degree of protection against electric shock: CF type application

Classification by degree of protection against liquid ingress: IPX0

Classification by safety when used with flammable anesthetic gases mixed with air or with flammable anesthetic gases mixed with oxygen or nitrous oxide: Not applicable

Power supply mode: Two AA 1.5V alkaline batteries (LR6)

Application section with protection against defibrillation discharge: CF type application

Signal output or input section: EGM input, pacing signal output, and data output

Non-permanent installation

Pollution degree: 2

Overvoltage category: I

Altitude: ⩽4000m

Other: Handheld device

5-3 EMC Declaration

Special EMC precautions must be taken. The device should be installed and used according to the EMC information provided.

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The device complies with YY 9706.102-2021 and is classified as Group 1, Class A according to GB 4824-2019.

Medical electrical equipment requires special EMC precautions and should be used in specified environments.

Portable and mobile RF communications equipment may affect device performance.

Use of accessories or cables other than those sold as internal components may increase emissions or decrease immunity of the Lifetech Cardio Dual-Chamber Temporary Pacemaker.

This device includes the following cable: Patient cable, length 1.86 m, shielded and dedicated.

Warning: Use of unauthorized accessories, transducers, or cables may increase emissions or decrease immunity of the device.

The device should not be used adjacent to or stacked with other equipment. If unavoidable, verify proper operation in the configuration.

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The user should ensure it is used in such an environment:		
Emission test	**Compliance**	**Electromagnetic environment -**

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		guide
RF Emissions GB 4824	Group 1	The device uses RF energy only for internal functions. Its RF emissions are low and unlikely to interfere with nearby electronic equipment
RF Emissions GB 4824	Class A	The device is suitable for use in all establishments except those directly connected to the public low-voltage power supply network.
Harmonic emissions GB 17625.1	Not applicable	
Voltage fluctuations/flicker emissions GB 17625.2	Not applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
The Lifetech Cardio Dual-Chamber Temporary Pacemaker is intended for use in the electromagnetic environment specified below. The purchaser or user should ensure that it is used in such an environment.			
Immunity test	**IEC 60601 test level**	**Compliance level**	**Assessment**
Electrostatic discharge GB/T 17626.2	±6 kV contact discharge ± 8 kV air discharge	±8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete, or ceramic tile. If covered with synthetic material, relative humidity should be at least 30%.
Electrical Fast Transient/Burst GB/T 17626.4	Not Applicable	Not Applicable	/
Surge GB/T 17626.5	Not Applicable	Not Applicable	/
Voltage Dips, Short	Not Applicable	Not Applicable	/

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Interruptions and Voltage Variations on Power Supply Input Lines GB/T 17626.11			
Power Frequency Magnetic Field (50 Hz/60 Hz) GB/T 17626.8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF GB/T 17626.6 Radiated RF GB/T 17626.3	3 V (rms) 150 kHz to 80 MHz (excluding ISM bands) 10 V (rms) 150 kHz to 80 MHz (for ISM bands) 10 V/m 80 MHz to 2.5 GHz	3 V (rms) 10 V (rms) 10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the 6108 Temporary Pacemaker, including cables, than the recommended separation distance calculated from the applicable formula. Recommended separation distance: d = [1.17]√‾Pd = [1.20]√‾Pd = [1.20],√‾P (80 MHz – 800 MHz)d = [2.30]

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			√P (800 MHz - 2.5 GHz) Where: P = max rated output power of the transmitter in watts (W)d = recommended separation distance in meters (m) The field strengths of fixed RF transmitters should be determined by an electromagnetic site survey and should be less than the applicable compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol. 
Note 1: At 80 MHz and 800 MHz, the separation distance formula for the higher frequency range applies. Note 2: These guidelines may not apply in all cases. Electromagnetic propagation is			

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affected by absorption and reflection from structures, objects and people.

a) The ISM frequency bands between 150 kHz and 80 MHz refer to: 6.765–6.795 MHz, 13.553–13.567 MHz, 26.957–27.283 MHz, and 40.66–40.70 MHz

b) For ISM bands from 150 kHz to 80 MHz and frequencies from 80 MHz to 2.5 GHz, compliance levels are designed to reduce the likelihood of interference from portable/mobile RF communications equipment accidentally brought into the patient area. A 10/3 factor is applied to the separation distance in these cases.

c) Fixed transmitters, such as base stations for radio (cellular/cordless), land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast, may not be predictable. A site survey should be considered to assess the electromagnetic environment.

If the measured field strength at the location exceeds the applicable RF compliance level, the device should be observed to verify normal operation.

If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device..

d) Across the 150 kHz–80 MHz frequency range, field strength should be less than 3 V/m.

Recommended separation distances between portable/mobile RF communications equipment and the Lifetech Cardio Dual-Chamber Temporary Pacemaker

The pacemaker is intended for use in an RF-controlled electromagnetic environment. The user can help prevent electromagnetic interference by maintaining the recommended minimum distance between the device and any RF transmitter according to its output power.

Maximum Rated Output Power of Transmitter (W)	Corresponding Separation Distances by Frequency (m)			
	150 kHz～80 MHz (excluding ISM bands) d=[1.17]√P	150 kHz～80 MHz (excluding ISM bands) d=[1.20] √P	80 MHz～800 MHz d=[1.20] √P	800 MHz～2.5 GHz d=[2.30] √P
0.01	0.12	0.12	0.12	0.23
0.1	0.37	0.38	0.38	0.73
1	1.17	1.20	1.20	2.30

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For transmitters not listed above, the recommended separation distance d (in meters) can be calculated using the appropriate formula, where P is the transmitter's maximum rated output power in watts (W).

Note 1: At 80 MHz and 800 MHz, use the formula for the higher frequency band.

Note 2: The ISM bands from 150 kHz to 80 MHz include: 6.765–6.795 MHz, 13.553–13.567 MHz, 26.957–27.283 MHz, and 40.66–40.70 MHz.

Note 3: A 10/3 factor is used in the calculation of separation distance in ISM bands (150 kHz–80 MHz) and 80 MHz–2.5 GHz, to reduce the risk of interference in patient areas.

Note 4: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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5-4 Product Models

The Lifetech Cardio Dual-Chamber Temporary Pacemaker models 9102, 9101, and 9100 differ in configuration as specified in section 1.3.
All three models share the same pacing functionality and safety features; configuration differences are shown in the table below.

Model / Other Features	9102	9101	9100
Intracardiac Electrogram (EGM)	√	√	——
Pacing Percentage	√	——	√
AV Conduction Time	√	——	——
Lead Impedance Monitoring Alerts	√	√	——
Lead Impedance Trend	√	√	——
Rapid Atrial Pacing (RAP)	√	√	√
Auto-Adjustment Function	√	√	√

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6 Support, Disclaimer, and Warranty

6-1 Support

Lifetech Cardio provides professional technical and after-sales support. Contact information is available on the front page.
Each device has a serial number on its back. Please provide this number when requesting support.

Troubleshooting
In the event of malfunction or accidental damage (e.g., fluid intrusion or dropping), the operator will not be harmed. Replace with another unit immediately for temporary pacing.
Power off the faulty device, restore its original condition if possible, wait 10 seconds and restart. If still unresponsive, reinstall the batteries and restart.

If issues persist, contact Lifetech Cardio after-sales support (see front page).

Device Repair
Do not attempt unauthorized repairs. Contact Lifetech Cardio for service.
Professional Training
The after-sales team offers training for operators who require it.

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6-2 Disclaimer

The device is intended for temporary pacing in cases of arrhythmia or bradycardia. Lifetech Cardio does not guarantee effectiveness in all patients.

Patient cable

All cables are subject to fatigue and damage due to repeated use.
Improper handling and repeated pressure may accelerate fatigue. Lifetech Cardio does not warrant the cable's durability under long-term, repeated use.
Despite careful design and testing, cables may still fail due to unexpected factors.
Warnings on product labeling provide further detail and form part of the disclaimer.
Lifetech Cardio assumes no liability for damage resulting from use, defect, or failure of patient cables.
The above disclaimers shall not be deemed to contradict applicable mandatory legal provisions.
If any part or provision of this warranty statement is found to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the remainder of the warranty statement shall not be affected, and all rights and obligations shall remain recognized and enforceable as if the invalidated portion had not been included.

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6-3 Warranty Information

Overview

In the absence of man-made damage, Lifetech Cardio provides the operator with a one-year warranty for the Dual-Chamber Temporary Pacemaker, starting from the date of delivery.

If the pacemaker experiences non-man-made damage during normal use, please contact Lifetech Cardio after-sales service. Contact information is provided on the front page.

Lifetech Cardio reserves the right to choose one of the following two methods of repair:

1) Repair or replace any defective component of the device;

2) Provide a replacement device with equivalent functionality.

Limitations

This warranty statement is valid only under the following conditions:

1. The device has not been tampered with, disassembled, or repaired by anyone other than Lifetech Cardio.

2. The failure is not caused by improper operation, mishandling, or accident by the user.

3. Lifetech Cardio after-sales service is contacted within 60 days of the failure.

Disclaimer

This warranty statement is limited to the express terms set forth herein and

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does not include any other express or implied warranties. Lifetech Cardio shall not be liable for any incidental or consequential events arising from or related to the use, defect, or failure of the device, regardless of whether such claims are based on warranty, contract, or tort.

The purpose of the above disclaimers and limitations is not to conflict with any mandatory provisions of applicable law.

If any part or provision of this warranty statement is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the remainder of the statement shall remain in full force and effect, and all rights and obligations shall continue to be recognized and performed.

As if the invalidated provision had never been included in the warranty statement.

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