iChemVELOCITY Wash Solution

REF 800-7217

Rx Only in the U.S.A.

INTENDED USE

The iChemVELOCITY Wash Solution is a buffer solution for routine rinse of the iChemVELOCITY or DxU 810c Series of Automated Urine Chemistry Analyzers. This reagent is intended for *in vitro* diagnostic use with the iChemVELOCITY or DxU 810c Series of Automated Urine Chemistry Analyzers.

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PRINCIPLE

The iChemVELOCITY Wash Solution is a ready-for-use solution specifically developed for the application in the iChemVELOCITY or DxU 810c Series of Automated Urine Chemistry Analyzers. Use of this Wash Solution properly provides the optimal operating condition for the System.

WARNING AND PRECAUTIONS

- The iChemVELOCITY Wash Solution has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Product is colorless, clear liquid with aromatic odor. Product that differs from these attributes should be discarded as they may have deteriorated.
- Avoid contact of solutions with eyes and skin. Refer to the Safety Data Sheet for first aid measures.
- Wear appropriate personal protective equipment when handling.
- Avoid long term exposure to skin or by inhalation. Refer to the Safety Data Sheet for first aid measures.
- Consult local or state authorities with regard to recommended method of disposal.

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- Water
- Stabilizers
- Preservatives

STORAGE AND STABILITY

- Store at 20°C to 28°C (68°F to 82.4°F).
- iChemVELOCITY Wash Solution is stable for twelve months from the date of manufacture.
- Open expiration: Unused solutions that remain in the original capped container are stable until the expiration date.
- Do not freeze
- Protect the solutions against light.
- Do not use solutions after the expiration date printed on the bottle label.

MATERIALS PROVIDED

- 2 x 7L iChemVELOCITY Wash Solution
- 2 x 1 µm Filter

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- iChemVELOCITY Wash Solution Tubing Assembly

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- iChemVELOCITY or DxU 810c Series of Automated Urine Chemistry Analyzers
- Personal Protective Equipment

PROCEDURE

Follow the directions in the iChemVELOCITY or DxU 810c Instructions for Use. The iChemVELOCITY or DxU 810c will automatically notify the operator when the Wash Solution liquid level is below minimal requirements.

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Rx Only in the U.S.A.

INTENDED USE

Iris System Cleanser is intended for *in vitro* diagnostic use with iQ200 or DxU 840/850m Series of Automated Urine Microscopy Analyzers and iChemVELOCITY or DxU 810c Automated Urine Chemistry Analyzer.

INTENDED USER

The device is intended to be used by a laboratory professional.

PRINCIPLE

Iris System Cleanser is a ready-for-use solution specifically developed as a wash solution for the iQ200 or DxU 840/850m Series of Automated Urine Microscopy Analyzers. For the iChemVELOCITY and DxU 810c Automated Urine Chemistry Analyzers, the Iris System Cleanser is used as a cleaning solution. A wash cycle must be performed at the beginning of each work day to prevent accumulation of residue in the fluidic system, sample lines and flow cell.

WARNING AND PRECAUTIONS

- Use only Iris System Cleanser for the Wash Cycle. Other similar solutions of hypochlorite should not be used because of the risk of particulate contaminants, trace oils or releasing compounds that could cause blockage or damage to the optical window.
- This material is harmful to eyes, skin and if inhaled or injested. Refer to the Safety Data Sheet for first aid measures.
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Do not reuse dispensed System Cleanser.
- Wear appropriate personal protective equipment when handling.
- Systemic allergic reactions are possible for sensitive individuals.
- Consult local or state authorities with regard to recommended method of disposal.
- Do not mix with other products- can release dangerous gas (chlorine).
- Product should present as a pale yellow, clear liquid; bleach odor. Product that differs from these attributes should be discarded.

GHS HAZARD CLASSIFICATION

Iris System Cleanser — DANGER

H314	Causes severe skin burns and eye damage.
H400	Very toxic to aquatic life.
H411	Toxic to aquatic life with long lasting effects.
P273	Avoid release to the environment.
P280	Wear protective gloves, protective clothing and eye/face protection.
P301+P330+P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
P303+P361+P353	IF ON SKIN (or hair): Rinse skin with water.

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P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor/physician.

Sodium Hydroxide 0.1 - 1%
Sodium Hypochlorite 1 - 5%

SDS Safety Data Sheet is available at beckmancoulter.com/techdocs

REACTIVE INGREDIENTS

Sodium hypochlorite, 4%

NON-REACTIVE INGREDIENTS

Purified Solvent

STORAGE AND STABILITY

- Store at 20°C to 28°C (68°F to 82.4°F).
- Iris System Cleanser is stable for eighteen months from date of manufacture.
- Open expiration: Unused solution that remains in the original capped container is stable until the expiration date.
- Do not freeze.
- Do not store in direct sunlight.
- Do not use solution after the expiration date printed on the bottle label.

MATERIALS PROVIDED

- 4 X 425 mL Iris System Cleanser
- 4 Flip-Top Caps

MATERIALS REQUIRED, BUT NOT PROVIDED

- Approved tubes listed in the iQ200 or DxU 840/850m Instructions for Use and the iChemVELOCITY or DxU 810c Instructions for Use. Do not use preservative tubes.
- Personal Protective Equipment.

PROCEDURE

- Follow the directions in the instrument IFU (Instructions for Use).

TECHNICAL SUPPORT

For additional information, or if damaged product is received, call Beckman Coulter Customer Service at 800-854-3633 (USA or Canada) or contact your local Beckman Coulter Representative.

A Glossary of Symbols (PN C27648) is available at beckmancoulter.com/techdocs.

NOTICE TO USER

For a patient/user/third party in the European Union and in countries with identical regulatory regime (*Regulation 2017/746/EU on In vitro Diagnostic Medical Devices*); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

May be covered by one or more pat. -see www.beckmancoulter.com/patents.

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Rx Only in the U.S.A.

INTENDED USE

Iris Diluent is intended for *in vitro* Diagnostic Use with the iQ200 or DxU 840/850m Series of Automated Urine Microscopy Analyzers.

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PRINCIPLE

Iris Diluent is a particle free aqueous salt solution which is used for diluting urine and body fluid samples. It is also used to rinse the system to prevent carryover of system cleanser and samples.

WARNING AND PRECAUTIONS

- Dispensed material should not be re-used.
- Iris Diluent has been determined to be nonhazardous under the guidelines issued by OSHA in 29 CFR 1910. 1200(d).
- Appropriate precautions in the collection, handling, storage and disposal of specimens should be exercised accordingly.
- Wear appropriate personal protective equipment when handling.
- Product should present as a clear, colorless, odorless liquid. Product that differs from these attributes should be discarded as they may have deteriorated.
- Consult local and state authorities with regard to recommended method of disposal

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- Water
- Salts
- Preservatives

STORAGE AND STABILITY

- Store at 20°C to 28°C (68°F to 82.4°F).
- Iris Diluent is stable for twelve months from date of manufacture.
- Open expiration: Unused solutions that remain in the original capped container are stable until the expiration date.
- Do not freeze.
- Do not store in direct sunlight.
- Do not use solutions after the expiration date printed on the bottle label.

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- Approved tubes listed in the iQ200 or DxU 840/850m Instructions for Use. Do not use preservative tubes.
- Personal Protective Equipment.

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PROCEDURE

Follow the directions in the iQ200 or DxU 840/850m Instructions for Use.

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