GRA-4214 Rev02, DCR2025-0372
Date of issue: 2024.06.16

Conger Plus

HYDROPHILIC GUIDEWIRE

INSTRUCTIONS FOR USE	ENGLISH
INSTRUCCIONES DE USO	ESPAÑOL
ISTRUZIONI PER L'USO	ITALIANO
INSTRUCTIONS D'UTILISATION	FRANÇAIS
GEBRAUCHSANWEISUNG	DEUTSCH
INSTRUKCJA UŻYCIA	POLSKI
INSTRUÇÕES PARA UTILIZAÇÃO	PORTUGUÊS
KULLANIM TALİMATLARI	TÜRKÇE
GEBRUIKSAANWIJZING	NEDERLANDS
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	ΕΛΛΗΝΙΚΑ
ИНСТРУКЦИИ ЗА УПОТРЕБА	БЪЛГАРСКИ
HASZNÁLATI UTASÍTÁS	MAGYAR
INSTRUCȚIUNI DE UTILIZARE	ROMÂNĂ
ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ	РУССКИЙ
UPUTE ZA UPORABU	HRVATSKI
NAVODILA ZA UPORABO	SLOVENŠČINA
BRUGSANVISNING	DANSK
KÄYTTÖOHJEET	SUOMI
BRUKSANVISNING	NORSK
BRUKSANVISNING	SVENSKA
NÁVOD NA POUŽITIE	SLOVENSKY
NÁVOD K POUŽITÍ	ČEŠTINA
KASUTUSJUHEND	EESTI
LIETOŠANAS	LATVIEŠ
NAUDOJIMO	ANGLŲ
ISTRUZZJONIJIE	INGLIŻ

BrosMed medical Co., Ltd

CE 0197

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Conger Plus HYDROPHILIC GUIDEWIRES

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGH OUT THESE INSTRUCTIONS. FAILURE TO DO MAY RESULT IN COMPLICATIONS.

1.0 Description:

Conger Plus Hydrophilic Guidewire is constructed from a high quality, steerable metallic alloy core wire with a polymer coating utilizing a sophisticated construction process to improve crossability of the guidewire. The hydrophilic coating is applied over the radiopaque polymer jacket to reduce the friction. Conger Plus Hydrophilic Guidewire is supplied sterile, non-pyrogenic and is intended for single use only.

1.1 Performance characteristics:

Outer diameter: 0.018 inch – 0.038 inch
Effective length: 45cm – 300 cm
Distal tip shape: angled, straight
Hardness: standard, stiff

2.0 Intended Use

The Conger Plus Hydrophilic Guidewire is used to facilitate the subsequent introductions of intravascular devices into peripheral vasculature during interventional procedures.

2.1 Indication

Conger Plus Hydrophilic Guidewire is used for peripheral vascular diseases that require interventional procedures.

2.2 Contraindication

Not intended for use in the heart or the central circulatory system.

2.3 Intended Patient Populations

This device is suitable for any patients with above disease or medical condition except for whom with contraindications.

2.4 Intended User

Intended users are the competent physicians who have been trained the interventional procedures.

2.5 Clinical benefit to be expected

High success rate of directing a catheter to the target location during interventional procedures.

2.6 Compatibility

The lumen size of the interventional device should be larger than the maximum outer diameter of the guidewire. The choice of guidewire size should be based on the recommendations of the interventional device label.

3.0 Warnings:

Carefully read all instructions prior to use. Failure to abide by the following warnings might result in damage to the vessel, shearing of the product and release of plastic fragments from the guidewire such pieces or fragments from the guidewire may have to be removed from the vessel.

● Do not manipulate or withdraw the guidewire through a metal entry needle or a metal dilator. Manipulation and/or withdrawal through a metal entry needle or a metal dilator may result in destruction and/or separation of the outer polyurethane coating requiring retrieval. A plastic entry needle is recommended when using this guidewire for initial placement.

● Do not use the guidewire with devices which contain metal parts such as atherectomy catheters, laser catheter, or metal introduction devices as they may cause the guidewire plastic coating to shear and/or sever the guidewire.

● Do not attempt to use the guidewire if it has been bent, kinked, or damaged. Use of a damaged guidewire may result in damage to the vessel or the release of fragments into the vessel.

● Do not reshape the guidewire by any means. Attempting to reshape the guidewire may cause damage, resulting in the release of the guidewire fragments into the vessel.

● When exchanging or withdrawing a catheter over the guidewire, secure and maintain the guidewire in place under fluoroscopy to avoid unexpected guidewire advancement; otherwise the guidewire's tip may damage the vessel.

● A retrieving device, such as a gripper or basket forceps, can only be used after the guidewire has been removed from the vessel. Using a retrieving device while the guidewire is in the vessel may cause the guidewire to break.

● Manipulate the guidewire slowly and carefully in the vessel while confirming the behavior, and location of the guidewire's tip under fluoroscopy. Excessive manipulation of the guidewire without fluoroscopic confirmation may result in vessel perforation. If any resistance is felt or the tip's behavior and/or location seems improper, stop manipulating the guidewire and/or the catheter, and determine the cause by fluoroscopy. Failure to exercise proper caution may result in bending, kinking and separation of the guidewire's tip, damage to the catheter, or damage to the vessel.

● Use the guidewire prior to the " Use-by date" (Expiration Date) specified on the package.

4.0Adverse Events:

Potential complications which may result from the improper use of a guidewire include, but are not limited to:

- Thrombus
- Air embolism
- Arterial or venous vessel wall damage
- Hematoma
- Infection
- Vessel spasm
- Intimal rupture
- Sterile inflammation or granulomas at the access site·
- Embolic stroke
- Additional surgical intervention

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5.0 Precautions:

● The guidewire is for single use only. Do not reuse reprocess or resterilize. Inspect the guidewire for damage prior to use, do not use any unit if package is opened or damaged. This device is sterilized by ethylene oxide gas. It is recommended that a plastic torque device be used to handle the guidewire. Use of a metal torque device may damage the guidewire surface coating.

● The surface of the guidewire is not lubricous unless it is wet. Before taking it out of its holder and inserting it through a catheter, fill the holder and the catheter with heparinized physiological saline solution. Avoid the guidewire exposure to any substance which reacts with coating, such as water and alcohol.

● **Storage conditions**: Keep dry. Keep away from sunlight. Store within temperature limit.

● After use, dispose and discard the product and packaging in accordance with hospital, administrative and /or local government policy.

6.0 Instruction for Use:

● Remove the guidewire and the dispenser together from the package and place them in the sterile area.

● Fill the dispenser with heparinized physiological saline

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solution using a syringe to fully soak the guidewire and activate the hydrophilic coating on the guidewire surface. [Note]: If the guidewire is not fully soaked, it may be damaged or difficult to take out the guidewire when it is extracted from the hoop dispenser.

● Remove the guidewire from the dispenser and inspect the guidewire prior to use, to verify that it is lubricated. If the guidewire cannot be easily removed from the dispenser, inject more heparinized physiological saline solution into the dispenser, and try again.

● Prior to use, Soak the target device (such as catheter) in the heparinized physiological saline solution to ensure smooth movement of the guidewire within it.

● Insert the guidewire into the cavity of the target device, and under X-ray fluoroscopy, send the guidewire into the blood vessel, and push the guidewire to guide the target device to the predetermined position.

[Note]: When inserting or exchanging other instruments again, it is necessary to clean the excessive blood on the guidewire surface with sterile gauze fully wetted with heparinized physiological saline solution, and soak the guidewire again with heparinized physiological saline solution.

[Note]: To prevent contrast agent crystallization / clotting, a continuous saline flush should be maintained between the interventional device and the guide wire during the procedure. The size of the syringe used to flush the lumen should be adapted to the length and diameter of the interventional device.

● After the operation, the puncture point should be pressed to stop bleeding, and pressure bandaged.

7.0 Disclaimer of Warranty:

Although the guidewire, hereafter referred to as "product", has been manufactured under carefully controlled conditions, BrosMed medical Co., Ltd and its affiliates have no control over conditions under which this product is used. BrosMed medical Co., Ltd and its affiliates, therefore, disclaims all warranties, both expressed and implied, with respect to the product,including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. BrosMed medical Co., Ltd and its affiliates shall not be liable to any personal or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind BrosMed medical Co., Ltd and its affiliates to any representation or warranty with respect to the product.

The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law.

If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

BrosMed Medical Co., Ltd
15th Building, SMEs Venture Park,
SongShan Lake Hi-Tech industrial Development Zone,
Dongguan , Guangdong, 523808, China
Office: +86 (769) 2289 2018
E-mail: sales@brosmed.com
www.brosmed.com

EU REP BrosMed Medical B.V.
Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Tel:+31 850 140 901
Mail:cs@brosmed.com

Explanation of Symbols:

Description	symbol
Catalog Number	REF
Lot Number	LOT
Contents (numeral represents quantity of units inside)	#
Sterilized using ethylene oxide	STERILE EO
Single sterile barrier system with protective packaging inside	
Do not re-use	
Caution	
Consult instructions for use	
Do not resterilize	
Do not use if package is damaged	
Manufacturer	
Authorized representative	EU REP
Keep dry	
Keep away from sunlight	
Temperature limit	10 30
Date of Manufacture	
Use-by date	
Medical Device	MD
Unique Device Identifier	UDI
CE Mark	CE 0197

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BrosMed Medical Co., Ltd

15th Building, SMEs Venture Park,
SongShan Lake Hi-Tech industrial Development Zone,
Dongguan , Guangdong, 523808, China
Office: +86 (769) 2289 2018
E-mail: sales@brosmed.com
www.brosmed.com

EU REP BrosMed Medical B.V.
Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Tel:+31 850 140 901
Mail:cs@brosmed.com

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Khons GUIDEWIRE

CAREFULLY READ ALL INSTRUCTIONS BEFORE USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

1.0 Device description

The Khons guidewire is a steerable guidewire available in a nominal outer diameter of 0.018 inch (0.46mm) and 0.014 inch (0.36mm) and nominal lengths of 80, 110, 150, 190, 200, 260, 310, 380, 400, and 450 cm. The tip flexibility of the guidewire has three specifications: soft, standard, and stiff; and the tip shape of the guidewire has two specifications: straight and angled. The distal shaft of the guidewire is welded with a radiopaque coil, wrapped in a polymer jacket, and coated with a hydrophilic coating on the surface. The proximal shaft of the guidewire is coated with polytetrafluoroethylene (PTFE). Khons guidewire is supplied sterile, non-pyrogenic, and intended for single use only.

1.1Performance characteristics

Outer diameter: 0.014 inch, 0.018 inch

Effective length: 80cm–450cm

Distal tip shape: angled, straight

Distal tip flexibility: soft, standard, stiff

2.0 How supplied

2.1 Content:

One guidewire with a torque device for each pouch.

2.2 Sterile:

Sterilized with ethylene oxide gas. Non-pyrogenic.

2.3 Storage Conditions:

Keep dry. Keep away from sunlight. Store within the temperature limit from 10°C to 30°C.

3.0 Intended use

The Khons Guidewire is used to facilitate the subsequent introductions of intravascular devices into peripheral vasculature during interventional procedures.

3.1 Indication

Khons Guidewire is used for peripheral vascular diseases that require interventional procedures.

3.2 Intended Patient Populations

This device is suitable for any patients undergoing peripheral vascular procedures.

3.3 Intended User

Intended users are the competent physicians who have been trained in the interventional procedures.

3.4 Clinical benefit to be expected

High success rate of directing intravascular devices to the target location during interventional procedures.

3.5 Contraindication

Not intended for use in the coronary arteries or the neurovasculature.

4.0 Warnings

- When the guidewire is in the body, it should be manipulated under fluoroscopy. Do not attempt to move the wire without observing the resultant tip response.
- Do not advance or withdraw the guidewire when excessive resistance is met until the cause of resistance is determined; otherwise, the guidewire may be bent, damaged, etc., which may lead to damage to the catheter and blood vessels
- When the guidewire torques inside the blood vessel, please turn it alternately clockwise and counterclockwise.
- The distal region of the guidewire should be hydrated by flushing saline prior to use. Otherwise, it may cause coating peeling and damage, which hinders the smooth movement of the guidewire.
- Avoid repeated bending at the same point to prevent damage to or separation of the guidewire.
- Do not re-sterilize, reprocess, and/or reuse the device, or it may damage the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness, or death. Re-sterilization, reprocessing, and/or reuse may also pose a risk of device contamination and/or lead to patient infection or cross-infection, including but not limited to the transmission of infectious diseases from one patient to another. Device contamination may cause injury, illness, or death to patients.

5.0 Precautions

- It is important to follow the instructions for use prior to using this product.
- The guidewire is intended for single use only. Do not re-sterilize and/or reuse the device.
- Do not use if the packaging is damaged.
- Please use the guidewire within the validity period indicated on the label. If the guidewire exceeds its expiration date, please dispose of it according to relevant procedures.
- Verify that the inner diameter of any diagnostic or therapeutic device to be used with the guidewire is compatible with the outer diameter of the guidewire prior to use.
- Be cautious when manipulating, advancing, and withdrawing these devices through the insertion tool, microcatheter, or other devices with sharp edges, or through tortuous or calcified blood vessels.
- The guidewire should be used immediately after opening. After use, dispose of it in accordance with hospital and/or local government policy.
- When medication or other intervention devices are used in combination with this device, to prevent damage to the guidewire, it is necessary to fully understand the characteristics of the medication or device, and confirm that it does not cause damage to the guidewire.
- Please advance or withdraw the guidewire slowly.
- The torque device is included to aid in the use of the guidewire. It is not intended to enter the patient's body at any time.
- Exercise care in handling the guidewire to reduce the chance of accidental damage. Do not expose the guidewire surface to organic solvents such as alcohol or medications, which might damage the guidewire coatings and/or cause the guidewire to lose lubricity.

6.0 Adverse events

Potential complications include, but are not limited to:

- Air embolism/Thromboembolism
- Allergic reaction (to contrast, device, or other)
- Amputation
- Arteriovenous (AV) fistula
- Hematoma/Hemorrhage
- Sepsis/Infection
- Vessel occlusion
- Vessel spasm, trauma, perforation, dissection
- Myocardial infarction or ischemia
- Embolization (plaque, thrombus, device, tissue, or other)
- Death
- Stroke/Cerebral vascular accident (CVA)/Transient ischemic attack
- Wire entrapment/Entanglement
- Foreign body/Wire fracture

Notice: any serious incident that has occurred ***in relation to the device*** should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

7.0 Instruction for use:

- Prior to use, please carefully check and confirm the sterilization package is not damaged and the guidewire is within its validity period. Remove the guidewire and the dispenser from the package.

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● Before removing the guidewire from the dispenser, hydrate the hydrophilic segment by flushing heparinized physiological saline through the dispenser using a syringe attached to the dispenser hub.
● Remove the proximal end of the guidewire from the dispenser and slowly push the guidewire. Repeat injection of heparinized saline if difficult removal of the guidewire occurs. After removing the guidewire from the dispenser, carefully inspect the wire prior to use to verify that it is lubricated.
● Flush any diagnostic or therapeutic device (such as a catheter) to be used with the guidewire with heparinized saline solution prior to use to ensure the guidewire can smoothly move within the target device.
● Carefully insert the distal section of the guidewire into the catheter and advance. Advance the guidewire until the distal tip is near the distal end of the catheter.
● If necessary, the guidewire can be used together with the provided torque device. Slip the torque device over the proximal end of the guidewire to the desired location. Ensure the torque device is in place by tightening the rotating knob. The torque device may be repositioned by loosening and retightening the rotating knob.
● During navigation in the vasculature, loosen the hemostatic connector, advance, and rotate the guidewire by rotating the torque device in either direction to facilitate vessel selection.
● When reinserting or exchanging other devices, it is necessary to clean the surface of the guidewire with sterile gauze that has been fully hydrated in heparinized saline solution and then hydrate the guidewire again with heparinized saline solution.
● Inspect the guidewire integrity after the operation.
● After the operation, the puncture point should be compressed to stop bleeding, and a tourniquet should be applied. After use, dispose of the guidewire in accordance with hospital and/or local government policy.

8.0 Disclaimer of warranty

Although the guidewire (hereafter referred to as the 'product') has been manufactured under carefully controlled conditions, BrosMed Medical Co., Ltd. and its affiliates have no control over the conditions in which the product is used. Therefore, BrosMed Medical Co., Ltd. and its affiliates disclaim all warranties, both express and implied, with respect to the product. This includes any implied warranty of merchantability or fitness for a particular purpose. BrosMed Medical Co., Ltd. and its affiliates shall not be liable for any medical expenses or any direct, incidental or consequential damages caused by the use, defect, failure or malfunction of the product, whether a claim for such damages is based on a warranty, contract or tort. No person has any authority to bind BrosMed Medical Co., Ltd. and its affiliates to any representation or warranty with respect to the product.
The exclusion and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court or competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

Manufacturer:
BrosMed Medical Co., Ltd.
15th building, SMEs Venture Park, Songshan Lake Hi-Tech Industrial Development Zone, Dongguan, Guangdong, 523808, China
Office: +86(769)2289 2018
E-mail: sales@brosmed.com
www.brosmed.com

EU REP EU Authorized Representative:
BrosMed Medical B.V.
Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Tel: + 31 850140901
E-mail: cs@brosmed.com

9.0 Explanation of symbols

Description	Symbol
Manufacturer	
Catalog number	REF
Batch code	LOT
Use-by date	
Date of manufacture	
Contents (numeral represents quantity of units inside)	#
Keep away from sunlight	
Keep dry	
Temperature limit	10°C 30°C
Medical device	MD
Unique device identifier	UDI
Do not re-use	
Caution	
Consult instructions for use	
Sterilized using ethylene oxide	STERILE EO
Single sterile barrier system with protective packaging inside	
Do not use if package is damaged	
Do not resterilize	
CE mark	CE 0197
Authorized representative	EU REP

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
Pregarde Guiding Sheath Instructions for Use English


Please read the instructions for use carefully before use.
Follow all items in this instruction, including the warnings and precautions, otherwise complications may occur.

1.0 Product Description

Pregarde Guiding Sheath consists of a sheath introducer and a dilator. The sheath introducer is mainly composed of a sheath and a haemostasis valve. The haemostasis valve has two types: Cross-Cut and Y-Shape. Both valve types are equipped with a three-way stopcock connected by a sidearm, as shown below. The sheath is composed of an outer layer, an intermediate stainless steel spring coil and an inner tube, with an outer surface coated with a hydrophilic coating and a radiopaque marking band at the distal end.

Y-Shape Valve (YSV)

Cross-Cut Valve (CCV)

1.1 Performance Characteristics:

Inner diameter: 4F, 5F, 6F, 7F, 8F
Effective length: 45cm, 65cm, 90cm, 110cm
Distal tip shape: ST, HS, LM, RDC, MP, CO

2.0 How Supplied

Contents:	One (1) sheath introducer and One (1) dilator.
Sterility:	Sterilized with ethylene oxide gas. Non-pyrogenic.
Transport condition:	Protected from heavy pressure, direct sunlight, rain, and snow.
Storage:	Keep away from sunlight, keep dry, and store within the temperature limit (10-30℃).

3.0 Intended Use/Indication

The Pregarde Guiding Sheath is intended to introduce therapeutic or diagnostic devices into the vasculature.

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3.4 Clinical Benefit

The intended clinical benefit of Guiding Sheath is to help make the introduction easier and protect the puncture site.

3.5 Compatibility

Before use, make sure the sheath size (Fr.) is appropriate for the access vessel and the interventional/ diagnostic device to be used.

4.0 Warnings

☆ This Guiding Sheath is for single use only. Do not repeatedly sterilize or use it.
☆ Do not use if the hydrophilic coating comes off.
☆ Do not use high-pressure syringes to inject contrast material.
☆ If you feel a significant increase in resistance when inserting the guide wire or sheath, discontinue insertion and withdraw the guide wire or sheath.
☆ Before removing the Guiding Sheath, if there is distortion of the vessels through the Guiding Sheath, it is necessary to insert the dilator to avoid possible damage.

5.0 Precautions

☆ Do not use if the product packaging has been opened or damaged before use.
☆ Check the Guiding Sheath before use, do not use it if it is damaged.
☆ Before use, ensure that all devices or catheters used with the product are suitable for the size of the sheath introducer and dilator, and must be able to move freely in the valve and sheath.
☆ The Guiding Sheath should be aspirated through a three-way stopcock, cleaned, and flushed with heparin saline before the device is removed through the Guiding Sheath or inserted.
☆ When inflating the balloon at or close to the tip of the Guiding Sheath, be sure not to inflate it inside the distal end of the Guiding Sheath.
☆ This product is valid for three years. Please use it before the expiry date indicated on the package.
☆ It should be safely disposed of after one-time use to avoid infection.

6.0 Adverse Events

The adverse events may occur, including but not limited to the following:
☆ Infection
☆ Aeroembolism
☆ Extravasation
☆ Hematoma
☆ Thrombosis
☆ Inner membrane rupture
☆ Distal embolization
☆ Arteriospasm
Notice: Any serious incident that has occurred ***in relation to the device*** should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

7.0 Instructions for Use

7.1 Upon removal from the package, ensure the inner diameter (ID) of the introducer is appropriate for the maximum diameter of the instrument or catheter to be introduced.
7.2 Remove the product from the packaging.
7.3 Fill the sheath introducer completely with heparin saline through the three-way stopcock to flush the sheath introducer. After cleaning, rotate the three-way stopcock to the closed position of the sidearm to prevent backflow of blood from the sidearm.
7.4 Flush the dilator with heparin saline.
7.5 Insert the dilator completely into the sheath introducer so that the dilator hub fits closely with the cap. For sheath introducers with Y-Shape Valve (YSV), tighten the Y-Shape Valve (YSV) around the dilator.
Note: Insert the dilator into the center of the haemostasis valve of the sheath introducer. Failure to insert the dilator in the center of the valve may cause valve damage, resulting in blood leakage.
7.5 The outer surface of the sheath was humidified with heparin saline to activate the hydrophilic coating.
Note: For best results, maintain a wetted condition of the device during placement.
7.6 Using the standard Seldinger technique, access the target vessel with the appropriate needle.

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7.7 Insert guide wire into vessel through the needle, then remove needle, leaving guide wire in place.
7.8 Insert the dilator and sheath introducer assembly together through the guide wire into the vessel.
Note: To prevent damage to the tip of the sheath introducer, do not withdraw the dilator until the dilator and sheath introducer assembly have been inserted into the vessel.
7.9 Release the dilator from the cap by gently pressing it. For sheath introducers with Y-Shape Valves (YSV), it is necessary to loosen the Y-Shape Valve (YSV).
7.10 Slowly pull out the guide wire and dilator, leaving the sheath introducer in the vessel.
7.11 Slowly aspirate with a syringe through the three-way stopcock to remove residual air in the sheath introducer. A suitable heparinized solution injected from a three-way stopcock after aspiration may be considered to prevent thrombosis.
7.12 The selected catheter or device is inserted into the center of the haemostasis valve and then passed through the sheath into the vessel.
Note: When withdrawing any endovascular device, aspirate through a three-way stopcock to remove fibrin or clots that may be present inside or at the tip of the sheath.
7.13 Aspirate with a syringe in the sheath introducer before withdrawing, and then slowly withdraw the sheath introducer to properly press the puncture site.
Note: If you encounter difficulties in withdrawing the sheath introducer, please investigate the cause first, and continue to withdraw after the cause is confirmed. Excessive withdrawal force may result in product separation or vascular injury.
Note: After use, dispose of and discard the product and packaging in accordance with hospital, administrative, and/or local government policy.

8.0 Disclaimer

The statement or description in the publications of BrosMed Medical Co., Ltd. (including these instructions for use) only serves as a general description of this product and does not constitute any explicit guarantee. BrosMed Medical Co., Ltd. accepts no liability for any accident, direct or indirect damage caused by the misuse or incorrect operation of this product.

Manufacturer:
BrosMed Medical Co., Ltd.
15th building, SMEs Venture Park, Songshan Lake
Hi-Tech Industrial Development Zone, Dongguan,
Guangdong, 523808, China
Office: +86 (769) 2289 2018
E-mail: sales@brosmed.com

EU Authorized Representative / International Sales Office:
BrosMed Medical B.V.
Mgr. Buckxstraat 8, 6134 AP Sittard, The Netherlands
Office: +31 850 140 901
E-mail: cs@brosmed.com

Explanation of Symbols:

Description	Symbol
Catalogue number	REF
Batch code	LOT
Date of manufacture	
Use-by date	
Contents (numeral represents quantity of units inside)	n
Manufacturer	
Authorized representative	EU REP
CE Mark	CE 0197
Medical device	MD
Unique device identifier	UDI
Sterilized using ethylene oxide	STERILE EO
Single sterile barrier system with protective packaging inside	
Do not re-use	
Do not use if the packaging is damaged	
Do not resterilize	
Caution	
Consult instructions for use	
Keep away from sunlight	
Keep dry	
Temperature limit	30℃ 10℃

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